

medisana®

PRODUCT RANGE 2022



HEALTH CONTROL THERAPY & HEALTHY HOME
WELLNESS & WARMTH SPORT BODY CARE

HEALTH CONTROL

Blood pressure monitor	3
Thermometer	8
Pulse oximeter	11
Body analysis scale	11
Personal scale	14
Blood glucose monitor	18
Particle filtering half mask	18

THERAPY & HEALTHY HOME

Inhaler	19
Lice comb	20
Hygrometer	20
Electro therapy	22
Air humidifier	22
Aroma diffuser	24
Light therapy	25

SPORT

Body trainer & Activity tracker	26
---------------------------------	----

WELLNESS & WARMTH

Massage seat cover	27
Massage mat	29
Neck massager	29
Massage cushion	31
Foot massage	32
Hand massage	35
Massage gun	36
Warmth & Heat pad	37
Heating blanket	40
Heating underblanket	41
Foot warmer	43
Warmth for the office	45
Outdoor head pad	46
Foot spa	47
Bath spa	49
Massage chair	50

BODY CARE

Hand & Foot care	51
Face & Skin care	52
Cellulite massager	53
Cosmetic mirror	54

ACCESSORIES

Accessories	56
-------------	----

NEW



BU 584 CONNECT

Upper arm blood pressure monitor

Connect device

- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*
- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- Motion detector for accurate measurement results
- Blue backlight

Product features

- Elegant classic design
- Indication of application errors
- Blood pressure classification – measurements classified using a traffic-light scale
- Irregular heart beat display
- 120 memory spaces for each of 2 users
- Big LCD display (86×78 mm) with blue backlight and extra large digits (40 mm)
- Incl. batteries (4×1.5 V AAA)
- Incl. storage bag
- Certified medical device



Size: approx. 13.7 × 9.85 × 5.25 cm
Weight: approx. 0.25 kg (without batteries)
Measuring method: Oscillometric
Measuring range: Blood pressure: 40–280 mmHg
 Pulse: 40–199 beats/minute
Batteries: 4×1.5 V (AAA)
Article No.: 51584
EAN: 40 15588 51584 2
Warranty: 3 years

NEW



BU 586 VOICE

Upper arm blood pressure monitor

With voice output

- Voice output in selectable languages: DE, GB, NL, FR, IT, TR
- Elegant classic design incl. big LCD display with extra large digits (40 mm)
- Motion detector to ensure accurate measurement results
- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff

Product features

- Indication of application errors
- Blood pressure classification – measurements classified using a traffic-light scale
- Irregular heartbeat display
- 120 memory slots for each of 2 users
- Extra large upper arm cuff (arm circumferences 23 – 43 cm)
- Incl. storage bag
- Incl. batteries (4 × AAA)
- Certified medical device



Size: approx. 13.7 × 9.85 × 5.25 cm
Weight: approx. 0.25 kg (without batteries)
Measuring method: Oscillometric
Measuring range: Blood pressure: 40–280 mmHg
 Pulse: 40–199 beats/minute
Batteries: 4×1.5 V (AAA)
Article No.: 51586
EAN: 40 15588 51586 6
Warranty: 3 years

Blood pressure monitor



BU 565

Upper arm blood pressure monitor

Design in pocket size

- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- Extra large and easy-to-read display
- Accurate blood pressure measurement on the upper arm with blood pressure classification
- Blue backlight

Product features

- Big LCD display: 7.3×6 cm with backlight
- 120 memory slots for each of 2 users
- Calculates average value of data of the last 3 results
- With indicator of irregular heartbeat
- Classification of measurements with traffic-light scale
- Certified medical device



Size: approx. 9.26 × 12.95 × 3.03 cm
Weight: approx. 0.188 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0–299 mmHg
 Pulse: 40–180 beats/minute
Batteries: 4×1.5 V (AAA)
Article No.: 51202
EAN: 40 15588 51202 5
Warranty: 3 years

* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer.
 Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.2

Blood pressure monitor



BU 535 VOICE (south)

Upper arm blood pressure monitor

With voice output

- Accurate blood pressure measurement on the upper arm
- Voice output in selectable languages: GB, FR, RU ES, GR, AR (Arabic)
- With indicator of irregular heartbeat
- 120 memory slots for each of 2 users
- Very large digits for easy reading

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Upper arm cuff for 22 – 36 cm arm circumference
- Incl. storage bag
- Certified medical device



Size: approx. 15 × 9.5 × 4.1 cm
Weight: approx. 0.23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40 – 260 mmHg
Pulse: 40 – 180 beats/minute
Batteries: 4 × 1.5 V (AAA)
Article No.: 51179
EAN: 40 15588 51179 0
Warranty: 3 years



BU 535 VOICE (north)

Upper arm blood pressure monitor

With voice output

- Accurate blood pressure measurement on the upper arm
- Voice output in selectable languages: DE, GB, NL, FR, IT, TR
- With indicator of irregular heartbeat
- 120 memory slots for each of 2 users
- Very large digits for easy reading

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Upper arm cuff for 22 – 36 cm arm circumference
- Incl. storage bag
- Certified medical device



Size: approx. 15 × 9.5 × 4.1 cm
Weight: approx. 0.23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40 – 260 mmHg
Pulse: 40 – 180 beats/minute
Batteries: 4 × 1.5 V (AAA)
Article No.: 51178
EAN: 40 15588 51178 3
Warranty: 3 years



BU 535

Upper arm blood pressure monitor

With XL-Display

- Accurate blood pressure measurement on the upper arm
- With indicator of irregular heartbeat
- 120 memory slots for each of 2 users
- Very large digits for easy reading

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Upper arm cuff for 22 – 36 cm arm circumference
- Incl. storage bag
- Certified medical device



Size: approx. 15 × 9.5 × 4.1 cm
Weight: approx. 0.23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40 – 260 mmHg
Pulse: 40 – 180 beats/minute
Batteries: 4 × 1.5 V (AAA)
Article No.: 51176
EAN: 40 15588 51176 9
Warranty: 3 years



BU 514

Upper arm blood pressure monitor

With inverse display

- Accurate blood pressure measurement on the upper arm
- With indicator of irregular heartbeat
- 90 memory slots for each of 2 users
- Very large digits for easy reading (up to 22 mm in height)

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of all stored data
- Adapter - or battery operation (adapter optional)
- Upper arm cuff for 22 – 36 cm arm circumference
- Incl. storage bag
- Certified medical device



Size: approx. 13×11×6 cm
Weight: approx. 0.3 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0 – 299 mmHg
Pulse: 40 – 199 beats/minute
Batteries: 4×1.5 V (AA)
Article No.: 51165
EAN: 40 15588 51165 2
Warranty: 3 years



BU 512

Upper arm blood pressure monitor

Handy upper arm cuff

- Accurate blood pressure measurement on the upper arm
- With indicator of irregular heartbeat
- 90 memory slots for each of 2 users
- Very large digits for easy reading (up to 22 mm in height)

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of all stored data
- Battery operation
- Upper arm cuff for 22 – 36 cm arm circumference
- Incl. storage bag
- Certified medical device



Size: approx. 13×11×6 cm
Weight: approx. 0.3 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40 – 299 mmHg
Pulse: 40 – 199 beats/minute
Batteries: 4×1.5 V (AA)
Article No.: 51162
EAN: 40 15588 51162 2
Warranty: 3 years



MTP PRO

Upper arm blood pressure monitor

Our top seller

- Accurate blood pressure measurement on the upper arm
- With indicator of irregular heartbeat
- 99 memory slots for each of 2 users
- Accurate measurement results with strongly fluctuating blood pressure through patented 3 MAM technology (weighted average of 3 subsequent measurements)
- 10 years of warranty

Product features

- Display items: systolic, diastolic, pulse, date, time
- Extra large, easy-to-read display
- Classification of measurements with traffic-light scale
- Calculates average value of all stored data
- One-button operation
- Battery operation (power adapter optionally available)
- Upper arm cuff for 22 - 32 cm arm circumference
- Incl. storage bag
- Incl. blood pressure pass
- Certified medical device



Size: approx. 12×16.3×6.6 cm
Weight: approx. 0.5 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30 – 280 mmHg
Pulse: 40 – 200 beats/minute
Batteries: 4×1.5 V (AA)
Article No.: 51090
EAN: 40 15588 51090 8
Warranty: 10 years

Blood pressure monitor



BU 570 CONNECT

Upper arm blood pressure monitor

Ultra slim design

- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*
- Extra large and easy-to-read display
- Accurate blood pressure measurement on the upper arm with blood pressure classification
- Blue backlight

Product features

- Big LCD display: 7.3×6 cm with backlight
- 120 memory slots for each of 2 users
- Calculates average value of data of the last 3 results
- With indicator of irregular heartbeat
- Classification of measurements with traffic-light scale
- Certified medical device



Size: approx. 9.26×12.95×3.03 cm
Weight: approx. 0.188 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0–299 mmHg
Pulse: 40–180 beats/minute
Batteries: 4×1.5 V (AAA)
Article No.: 51206
EAN: 40 15588 51206 3
Warranty: 3 years



BU 542 CONNECT

Upper arm blood pressure monitor

With inverse display

- Accurate blood pressure measurement on the upper arm
- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*
- With indicator of irregular heartbeat
- 250 memory slots for each of 2 users
- Very large digits for easy reading

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Extra large upper arm cuff (arm circumferences 22 – 42 cm)
- Incl. storage bag
- Certified medical device



Size: approx. 15×9.5×4.1 cm
Weight: approx. 0.25 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40–260 mmHg
Pulse: 40–180 beats/minute
Batteries: 4×1.5 V (AA)
Article No.: 51184
EAN: 40 15588 51184 4
Warranty: 3 years



BU 540 CONNECT

Upper arm blood pressure monitor

Top price-quality ratio

- Accurate blood pressure measurement on the upper arm
- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*
- With indicator of irregular heartbeat
- 250 memory slots for each of 2 users
- Very large digits for easy reading

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Extra large upper arm cuff (arm circumferences 22 – 42 cm)
- Incl. storage bag
- Certified medical device



Size: approx. 15×9.5×4.1 cm
Weight: approx. 0.23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40–260 mmHg
Pulse: 40–180 beats/minute
Batteries: 4×1.5 V (AA)
Article No.: 51182
EAN: 40 15588 51182 0
Warranty: 3 years



ecomed BU-92E

Upper arm blood pressure monitor

Top price-quality ratio

- Accurate blood pressure measurement on the upper arm
- With indicator of irregular heartbeat
- 60 memory slots for each of 2 users

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Cuff for 22-30 cm arm circumferences (optional: 30-42 cm, Art. 51126)
- Certified medical device



Size: approx. 13.6 × 10 × 6.5 cm
Weight: approx. 0.23 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 40–260 mmHg
Pulse: 40–180 beats/minute
Batteries: 4 × 1.5 V (AA)
Article No.: 23205
EAN: 40 15588 23205 3
Warranty: 2 years



BW 335

Wrist blood pressure monitor

Flat design

- Accurate blood pressure measurement on the wrist
- With indicator of irregular heartbeat
- 120 memory slots for each of 2 users
- Very large digits for easy reading
- Ultra slim design

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Wrist cuff for 14 – 19.5 cm wrist circumferences
- Incl. storage box
- Certified medical device



Size: approx. 8.9 × 6 × 3.1 cm
Weight: approx. 0.07 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 0 – 299 mmHg
Pulse: 40 – 180 beats/minute
Batteries: 2 × 1.5 V (AAA)
Article No.: 51077
EAN: 40 15588 51077 9
Warranty: 3 years

Blood pressure monitor



BW 320

Wrist blood pressure monitor

Top price-quality ratio

- Accurate blood pressure measurement on the wrist
- With indicator of irregular heartbeat
- 120 memory slots for each of 2 users
- Extra large, easy-to-read display

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Battery level indicator
- Extra large wrist cuff for 12.5 – 21.5 cm wrist circumferences
- Certified medical device



Size: approx. 7.5 × 6.4 × 2.9 cm
Weight: approx. 0.117 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30 – 280 mmHg
Pulse: 40 – 200 beats/minute
Batteries: 2 × 1.5 V (AAA)
Article No.: 51074
EAN: 40 15588 51074 8
Warranty: 3 years



BW 315

Wrist blood pressure monitor

Our bestseller

- Accurate blood pressure measurement on the wrist
- With indicator of irregular heartbeat
- 120 memory slots for each of 2 users
- Extra large, easy-to-read display

Product features

- Display items: systolic, diastolic, pulse, date, time
- Classification of measurements with traffic-light scale
- Calculates average value of data of the last 3 results
- Battery level indicator
- Extra large wrist cuff for 12.5 – 21.5 cm wrist circumferences
- Certified medical device



Size: approx. 6.3×6.3×2.9 cm
Weight: approx. 0.11 kg
Measuring method: Oscillometric
Measuring range: Blood pressure: 30 – 280 mmHg
Pulse: 40 – 200 beats/minute
Batteries: 2×1.5 V (AAA)
Article No.: 51072
EAN: 40 15588 51072 4
Warranty: 3 years



TM 750 CONNECT

Infrared multifunctional thermometer

Connect device

- Precise measurement of the body temperature in the ear and on the forehead in seconds
- Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100°C)
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock® Online*

Product features

- Automatic storage of the last 30 measurement results
- Visual fever alarm with colour change – green: normal temperature, red: raised temperature
- Acoustic signal at the beginning and end of temperature reading
- Display of date and time
- Incl. storage box
- Certified medical device



Size: approx. 13.4×4.2×4.7 cm
Weight: approx. 0.065 kg
Measuring range: Body: 34 – 43°C
Object: 0 – 100°C
Measuring accuracy: Forehead: ± 0.3°C (34 – 43°C)
Ear: ± 0.3°C (35.5 – 42°C)
Object: ≤30°C ±1.5°C / >30°C ±5%
Batteries: 2×1.5 V (AAA)
Article No.: 76145
EAN: 40 15588 76145 4
Warranty: 3 years



TM 750

Infrared multifunctional thermometer

Visual fever alarm with colour change

- Precise measurement of the body temperature in the ear and on the forehead in seconds
- Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100°C)

Product features

- Automatic storage of the last 30 measurement results
- Visual fever alarm with colour change – green: normal temperature, red: raised temperature
- Acoustic signal at the beginning and end of temperature reading
- Display of date and time
- Incl. storage box
- Certified medical device



Size: approx. 13.4×4.2×4.7 cm
Weight: approx. 0.065 kg
Measuring range: Body: 34 – 43°C
Object: 0 – 100°C
Measuring accuracy: Forehead: ± 0.3°C (34 – 43°C)
Ear: ± 0.3°C (35.5 – 42°C)
Object: ≤30°C ±1.5°C / >30°C ±5%
Batteries: 2×1.5 V (AAA)
Article No.: 76140
EAN: 40 15588 76140 9
Warranty: 3 years



TM 760

Infrared multifunctional thermometer

With LED tracking light

- Precise non-contact measurement of the body temperature
- Additional function: measurement of ambient temperature and surfaces (0 – 100°C)
- Visual and acoustic fever alarm with colour change – green: normal temperature, red: raised temperature
- Extra large, easy-to-read display

Product features

- LED Tracking light to indicate the correct point of measurement
- Acoustic signal at the beginning and end of temperature reading (can be switched off to avoid disturbance during sleep)
- Automatic storage of the last 30 measurement results
- Display of date and time
- Conversion from °C to °F
- Ergonomic and light-weight design
- Incl. storage bag
- Certified medical device



Size: approx. 15.8 × 4.2 × 4.3 cm
Weight: approx. 0.067 kg
Measuring range: Body: 34 – 43°C
 Object: 0 – 100°C
Measuring accuracy: Forehead: ± 0.3°C (34 – 42°C)
 Object: ± 1°C (0 – 100°C)
Batteries: 1 × CR2032
Article No.: 76121
EAN: 40 15588 76121 8
Warranty: 3 years



TM A79

Infrared body thermometer

Two modes – body and surfaces

- Precise and fast measurement of the body temperature on the forehead within seconds
- Additional function: measurement of ambient temperature, liquids and surfaces (0 – 100°C)
- Visual and acoustic fever alarm with colour change (green: normal temperature, yellow: slightly raised temperature, red: raised temperature)

Product features

- Extra large, easy-to-read LCD display
- Two modes integrated: Body mode (34 – 43°C) and surface mode (0 – 100°C)
- Automatic storage of the last 50 measurement results
- Automatic switch-off after 18 seconds
- Conversion from °C to °F
- Certified medical device



Size: approx. 13.8 × 9.5 × 4.0 cm
Weight: approx. 0.09 kg (without batteries)
Measuring range: Body: 34 – 43°C
 Object: 0 – 100°C
Measuring accuracy: Body: ± 0.3°C (34.0 – 42.0°C)
 Object: ± 2°C
Batteries: 2 × 1.5 V (AAA)
Article No.: 99663
EAN: 40 15588 99663 4
Warranty: 3 years



TM-65E

Infrared thermometer

Top price-quality ratio

- Precise and fast measurement of the body temperature
- Visual fever alarm
- Additional function: measurement of ambient temperature, liquids and surfaces

Product features

- Measurement in a few seconds
- Extra large, easy-to-read LCD display
- Automatic storage of the last 30 measurement results
- Automatic switch-off after 30 seconds
- Certified medical device



Size: approx. 12.3 × 5 × 5 cm
Weight: approx. 0.06 kg (without batteries)
Measuring range: Body: 34 – 43°C
 Object: 0 – 100°C
Measuring accuracy: Forehead: ± 0.3°C (34 – 42°C)
 Object: ± 2°C
Batteries: 2 × 1.5 V (AAA)
Article No.: 23400
EAN: 40 15588 23400 2
Warranty: 2 years

Thermometer



TM 705
Thermometer

With fever indicator

- Precise measurement of the body temperature: oral, axillary, rectal
- Flexible tip
- Measurement in 10 seconds
- Water resistant

Product features

- Acoustic fever alarm
- Acoustic signal at the beginning and end of temperature reading
- Automatic storage of the last measurement result
- Automatic switch-off after 10 minutes
- Digital LCD display
- Conversion from °C to °F
- Easy battery replacement
- Incl. storage box
- Certified medical device



Size: approx. 13.9×2.2×1.3 cm
Weight: approx. 0.013 kg
Measuring range: 32 – 42.9°C
Scaling: 0.1°C
Batteries: 1×1.5 V (AG3/LR41)
Article No.: 77042
EAN: 40 15588 77042 5
Warranty: 3 years



TM 700
Thermometer

Our gentle classic

- Precise measurement of the body temperature: oral, axillary, rectal
- Flexible tip
- Measurement in 10 seconds
- Water resistant

Product features

- Acoustic fever alarm
- Acoustic signal at the beginning and end of temperature reading
- Automatic storage of the last measurement result
- Automatic switch-off after 10 minutes
- Digital LCD display
- Conversion from °C to °F
- Easy battery replacement
- Incl. storage box
- Certified medical device



Size: approx. 13.9×2.2×1.3 cm
Weight: approx. 0.013 kg
Measuring range: 32 – 42.9°C
Scaling: 0.1°C
Batteries: 1×1.5 V (AG3/LR41)
Article No.: 77040
EAN: 40 15588 77040 1
Warranty: 3 years



FTC
Thermometer

Top price-quality ratio


- Precise measurement of the body temperature: oral, axillary, rectal
- Acoustic signal at the beginning and end of temperature reading
- Water resistant

Product features

- Acoustic fever alarm
- Automatic storage of the last measurement result
- Automatic switch-off after 8 minutes
- Digital LCD display
- Easy battery replacement
- Incl. storage box
- Certified medical device



Size: approx. 12.8×1×1.9 cm
Weight: approx. 0.0098 kg
Measuring range: Body 32 – 42°C
Scaling: 0.1°C
Batteries: 1×1.5 V (AG3/LR41)
Article No.: 77030
EAN: 40 15588 77030 2
Warranty: 3 years



PM 100 CONNECT

Pulse oximeter

Now as connect device

- Pulse oximeter to measure the saturation of oxygen in blood (SpO2%) and the heart rate (pulse)
- View measurements on OLED display
- Simple "one touch" operation – ideal for on the go
- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*

Product features

- Oxygen saturation (SpO2) and pulse rate can be displayed as bars or waves
- Automatic switch-off after 8 seconds
- Low battery display
- Incl. carrying strap
- Single measurements or long-term measurements selectable
- Certified medical device



Size: approx. 5.8×3.4×3.5 cm
Weight: approx. 0.055 kg
Measuring range: SpO2: 70% ~ 100%
Pulse: 30 – 235 beats/min
Batteries: 2×1.5 V (AAA)
Article No.: 79456
EAN: 40 15588 79456 8
Warranty: 3 years



PM 100

Pulse oximeter

Simple "one touch" operation


- Pulse oximeter to measure the saturation of oxygen in blood (SpO2 %) and the heart rate (pulse)
- View measurements on OLED display
- Display brightness adjustable in 10 levels

Product features

- Oxygen saturation (SpO2) and pulse rate can be displayed as bars or waves
- Display of the perfusion index PI (pulse strength)
- 7 different views
- Automatic switch-off after 8 seconds
- Low battery display
- Incl. carrying strap
- Certified medical device



Size: approx. 5.8×3.4×3.5 cm
Weight: approx. 0.055 kg
Measuring range: SpO2: 70% ~ 100%
Pulse: 30 – 235 beats/min
Batteries: 2×1.5 V (AAA)
Article No.: 79455
EAN: 40 15588 79455 1
Warranty: 3 years



BS 650 CONNECT

WiFi & Bluetooth Body analysis scale

With high quality ITO-electrodes

- Convenient WiFi or Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*
- Measurement of weight, body fat, body water, muscle mass and bone mass
- Heart rate measurement with every weight measurement (pulse)

Product features

- One-time connection of the scale to the WLAN network via the VitaDock+ App ensures automatic, continuous synchronisation of new measured values
- Determination of the BMI and BMR value
- Convenient "step-on" switch and automatic switch off
- High-quality and invisible ITO electrodes for higher conductivity of body analysis
- Large easy-to-read LCD display with white luminous digits
- Automatic detection of up to 8 users
- Switchable between kg, lb and st
- Body analysis not available for persons with medical implants (e.g. pacemakers)



Size: approx. 31×31×2.5 cm
Weight: approx. 1.80 kg
Display size: 7.25×4.2 cm
Weight display: up to 180 kg, 396 lb, 28st 4lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4×1.5 V (AAA)
Article No.: 40500
EAN: 40 15588 40500 6
Warranty: 3 years

*iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer.
Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.2

Body analysis scale



BS 600 CONNECT

**WiFi & Bluetooth
Body analysis scale**

Innovation in modern design

- Convenient WiFi or Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*
- Heart rate measurement with every weight measurement (pulse)

Product features

- One-time connection of the scale to the WiFi network via the VitaDock+ App ensures automatic, continuous synchronisation of new measured values
- Measurement of weight, body fat, body water, muscle mass and bone mass
- Determination of the BMI and BMR value
- Convenient "step-on" switch and automatic switch off
- Large easy-to-read LCD display with white luminous digits
- Automatic detection of up to 8 users
- Switchable between kg, lb and st
- Body analysis not available for persons with medical implants (e.g. pacemakers)



Size: approx. 31×31×2.5 cm
Weight: approx. 1.80 kg
Display size: 7.25×4.2 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4×1.5 V (AAA)
Article No.: 40501
EAN: 40 15588 40501 3
Warranty: 3 years



BS 550 CONNECT

Bamboo body analysis scale

Made from genuine bamboo

- Body analysis scale with high-quality bamboo surface for a pleasant feeling on the feet
- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock® Online*

Product features

- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Large, easy-to-read LCD display
- Convenient "step-on" switch and automatic switch-off
- 240 memory slots for measured values
- Body analysis not available for persons with medical implants (e.g. pacemakers)



Size: approx. 30×30×3.3 cm
Weight: approx. 1.2 kg
Display size: 7.25×4.2 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 0.2 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40550
EAN: 40 15588 40550 1
Warranty: 3 years



BS 450 CONNECT

Body analysis scale

Connect device with LED display

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock® Online*

Product features

- Ultra high-quality ITO electrodes
- Convenient "step-on" activation
- Invisible LED display that only lights up when in use
- Switchable between kg, lb and st
- 4 high-precision strain gauge sensors for accurate results
- Automatic recognition of up to 8 users
- 30 memory slots for each of the 8 users
- Body analysis not available for persons with medical implants (e.g. pacemakers)
- Sensor touch buttons



Size: approx. 31.9×31.9×2.3 cm
Weight: approx. 1.7 kg
Display size: 7.4×2.9 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4×1.5 V (AAA)
Article No.: 40510
EAN: 40 15588 40510 5
Warranty: 3 years



BS 445 CONNECT

Body analysis scale

Connect device with LCD display

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock® Online*

Product features

- High-quality stainless steel electrodes
- Convenient "step-on" activation
- Large, easy-to-read LCD display
- Switchable between kg, lb and st
- 30 memory slots for each of 8 users
- Body analysis not available for persons with medical implants (e.g. pacemakers)



Size: approx. 30×30×2.3 cm
Weight: approx. 1.6 kg
Display size: 7.3×4.8 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3 × 1.5 V (AAA)
Article No.: 40441
EAN: 40 15588 40441 2
Warranty: 3 years



BS 444 CONNECT

Body analysis scale

Our bestseller

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock® Online*

Product features

- High-quality stainless steel electrodes
- Convenient "step-on" activation
- Large easy-to-read LCD display
- Switchable between kg, lb and st
- 4 high-precision strain gauge sensors for accurate results
- 30 memory slots for each of 8 users
- Body analysis not available for persons with medical implants (e.g. pacemakers)



Size: approx. 30×30×2.3 cm
Weight: approx. 1.6 kg
Display size: 7.3×4.8 cm
Weight display: up to 180 kg, 396 lb, 28 st 5 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3 × 1.5 V (AAA)
Article No.: 40444
EAN: 40 15588 40444 3
Warranty: 3 years

Body analysis scale



BS 465 (white)

Body analysis scale

With colour indication

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- LCD display shows by colour the difference between the current weight and the target weight
- Green: under target weight
Blue: target weight
Red: above target weight

Product features

- Automatic recognition of up to 8 users
- 9 memory slots for 9 users
- Individual adjustment of the target weight possible
- 4 high-precision strain gauge sensors for accurate readings
- 6 mm tempered safety glass platform
- Athlete mode for athletes
- Convenient "step-on" activation
- Automatic switch-off after 8 seconds
- Integrated calorie requirement analysis
- Switchable between kg, lb and st
- Body analysis not available for persons with medical implants (e.g. pacemakers)



Size: approx. 30×30×3 cm
Weight: approx. 1.7 kg
Display size: 4.4×9.1 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 4 × 1.5 V (AAA)
Article No.: 40483
EAN: 40 15588 40483 2
Warranty: 3 years

*iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer.
Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.2



ecomed BS-70E

Body analysis scale

Top price-quality ratio

- Measurement of weight, body fat, body water, muscle mass and bone mass
- Determination of BMI and BMR
- Convenient "step-on" activation
- High-quality stainless steel electrode

Product features

- Tread area of 6 mm safety glass
- 10 memory spaces
- Automatic switch-off after 10 seconds
- Incl. batteries



Size: approx. 30×30×1.75 cm
Weight: approx. 1.6 kg
Display size: 8.0×3.55 cm
Weight display: up to 180 kg
Graduation: 100 g, 0.2 lb
Batteries: 1×3 V (CR2032)
Article No.: 23501
EAN: 40 15588 23501 6
Warranty: 3 years



PS 470

XL glass personal scale

Particularly resilient

- Extra large, easy-to-read display (71 mm high)
- Tread area made of extra white safety glass with rounded corners
- 4 high-precision strain gauge sensors for accurate results

Product features

- Energy-saving automatic switch-off
- Convenient "step-on" activation
- Switchable between kg, lb and st
- Overload indication
- Low battery display



Size: approx. 34×34×2.5 cm
Weight: approx. 2 kg
Display size: 9.3×7.7 cm
Weight display: up to 250 kg, 550 lb, 39 st 6 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 1×3 V (CR2032)
Article No.: 40547
EAN: 40 15588 40547 1
Warranty: 3 years



PS 450

Bamboo personal scale

From genuine bamboo

- Personal scale with high-quality bamboo surface for a comfortable feeling
- Easy to read LCD display
- Energy-saving automatic switch-off

Product features

- Convenient "step-on" activation
- Extra flat design
- Incl. batteries



Size: approx. 30×30×3.3 cm
Weight: approx. 1.2 kg
Display size: 8.3×4.1 cm
Weight display: up to 250 kg, 550 lb, 39 st 6 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40548
EAN: 40 15588 40548 8
Warranty: 3 years



PS 440

Bamboo personal scale

Magic LED display


- Decorative personal scale made of bamboo
- Invisible LED display with luminous numbers - display only during weighing
- Energy-saving automatic switch-off

Product features

- 4 high-precision strain gauge sensors for accurate results
- Convenient "step-on" activation
- Switchable between kg, lb and st
- Overload indication
- Low battery display



Size: approx. 30×30×3.2 cm
Weight: approx. 1 kg
Display size: 7×5 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40544
EAN: 40 15588 40544 0
Warranty: 3 years



PS 435

Design glass scale

Elegant design

- Invisible LED display (only appears during the weighing)
- Large easy-to-read LED display (28 mm)
- Treads of the standing area made of 6 mm safety glass with rounded corners
- 4 high-precision strain gauge sensors for accurate readings

Product features

- Convenient "step-on" switch
- Automatic switch-off after 10 seconds
- Extra slim design
- Switchable between kg, lb and st
- Overload indication
- Low battery display



Size: approx. 30×30×2.1 cm
Weight: approx. 1.5 kg
Display size: 7×5 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40434
EAN: 40 15588 40434 4
Warranty: 3 years



PS 425

Personal scale

Our talking scale

- Voice output in 4 languages: DE, GB, FR, ES
- Easy-to-read LCD display
- 4 high-precision strain gauge sensors for accurate results
- Energy-saving automatic switch-off

Product features

- Tread area made of 6 mm safety glass with rounded corners
- Convenient "step-on" activation
- Switchable between kg and lb
- Overload indication
- Low battery display



Size: approx. 31.9 × 31.9 × 2.8 cm
Weight: approx. 1.8 kg
Display size: 7.5×3.6 cm
Weight display: up to 180 kg, 396 lb, 28 st 4 lb
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 3×1.5 V (AAA)
Article No.: 40448
EAN: 40 15588 40448 1
Warranty: 3 years

Personal scale



PS 400

Glass personal scale

Elegant design

- Four high-precision strain gauge sensors for accurate results
- Tread area of 6 mm safety glass
- Extra slim design

Product features

- Convenient "step-on" switch
- Large easy-to-read LCD display
- Switchable between kg, lb and st
- Automatic switch-off after 12 seconds



Size: approx. 30×30×1.7 cm
Weight: approx. 1.5 kg
Display Size: 7.2×3.0 cm
Weight display: up to 150 kg, 330 lb, 24 st
Graduation: 100 g, 0.2 lb, 1/4 lb
Batteries: 1×3 V (CR2032)
Article No.: 40455
EAN: 40 15588 40455 9
Warranty: 3 years



ecomed PS-72E

Personal scale

Top price-quality ratio

- Convenient "step-on" switch
- Tread area of 6 mm safety glass

Product features

- Automatic switch-off after 10 seconds
- Incl. batteries



Size: approx. 30×30×2.0 cm
Weight: approx. 1.5 kg
Display Size: 6.9×3.4 cm
Weight display: up to 150 kg
Graduation: 100 g
Batteries: 1×3 V (CR2032)
Article No.: 23511
EAN: 40 15588 23511 5
Warranty: 2 years





PSD

Analogue personal scale

Our bestseller

- Trendy retro design
- Battery-free, mechanical operation
- Large analogue full-view scale

Product features

- Non-slip tread area
- Solid metal casing



Size: approx. 28.5×42.5×8.5 cm
Weight: approx. 3.6 kg
Weight display: up to 150 kg
Graduation: 500 g
Article No.: 40461
EAN: 40 15588 40461 0
Warranty: 3 years



PS 100 (red)

Analogue personal scale

Retro style

- Trendy retro design
- Large analogue full-view scale
- Non-slip tread area

Product features

- Battery-free, mechanical operation
- Solid metal casing



Size: approx. 29×43×8.5 cm
Weight: approx. 3.2 kg
Weight display: up to 150 kg
Graduation: 500 g
Article No.: 40498
EAN: 40 15588 40498 6
Warranty: 3 years



PS 100 (black)

Analogue personal scale

Retro style

- Trendy retro design
- Large analogue full-view scale
- Non-slip tread area

Product features

- Battery-free, mechanical operation
- Solid metal casing



Size: approx. 29×43×8.5 cm
Weight: approx. 3.2 kg
Weight display: up to 150 kg
Graduation: 500 g
Article No.: 40499
EAN: 40 15588 40499 3
Warranty: 3 years

Blood glucose monitor | Mask



MEDITOUCH 2 CONNECT Blood glucose monitor

Connect device

- Precise measurement accuracy according to the new ISO 15197: 2015 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock® Online*

Product features

- Pre and post-meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0.6 µl amount of blood
- Extra large, easy-to-read display
- Incl. diabetic diary
- Certified medical device



Starter set incl. 10 test strips, lancing device, 10 lancets, AST cap, control solution and bag



Size: approx. 5×9.8×1.1 cm
Weight: approx. 0.047 kg
Batteries: 2×3 V (CR 2032)
dual Article No.: 79048
EAN: 40 15588 79048 5
Warranty: 3 years



MEDITOUCH 2 Blood glucose monitor

Our classic

- Precise measurement accuracy according to the new ISO 15197: 2015 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots

Product features

- Calculates average values for 7/14/30/90 days
- Pre and post-meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0.6 µl amount of blood
- Extra large, easy-to-read display
- Certified medical device



Starter set incl. 10 test strips, lancing device, 10 lancets, AST cap, control solution and bag



mg/dL
Article No.: 79030
EAN: 40 15588 79030 0



mmol/L
Article No.: 79034
EAN: 40 15588 79034 8

Size: approx. 5×9.8×1.1 cm
Weight: approx. 0.03 kg
Batteries: 2×3 V (CR 2032)
Warranty: 3 years



RM 100 Particle filtering half mask

With additional clip for perfect support

- Protection class FFP2 NR: Protection against harmful substances such as solid and liquid dusts, smoke and aerosols
- Efficient 3-layer filter system (Spun bond, melt bond, non-woven – high level filtering of 95% of all particles in the air)
- Hygienic individual packaging prevents contamination of the mask before use

Product features

- Tested and CE-Certified according to standard EN 149:2001 + A1:2009; regulation (EU) 2016/425 for personal protective equipment
- Effective prevention of dust, bacteria, allergens and germs
- One size fits all
- Skin-friendly material
- Excellent fitting: V-shape with nose contour design
- 3D design for comfortable breathing
- Ear loop with elastic band and additional clip for perfect fit
- Mask for single use
- Single foiled PE bags



10 pack
Article No.: 33333
EAN: 40 15588 33333 0



5 pack
Article No.: 33335
EAN: 40 15588 33335 4



Colour: black
Article No.: 33334
EAN: 40 15588 33334 7

Size: 12×15.5 cm (PE-bags)
Weight: 0.01 kg

NEW



IN 155

Inhaler

Child-friendly Inhaler

- Nebulization with compressed air technology
- Reliable relief of respiratory complaints
- Efficient inhalation through high nebulization performance (approx. 0.35 ml/min.)
- Includes stickers for individual design and cards to colour in
- Motif mask for babies and children enables easy inhalation

Product features

- For the effective treatment of the upper and lower respiratory tract (e.g. asthma or colds) in children and adults
- 2-in-1: Inhaler and nasal wash
- Storage bag included
- Extensive accessories for children and adults
- Certified medical device

Incl. nebulizer bottle, 2 air tubes (approx. 180 cm), motif mask for babies and children, adult mask, angled mouth piece, nasal piece, nasal washer, 5 replacement filters, storage bag, stickers, 3 postcards



Size: approx. 18 × 14 × 9.4 cm
Weight: approx. 1.45 kg
Power supply: 230V~ 50 Hz, 1A
Filling volume: max. 8 ml
Nebulizing output: approx. 0.35 ml/min
Particle size: approx. 3 µm
Article No.: 54555
EAN: 40 15588 54555 9
Warranty: 3 years

NEW



IN 165

Inhaler

Perfect for travel

- High efficiency through innovative micro-compressor technology
- Reliable relief of respiratory complaint
- Powered by rechargeable battery: Ideal for on the go
- Includes stickers for individual design and cards to colour in
- Motif mask for babies and children enables easy inhalation

Product features

- For the targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Low-noise nebulization
- Automatic switch-off with an empty solution tank
- Storage bag included
- Extensive accessories for children and adults
- Certified medical device

Incl. nebulizer bottle, 2 air tubes (approx. 180 cm), motif mask for babies and children, adult mask, angled mouth piece, nasal piece, nasal washer, 5 replacement filters, storage bag, stickers, 3 postcards



Size: approx. 10.8 × 7.2 × 4.4 cm
Weight: approx. 0.23 kg
Power supply: Input: 100 – 240 V~ 50/60 Hz, 0.5 A; Output: DC 5 V/2 A
Batteries: Lithium-ion battery; 3.7 V, 2,570 mAh
Filling volume: max. 6 ml
Nebulizing output: 0.25 ml/min
Particle size: approx. 2.9 µm
Article No.: 54556
EAN: 40 15588 54556 6
Warranty: 3 years



IN 510

Inhaler

With comprehensive accessories

- Nebulization with compressed air technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories

Product features

- Durable and modern design
- Easy to use
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Storage bag included
- Certified medical device

Incl. nebulizer bottle, air tube (approx. 180 cm), adult mask, paediatric mask, angled mouth piece, nose piece, 5 replacement filters, storage bag



Size: approx. 18 × 14 × 9.4 cm
Weight: approx. 1.45 kg
Power supply: 230V~ 50 Hz, 1A
Filling volume: max. 8 ml
Nebulizing output: approx. 0.35 ml/min
Particle size: approx. 3 µm
Article No.: 54547
EAN: 40 15588 54547 4
Warranty: 3 years

Inhaler



IN 600

Inhaler

With rechargeable battery

- High efficiency through innovative micro compressor technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Powered by rechargeable battery and micro USB: ideal for on the go
- Comprehensive set of accessories

Product features

- 4 LED battery charge indicator
- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulization
- Compact, light-weighted and easy to use
- Automatic switch-off with an empty solution tank
- Incl. storage bag
- Certified medical device



Incl. nebulizer bottle, air tube (approx. 80 cm), adult mask, paediatric mask, angled mouth piece, nose piece, 5 replacement filters, storage bag



Size: approx. 10.8 × 7.2 × 4.4 cm
Weight: approx. 0.23 kg
Power supply: Input: 100 – 240 V~ 50/60 Hz, 0.5 A; Output: DC 5 V/2 A
Batteries: Lithium-ion battery; 3.7 V, 2,570 mAh
Filling volume: max. 6 ml
Nebulizing output: approx. 0.25 ml/min
Particle size: approx. 2.9 µm
Article No.: 54542
EAN: 40 15588 54542 9
Warranty: 3 years



IN 530

Inhaler

Perfect for travelling

- High efficiency through innovative micro compressor technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Powered by rechargeable battery: ideal for travelling
- Comprehensive set of accessories

Product features

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulization
- Compact, light-weighted and easy to use
- Automatic switch-off with an empty solution tank
- Incl. storage bag
- Certified medical device



Incl. nebulizer bottle, air tube (approx. 180 cm), adult mask, paediatric mask, angled mouth piece, nose adapter, 5 replacement filters



Size: approx. 4.5 × 4.3 × 11.5 cm
Weight: approx. 0.19 kg
Power supply: Input: AC 100-240 V/50~60 Hz
Power output: DC 5 V/2 A
approx. 2.9 µm
Particle size:
Nebulization output (average): 0.25 ml/min
Tank capacity: max. 6 ml
Charging cable: Micro USB
Article No.: 54540
EAN: 40 15588 54540 5
Warranty: 3 years



IN 525

Ultrasonic inhaler

Small inhaler - perfect for travelling

- High efficiency through innovative micro-membrane nebulization (mesh technology) with ultrafine droplets
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories

Product features

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulization
- Low power consumption
- Automatic switch-off with an empty solution tank
- Incl. storage bag
- Certified medical device

Incl. mouthpiece and masks for children and adults, storage case



Size: 7.2 × 4 × 10.3 cm
Weight: approx. 0.1 kg (without batteries)
Batteries: 2 × 1.5 V (AA)
Filling volume: max. 8 ml
Nebulizing output: ≥ 0.2 ml/min
Particle size: approx. <5 µm
Article No.: 54115
EAN: 40 15588 54115 5
Warranty: 3 years



LC 860

Electrical lice comb

Hair-gentle treatment

- Simple combing frees hair from head lice and nits
- No side effects, no use of chemical agents
- Also suitable for dogs and cats

Product features

- The lice are destroyed by an electrical current upon contact with the teeth of the comb
- With integrated cleaning brush
- Batteries included (1×1.5 V AA)



Size: approx. 19.3 × 6.6 × 2.9 cm
Weight: approx. 0.095 kg
Batteries: 1 × 1.5 V (AA)
Article No.: 41017
EAN: 40 15588 41017 8
Warranty: 3 years



LC 875

Electrical lice comb

With LED-light


- Integrated LED light for an easier detection and removal of lice
- Simple combing frees hair from head lice and nits
- No side effects, no use of chemical agents
- Also suitable for dogs and cats

Product features

- The lice are destroyed by an electrical current upon contact with the teeth of the comb
- With integrated cleaning brush
- Batteries included (1×1.5 V AA)



Size: approx. 20.9 × 6.6 × 3.1 cm
Weight: approx. 0.1 kg
Batteries: 1 × 1.5 V (AA)
Article No.: 41020
EAN: 40 15588 41020 8
Warranty: 3 years



HG 100

Hygrometer

Including magnetic mount

- Displays humidity and room temperature for a healthy climate indoors
- Saves the highest and lowest values
- For hanging on the wall (magnetic mount) or standing upright

Product features

- Temperature measurement range: -10°C to +50°C (+14 °F to +122 °F)
- Humidity measurement range: 20% - 90%
- Comfort zone indicator (smiley display)
- Time display
- Alarm function
- Incl. batteries (1.5 V AAA)



Size: approx. 8.2 × 7 × 2 cm
Weight: approx. 0.07 kg
Batteries: 1 × 1.5 V (AAA)
Article No.: 60079
EAN: 40 15588 60079 1
Warranty: 3 years



TT 200

3-in-1 Electro therapy device

TENS, EMS and Massage all in one

- TENS, EMS and massage with multiple programs
- 22 pre-installed programs (TENS, EMS and massage) with adjustable times from 5 to 90 minutes
- Large, easily readable display

Product features

- TENS therapy for 10 body parts: Neck, shoulder, arm, hand, back, abdomen, hip, leg, foot & joint
- EMS therapy for 9 body parts: Neck, Shoulder, Arm, Hand, Back, Abdomen, Hip, Leg & Foot
- 3 massage programs included for all body parts: Knead, Rub & Tap
- Backlight - ideal in bright or low light
- 40 intensity levels



west
Article No.: 88347
EAN: 40 15588 88347 7



east
Article No.: 88357
EAN: 40 15588 88357 6

Size: approx. 10.9×5.5×2.3 cm
Weight: approx. 0.076 kg (without batteries)
Batteries: 3×1.5 V (AAA)
Power supply: DC 4.5 V
Power output: max. 120 mA
Frequency: 2 – 120 Hz
Warranty: 3 years



PT 100

PainShield ultrasonic therapy

Promotes blood circulation

- Therapeutic, portable pain therapy device for soft tissue and musculoskeletal system
- Improves blood circulation and accelerates the healing process of pain - without side effects
- Easy to use anywhere and anytime – at home, at work or even while you sleep (no clinical setting necessary)

Product features

- Extraordinarily advanced low intensity, low frequency ultrasound technology for safe use
- Treatment of large, hard to reach (body) surfaces or painful areas
- No side effects
- No ultrasound gel required
- Rechargeable lithium-ion battery



Size: approx. 13.7×4.7×2.1 cm
Weight: approx. 0.8 kg
Batteries: Lithium-ion battery, 3.7 V, 1,250 mAh
Power supply: Input: 100 – 240 V~ 50/60 Hz, 500 mA; Power output: DC 5 V/2.4 A
Article No.: 88346
EAN: 40 15588 88346 0
Warranty: 3 years



UHW

Air humidifier

Comprehensive set of accessories

- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Steplessly adjustable misting levels

Product features

- Automatic switch-off when the tank is empty
- Low-noise and energy-saving
- Water tank 4.200 ml
- For rooms up to approx. 30 m²



Size: approx. 26.5×33×20 cm
Weight: approx. 1.46 kg
Power supply: 220-240 V~ 50/60 Hz
Power output: 30 watts
Nebulizing output: max. 300 ml/h
Operating time: variable depending on the setting
Article No.: 60065
EAN: 40 15588 60065 4
Warranty: 3 years



AH 660
Air humidifier

Including water filter

- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Steplessly adjustable misting levels

Product features

- Removable tank
- Automatic switch-off when the tank is empty
- Low-noise and energy-saving
- Water tank 4.500 ml
- For rooms up to approx. 30 m²
- Incl. water filter



Size: approx. 24×17×36 cm
Weight: approx. 1.3 kg
Power supply: 220 – 240V~ 50/60 Hz
Power output: 30 watts
Nebulizing output: max. 300 ml/h
Operating time: approx. 15 hours
Article No.: 60067
EAN: 40 15588 60067 8
Warranty: 3 years



AH 661
Air humidifier

Germ-reducing pre-heating

- Air humidifier with ultrasound technology to improve the breathing air quality
- Hygienically clean through an integrated pre-heating function to reduce germs and bacteria
- Integrated aroma dispenser for scented oils

Product features

- 3 mist intensities selectable
- High efficiency through micro-fine atomisation
- Automatic switch-off when the tank is empty
- Water tank 3.500 ml
- Touch control buttons
- Low-noise and energy-saving
- Timer function (2/4/6 hrs)



Size: approx. 18.5×18.5×28.1 cm
Weight: approx. 1.1 kg
Power supply: 220 – 240V~ 50/60 Hz
Power output: 95 watts
Nebulizing output: cold 250 ml/h, warm 350 ml/h
Operating time: approx. 10 – 14 hours
Article No.: 60052
EAN: 40 15588 60052 4
Warranty: 3 years



AH 662
Air humidifier

Aroma dispenser for scented oils

- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Integrated aroma dispenser for scented oils

Product features

- Automatic switch-off when the tank is empty
- Low-noise and energy-saving
- Water level indicator (full=blue, empty=red)
- Water tank 900 ml
- Incl. cleaning brush



Size: approx. 15×15×22 cm
Weight: approx. 0.45 kg
Power supply: 100 – 240V~ 50/60 Hz
Power output: 12 watts
Nebulizing output: max. 60 ml/h
Operating time: approx. 15 hours
Article No.: 60077
EAN: 40 15588 60077 7
Warranty: 3 years

Aroma diffuser



AD 620 Aroma diffuser

Our bestseller

- Aroma diffuser for scented oils
- By using ultrasound technology, fine mist of fragrance is created and is spread throughout your entire room
- Wellness light with a colour change in 5 colours
- Timer function: 2 and 4 hours

Product features

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Incl. refill cup



Size: 13.2 × 13.2 × 14.5 cm
Weight: approx. 0.25 kg
Power supply: 100-240 V~ 50/60 Hz,
DC 24 V/500 mA~ 12 W
Water tank: 300 ml
Nebulizing output: max. 30 ml/h
Operating time: approx. 10 hours
Article No.: 60082
EAN: 40 15588 60082 1
Warranty: 3 years



AD 625 Aroma diffuser

Made of genuine bamboo

- Decorative bamboo aroma diffuser for scented oils
- Creates a fine mist with ultrasonic technology, spreading the fragrance throughout the entire room
- Wellness light with a colour change in 6 colours
- Stylish living accessory without visible control elements

Product features

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving



Size: approx. ø 16.6 × 7.2 cm
Weight: approx. 0.68 kg
Power supply: Mains adaptor
Water tank: 100 ml
Nebulizing output: max. 16 ml/h
Operating time: approx. 6 hours
Article No.: 60083
EAN: 40 15588 60083 8
Warranty: 3 years



AD 635 Aroma diffuser

Made of frosted glass

- Aroma diffuser for scented oils with high-quality housing made of frosted glass
- By using ultrasound technology, fine mist of fragrance is created and is spread throughout your entire room
- Wellness light with a colour change in 6 colours
- Stylish accessory for the living room without visible switches

Product features

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving



Size: ø 15 × 17.1 cm
Weight: approx. 0.6 kg
Power supply: Mains adaptor
Water tank: 100 ml
Nebulizing output: max. 16 ml/h
Operating time: approx. 6 hours
Article No.: 60085
EAN: 4015588 60085 2
Warranty: 3 years



LT 460
Daylight lamp

As flat as a tablet

- Energising daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Uniform light distribution with LED technology
- 2 light intensities selectable

Product features

- Light intensity: up to 10.000 Lux
- Comfortable one-button operation
- Can be used horizontally and vertically
- Ideal for home and for travelling
- LED technology
- Incl. storage bag
- Certified medical device



Size: approx. 26×15×1.7 cm
Weight: approx. 0.5 kg
Light area: 20.5×12 cm
Power supply: Input: 100 – 240 V~ 50/60 Hz, 500 mA; Output: 12 V-1 A
Power output: 12 watts
Cable length: approx. 180 cm
Article No.: 45220
EAN: 40 15588 45220 8
Warranty: 3 years



LT 500
Daylight lamp

With wellness light in various colours

- Energising daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Uniform light distribution through LED technology
- Colour-changing wellness light with 4 colours

Product features

- Light intensity: up to 10.000 Lux
- Long-lasting, energy-saving LED lighting
- 2 light intensities can be selected
- Particularly silent
- Comfortable two-button operation
- Certified medical device



Size: approx. 25×25×8 cm
Weight: approx. 0.6 kg
Light area: ø 25 cm
Power supply: 220 – 240 V~ 50 Hz 0.2 A
Power output: 23 watts
Cable length: approx. 180 cm
Article No.: 45226
EAN: 40 15588 45226 0
Warranty: 3 years



IR 850
Infrared lamp

With illuminated LED button

- Intense infrared light to relax muscles or treat colds
- Three-stage electronic timer (5, 10 or 15 minutes) with LED indicator
- Cool down function with auto off

Product features

- Ceramic glass to block UV light
- Adjustable tilt angle (-15° to +35°)
- Recommended treatment distance: 40 cm
- Illuminated LED on/off button
- With hot warning indicator
- Extra long cable (200 cm)
- Incl. lamp (halogen)
- Certified medical device



Size: approx. 20×18×28 cm
Weight: approx. 1.2 kg
Power supply: 220-240 V~ 50-60 Hz
Power output: 300 watts
Article No.: 88303
EAN: 40 15588 88303 3
Warranty: 3 years

Body trainer & Activity tracker



VIFIT RUN Activity tracker

With heart rate measurement

- Tracks steps, distance, burned calories, activity time, sleep and heart rate function
- Individual daily target can be set – with reminder function
- Displays incoming calls or messages – SMS, WhatsApp, Facebook Messenger
- Dynamic heart rate measurement

Product features

- 5 heart rate training zones: Warm up, Fat burn, Aerobics, Cardio, High
- Sleep mode: tracks activity during sleep and hours slept
- Displays date, time and battery status
- Alarm clock function with vibration alarm
- Display of the current day's weather
- Stores measurements of the last 30 days
- Bluetooth® data transfer to the VitaDock+ App for iOS and Android and to VitaDock Online*



Size: approx. 4.15×2×1.18 cm
Weight: approx. 0.032 kg
Batteries: Lithium-ion battery
Article No.: 79492
EAN: 40 15588 79492 6
Warranty: 3 years



EMS Body trainer

Targeted muscle building

- Helps to build muscles
- Auto programs: massage, acupuncture, tapping, combination
- Set including 3 controllers with 1 pad for the abdomen and 2 pads for arms or legs
- Automatic switch-off after 23 minutes

Product features

- Intensity adjustable in 15 steps
- PU leather pads



Including 2 pads for arms and legs



Size: 17×18.5 cm (abdomen pad)
5×20 cm (arm pad)
Heating zone: 30 g (abdomen pad)
13 g (arm pad)
Batteries: 3×CR2032
Article No.: 79524
EAN: 40 15588 79524 4
Warranty: 3 years





MCG 820
COMFORT **Shiatsu acupressure**
massage seat cover

Soft-Touch technology

- Gel massage heads provide a soft and natural massage experience
- Soothing shiatsu massage for the back
- Soft sprung massage heads gently massage the neck area

Product features

- 3 selectable massage zones for shiatsu massage: total, upper or lower back or point-specific spot massage
- Acupressure massage in the waist and thigh area with 3 selectable intensity levels
- Separate, height-adjustable neck massage
- Additional red light and heat function
- Handy cable remote control
- Removable and washable neck cover
- Easy to clean material
- Convenient Velcro for fastening at the back of the chair



silver



black



Size: approx. 124.5 × 19 × 43.5 cm
Weight: approx. 7 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 40 watts
Article No.: 88991 (silver)
 88995 (black)
EAN: 40 15588 88991 2 (silver)
 40 15588 88995 0 (black)
Warranty: 3 years



MCG 810
COMFORT **Shiatsu massage**
seat cover

Soft-Touch technology

- Gel massage heads provide a soft and natural massage experience
- Soothing shiatsu massage for the back
- Soft sprung massage heads gently massage the neck area

Product features

- 3 selectable massage zones for shiatsu massage: total, upper or lower back or point-specific spot massage
- Separate, height-adjustable neck massage
- Separate vibration massage in the seat area with 3 selectable intensity levels
- Additional red light and heat function
- Handy cable remote control
- Removable and washable neck cover
- Easy to clean material
- Convenient Velcro for fastening at the back of the chair



Size: approx. 124.5 × 19 × 43.5 cm
Weight: approx. 5.5 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 48 watts
Article No.: 88990
EAN: 40 15588 88990 5
Warranty: 3 years

Massage seat cover



MC 825
Shiatsu acupressure massage
seat cover

Acupressure for waist and thigh area

- Targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist and thigh area
- Separate, height-adjustable neck massage


Product features

- 3 selectable massage zones for shiatsu massage: total, upper or lower back or point-specific spot massage
- 3 intensity levels of acupressure massage selectable
- Intensive shiatsu massage in the back area
- Additional red light and heat function
- Auto-off after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient Velcro for fastening at the back of the chair



Size: approx. 110 × 42.5 × 18 cm
Weight: approx. 5.9 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 40 watts
Article No.: 88939
EAN: 40 15588 88939 4
Warranty: 3 years

Massage seat cover



MC 826
Shiatsu acupressure massage
seat cover

Comprehensive relaxation with tapping and acupressure massage

- Targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist and thigh area with 3 selectable intensity levels
- 2 massage programs for the back selectable (shiatsu and tapping massage)

Product features

- Relaxing massage for the total, upper or lower back
- Separate, automatic height-adjustable neck massage with switchable rotation direction
- Separate vibration massage in the seat with 3 selectable intensity levels
- Targeted spot massage for tapping massage
- For tapping massage 3 selectable programs and 3 intensity levels
- Additional red light and heat function
- Auto-off after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient Velcro for fastening at the back of the chair



Size: approx. 120 × 26 × 46 cm
Weight: approx. 8.5 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: approx. 60 watts
Article No.: 88926
EAN: 40 15588 88926 4
Warranty: 3 years



MCN (new)
Shiatsu massage
seat cover

New Generation

- Targeted, intensive massage for tangible sense of well-being
- Separate neck massage
- Shiatsu massage for the entire back or as point-specific spot massage

Product features

- 3 selectable massage zones: total, upper or lower back
- Additional red light and heat function
- Auto-off after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient Velcro for fastening at the back of the chair



Size: approx. 118 × 40 × 9 cm
Weight: approx. 5.5 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 48 watts
Article No.: 88930
EAN: 40 15588 88930 1
Warranty: 3 years



MCH
Vibration massage
seat cover

With additional heat function

- 4 massage zones in the back and seat area
- 3 selectable massage intensities
- 5 selectable massage programs
- Additional heat function

Product features

- 5 vibration motors
- Auto-off after 15 minutes
- Handy cable remote control
- Easy to clean material
- Very handy and easy to store
- Also ideal for office use or in the car
- Incl. car power adapter



Size: approx. 107 × 47 × 3.5 cm
Weight: approx. 1.1 kg
Power supply: 100 – 240 V~ 50/60 Hz, 12 V DC 800 mA
Power output: 9.6 watts
Article No.: 88935
EAN: 40 15588 88935 6
Warranty: 3 years



MM 825
Massage mat

Massage for the entire body


- 4 selectable massage zones: upper and lower back, upper and lower leg incl. 2 intensity levels for each massage zone
- Comfortable vibration massage for the whole body
- Soothing heat function in the back area

Product features

- 2 intensity levels for each massage zone
- 5 selectable massage programs
- Additional heat function
- Auto-off after 15 minutes
- Handy cable remote control
- Extra soft and cosy fleece cover
- Cushion in the head area of the mat for lying comfortably



Size: approx. 180 × 62 × 10 cm
Weight: approx. 1.8 kg
Power supply: 100 – 240 V~ 50/60 Hz, 12 V DC 1200 mA
Power output: 14.4 watts
Article No.: 88955
EAN: 40 15588 88955 4
Warranty: 3 years



NMG 850
COMFORT Shiatsu neck massager

Soft-Touch technology

- Gel massage heads provide a soft and natural massage experience
- Relaxing Shiatsu massage for the neck and shoulder
- 3 massage speed levels selectable

Product features

- Massage intensity can be individually adjusted by pulling the holding straps
- Multi-functional usability – also applicable on shoulders, back, abdomen, upper and lower legs
- Easy operation through integrated control panel
- Auto-off after 15 minutes
- Easy to clean material
- Additional red light and heat function
- Removable and washable neck cover



Size: approx. 39.5 × 18 × 11 cm
Weight: approx. 1.69 kg
Power supply: 100 – 240 V~ 50/60 Hz, 12 V DC 2000 mA
Power output: 24 watts
Article No.: 88993
EAN: 40 15588 88993 6
Warranty: 3 years



NM 860
Neck massager

Intensive shiatsu massage

- 4 rotating massage heads
- Massage intensity individually adjustable by pulling the holding straps

Product features

- Flexible use on shoulder, back, abdomen, upper and lower leg
- 2 rotation speeds selectable
- Additional light and heat function in the upper massage heads
- Easy wireless operation through integrated control panel
- Removable and washable neck cover



Size: approx. 58 × 34 × 20 cm
Weight: approx. 1.55 kg
Power supply: 100 – 240 V~ 50/60 Hz, 12 V DC 1500 mA
Power output: 18 watts
Article No.: 88961
EAN: 40 15588 88961 4
Warranty: 3 years

Neck massager



NM 868

Vibration neck massager

Improved heating function

- 2-level vibration massage for neck and shoulders
- Ideal for home and on the go
- Additional heat function
- Removable and washable neck cover

Product features

- Functional design
- Intensity can be individually adjusted by pulling on the holding straps
- Optional mains or battery operation
- Easy wireless operation through integrated control panel



Size: approx. 38 × 33.5 × 17 cm
Weight: approx. 1.6 kg
Power supply: 140 – 220 V~ 50/60 Hz, 5.9 V DC 1200 mA
Power output: 7 watts
Article No.: 88962
EAN: 40 15588 88962 2
Warranty: 3 years



NM 885

Shiatsu neck massager

With washable neck cover

- Relaxing shiatsu massage for neck and shoulders
- Ergonomic fit
- Additional heat function
- Removable and washable neck cover
- Multi functional usability - also applicable on shoulders, back, abdomen, upper and lower legs

Product features

- 3 massage speed levels selectable
- Easy care polyester cover and pleasant foam padding
- Intensity can be individually adjusted by pulling on the holding straps
- Optional mains or battery operation
- LED battery level indicator
- Integrated control panel



Size: approx. 40 × 18.5 × 21 cm
Weight: approx. 1.6 kg
Power supply: 220 – 240 V~ 50/60 Hz, 12.6 V DC 1500 mA
Power output: 19 watts
Article No.: 88958
EAN: 40 15588 88958 5
Warranty: 3 years



NM 890

Shiatsu neck massager

Like a real kneading massage

- Intensive shiatsu massage for the neck
- Specifically adapted to the ergonomics of the neck area
- 2 massage modes selectable: Clockwise and anti-clockwise rotation
- 3 levels of rotation speed selectable

Product features

- Intensity can be individually adjusted by pulling the holding straps
- Also applicable on shoulder, back, abdomen, upper and lower leg
- Easy wireless operation through integrated control panel
- Additional heat function
- Removable and washable neck cover



Size: approx. 38 × 18.5 × 29.5 cm
Weight: approx. 2.1 kg
Power supply: 100 – 240 V~ 50/60 Hz, 12 V DC 2500 mA
Power output: 30 watts
Article No.: 88975
EAN: 40 15588 88975 2
Warranty: 3 years



NM 870
Neck massage cushion

Relaxed travelling companion

- Pleasant vibration massage for the neck and shoulders
- Ideal for travelling and at home
- Extra soft and fluffy fleece cover

Product features

- Neck massager made of comfortable moulded foam – fits perfectly to the contours of the head and neck
- With integrated “One Touch” control panel
- With integrated storage bag for space-saving storing



Size: approx. 28 × 26 × 9 cm
Weight: approx. 0.425 kg
Batteries: 2 × 1.5 V (AA)
Power output: 3.6 watts
Article No.: 88946
EAN: 40 15588 88946 2
Warranty: 3 years

NEW



CL 300
Contour Shiatsu massage cushion

With Ergonomic Flex Technology

- “Ergonomic Flex Technology” - Soft suspension of the cushion ensures individual adaptation of the massage heads to the body contour
- Soothing shiatsu massage for neck, shoulder, back and legs
- Multifunctional usability on several body parts due to flexible design

Product features

- 4 rotating massage heads, with automatically changing direction of rotation
- Additional red light and heat function
- With integrated “One Touch” control panel
- Easy to clean material
- Small and handy design
- Convenient holding straps on the back for better fixation and handling during the massage



Size: approx. 35 × 24 × 11 cm
Weight: approx. 1.9 kg
Power supply: 100 – 240 V~ 50/60 Hz, 24 V, 1.5 A
Power output: 30 watts
Article No.: 88700
EAN: 40 15588 88700 0
Warranty: 3 years

Neck cushion | Massage cushion



MCG 800
COMFORT Shiatsu massage cushion

Soft-Touch technology

- Gel massage heads provide a soft and natural massage experience
- Soothing shiatsu massage for neck, shoulder, back and legs
- Multifunctional usability

Product features

- 4 rotating gel massage heads with selectable rotating direction
- Additional red light and heat function
- Particularly flat design for flexible, convenient operation
- Easy to clean material
- Convenient Velcro for fastening at the back of the chair
- Handy cable remote control



Size: approx. 34.5 × 22 × 12.5 cm
Weight: approx. 1.8 kg
Power supply: 100 – 240 V~ 50/60 Hz, 12 V DC 1500 mA
Power output: 18 watts
Article No.: 88992
EAN: 40 15588 88992 9
Warranty: 3 years

Massage cushion

NEW



MC 250
Shiatsu massage cushion

Soothing Shiatsu massage

- Soothing Shiatsu massage for back, shoulder, legs and neck
- 4 rotating massage heads
- Simple one-button operation
- Comfortable and soft fabric

Product features

- Integrated control panel
- Convenient Velcro for fastening at the back of the chair
- Easy to clean material
- Small and handy design
- Additional red light and heat function



Size: approx. 34 × 11 × 22 cm
Weight: approx. 1.5 kg
Power supply: 12 V DC 2A
Power output: 24 watts
Article No.: 88957
EAN: 4015588 88957 8
Warranty: 3 years

NEW



MC 850
Shiatsu massage cushion

With flexible cushion

- Soothing Shiatsu massage for back, shoulder, legs and neck
- 2 massage speed levels
- Additional red light and heat function
- With flexible flap for softer massage

Product features

- 4 rotating massage heads with automatically changing direction of rotation
- Convenient Velcro for fastening at the back of the chair
- Integrated two button control panel for simple operation
- Handy and portable thanks to side handle
- Easy to clean material
- Comfortable and soft fabric



Size: approx. 34 × 10.4 × 22.4 cm
Weight: approx. 1.3 kg
Power supply: 12 V DC 2A
Power output: 24 watts
Article No.: 88947
EAN: 40 15588 88947 9
Warranty: 3 years



**FM 880**COMFORT **Shiatsu foot massager****Soft-Touch technology**

- Gel massage heads provide a soft and natural massage experience
- Deep kneading shiatsu massage of the foot and heels
- Upper foot compression massage with 3 selectable intensity levels

Product features

- Promotes blood circulation, relaxes tense muscles and invigorates tired feet
- Additional red light and heat function
- Comfortable to use with the foot
- Easy to clean material



Size: approx. 34.5 × 32 × 15 cm
Weight: approx. 2.4 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 30 watts
Article No.: 88994
EAN: 40 15588 88994 3
Warranty: 3 years

NEW**FM 900****Foot reflexology massage****Premium foot massager**

- Intensive foot reflexology massage
- 2-in-1: Compression massage for the instep and acupressure point massage for the soles of the feet
- Inspired by traditional Chinese medicine therapy (TCM), gentle acupressure massage heads promote the flow of energy

Product features

- Promotes the blood circulation and invigorates tired feet
- Additional heat function
- 2 speed settings for acupressure point massage
- 3 intensity level for compression massage
- Easy operation thanks to integrated press buttons
- Easy to clean material



Size: approx. 35.6 × 33.4 × 16 cm
Weight: approx. 3.5 kg
Power supply: 24 V 1.5 A
Power output: 30 watts
Article No.: 88394
EAN: 40 15588 88394 1
Warranty: 3 years

**FM 883****Shiatsu foot massager****Additional red light and heat function**

- Soothing shiatsu massage
- Suitable for feet and back
- 2 selectable massage intensity levels

Product features

- 6 x 3 rotating massage heads
- Two height adjustable feet to adjust the massage position
- Handy cable rewinds under the device
- Simple one-button operation



Size: approx. 40.5 × 36.5 × 12 cm
Weight: approx. 3.3 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 30 watts
Article No.: 88396
EAN: 40 15588 88396 5
Warranty: 3 years

Foot massage

NEW



FM 890 Shiatsu foot massager

Intensive foot reflexology massage

- 2 different massage programs selectable
- Shiatsu massage of the foot soles and heels
- Upper foot compression massage
- Comfortable size – suitable up to shoe size 11 (UK)

Product features

- Additional red light and heat function
- 2 intensity levels for the compression massage selectable
- Promotes the blood circulation and invigorates tired feet
- Removable and washable inner lining
- Easy to clean material



Size:	approx. 41.7 × 35.7 × 22.4 cm
Weight:	approx. 4.1 kg
Power supply:	220 – 240 V~ 50/60 Hz
Power output:	30 watts
Article No.:	88399
EAN:	40 15588 88399 9
Warranty:	3 years



NEW



HM 650
Hand massager

Deep-acting tapping massage

- Vitalizing and deep-acting thanks to powerful tapping massage
- 6 different massage attachments for individual needs
- Effectively reduces muscle tension, relieves pain and promotes circulation

Product features

- Direct and pleasant transmission of the vibration to the tissue and the muscles
- Flexible application on the whole body, such as neck, shoulders, back, legs etc.
- Continuously adjustable speed control
- Long ergonomic handle with non-slip material for comfortable self-massage
- Auto-off after 15 minutes
- Cordless use thanks to rechargeable battery



Incl. 6 attachments



Size: approx. 41.5 × 7 × 10 cm
Weight: approx. 1.1 kg
Power supply: 9 V / 1A
Power output: 40 watts / 3,300 rpm
Charging time: 135 minutes
Operating time: 135 minutes
Article No.: 88264
EAN: 40 15588 88264 7
Warranty: 3 years

NEW



HM 630
Handheld roller massager

Vibration and rolling massage

- Invigorating vibration and rolling massage for the whole body
- Vibration massage pleasantly promotes blood circulation
- Additional massage knob for deep-seated tensions and trigger point therapy
- 3 intensity levels

Product features

- Flexible application on legs, shoulders, neck, arms, feet and other body parts
- Cordless and comfortable use thanks to rechargeable battery
- Convenient "one-button" operation
- Auto-off after 5 minutes
- Long battery life of about 2.5 h
- Splash-proof
- Ergonomic design



Size: approx. 22 × 8.5 × 8 cm
Weight: approx. 290 g
Batteries: Lithium-ion battery (3.7 V 1,000 mAh)
Waterproof: IPX6
Charging : approx. 3 h
Article No.: 88267
EAN: 40 15588 88267 8
Warranty: 3 years

NEW



HM 300
Hand massager

Comfortable vibration massage

- Comfortable vibration massage with 3 acupressure points
- 2-in-1 massage pad attachment:
 - Natural loofah material for intense exfoliation
 - Soft microfiber fabric for a pleasant massage

Product features

- Flexible application on legs, shoulders, neck, arms, feet and other body parts
- Light and handy for convenient self-massage
- Convenient "one-button" operation
- Splash-proof, suitable for use on wet skin after shower or bath



Size: approx. 8 × 10.1 × 8.9 cm
Weight: approx. 0.2 kg
Batteries: 2 x 1.5 V (AAA)
Article No.: 88269
EAN: 40 15588 88269 2
Warranty: 3 years

Hand massager | Massage gun



ITM
Hand massager

Tapping, acupressure and spot massage

- Vigorous tapping massage
- Continuously adjustment of massage intensity
- Additional red light and heat function
- Various massage attachments: acupressure and spot massage

Product features

- Easy to use thanks to ergonomic shape
- Non-slip grip
- Extra-long cable



Including 2 attachments for finger massage and 2 attachments for acupressure



Size: approx. 44 × 14 × 13 cm
Weight: approx. 1.7 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 35 watts
Article No.: 88275
EAN: 40 15588 88275 3
Warranty: 3 years



NEW

MG 500
Massage gun pro

Powerful deep tissue massage

- Deep tissue massage for muscle activation and promotion of muscle regeneration
- For self-massage at different body areas, such as shoulder muscles, pectoral muscles, dorsal muscles, gluteus muscles, thigh muscles, calf muscles
- Helps soothing muscle fatigue and soreness, stimulates blood circulation and reduces muscle recovery time

Product features

- 3 powerful vibration levels for deep muscle stimulation
- Long battery life (approx. 1.5 hours in working mode)
- Auto-off after 10 minutes
- Small travel case for easy and convenient storage
- 4 massage head attachments in different shapes for multiple massage effects:

Round: For use on large muscle groups to activate muscles and promote blood circulation

Bullet: For targeting specific areas such as trigger points or smaller muscle groups

Spinal: For relaxing muscles around the spine

Pebble: For the treatment of joint areas, deep tissue, fascia and trigger point areas



Size: approx. 23.5 × 6 × 19 cm
Weight: approx. 692 g
Batteries: Lithium-ion battery (7.4 V 2,600 mAh)
Power supply: 100 – 240 V~ 50-60 Hz, DC 9 V 1A
Charging time: approx. 2-3 h
Article No.: 88985
EAN: 40 15588 88985 1
Warranty: 3 years



NEW

MG 150
Massage gun mini

Powerful deep tissue massage

- Deep tissue massage for muscle activation and promotion of muscle regeneration
- For self-massage at different body areas, such as shoulder muscles, pectoral muscles, dorsal muscles, gluteus muscles, thigh muscles, calf muscles

Product features

- 5 powerful vibration levels for deep muscle stimulation
- Ergonomic silicone handle for better grip
- USB-C charging port
- Auto-off after 10 minutes
- 6 massage head attachments in different shapes for multiple massage effects



Size: approx. 15.5 × 13.8 × 4.5 cm
Weight: approx. 372 g
Batteries: Lithium-Ion-battery (7.4 V 2,600 mAh)
Power supply: 24 watts
Charging time: approx. 2-3 h
Article No.: 88986
EAN: 40 15588 88986 8
Warranty: 3 years



HS 200

Electric cherry stone heat pad

With natural cherry stones

- Electric heat pad with cover filled with natural cherry stones
- 3-in-1: Electric heat pad, cherry stone cushion and combined application
- Cherry stones increase comfort due to better body contouring
- 4 temperature settings with heating up indicator light

Product features

- Improved heating experience thanks to long-lasting heat output
- Cherry stone cover can also be heated separately in a microwave
- APS technology - Active Protect System: visual and acoustic temperature control
- Cover made of 100% cotton
- Cherry stone cover washable with damp cloth
- Auto-off after 90 minutes



Size: approx. 30 × 40 cm
Weight: approx. 0.9 kg
Power supply: 230 V~ 50 Hz
Power output: 100 watts
Article No.: 61168
EAN: 40 15588 61168 1
Warranty: 3 years



HP 650

Heat pad XL

New XL-Size

- Heat pad in comfortable XL size (45 x 35 cm)
- Washable thanks to detachable controller
- Particularly soft, pleasant material
- 4 temperature settings

Product features

- Safety system with overheating protection
- Quick heating with 100 watts
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Oeko-Tex® Standard 100
- Auto-off after 90 minutes



Size: approx. 45 × 35 cm
Weight: approx. 0.33 kg
Power supply: 230 V~ 50 Hz
Power output: 100 watts
Article No.: 61150
EAN: 40 15588 61150 6
Warranty: 3 years



HP 605

Heat pad

With removable cover

- Removable cover of pleasantly soft upper material
- Inner heat pad washable thanks to detachable controller
- 4 temperature settings

Product features

- Particularly soft, pleasant material
- Safety system with overheating protection
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Quick heating with 100 watts
- Cover machine-washable at 30°C
- Oeko-Tex® Standard 100
- Auto-off after 90 minutes



Size: approx. 40 × 30 cm
Weight: approx. 0.35 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 100 watts
Article No.: 61148
EAN: 40 15588 61148 3
Warranty: 3 years

Neck & back heat pad



HP 405

Heat pad

The classic with 3 temperature settings

- Pleasantly soft upper material
- Washable thanks to detachable controller

Product features

- Safety system with overheating protection
- Quick heating with 100 watts
- Handy controller
- Auto-off after 90 minutes
- Oeko-Tex® Standard 100



Size: approx. 30 × 40 cm
Weight: approx. 0.35 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 100 watts
Article No.: 61003
EAN: 40 15588 61003 5
Warranty: 3 years



HP 460

Heat pad for back and neck

Super fluffy, extra soft

- Perfect fit for the back and neck
- Washable thanks to detachable controller
- 3 temperature settings

Product features

- Quick heating with 100 watts
- Safety system with overheating protection
- Pleasantly soft upper material
- Handy controller
- Auto-off after 90 minutes
- Oeko-Tex® Standard 100



Size: approx. 63 × 42 cm
Weight: approx. 0.5 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 100 watts
Article No.: 61012
EAN: 40 15588 61012 1
Warranty: 3 years





HP 622

Head pad for the shoulders and back

With 6 temperature settings

- Flexible heating wire and matrix sensor system for constant heat distribution
- Ergonomic shape for pleasant warmth at shoulder and neck

Product features

- Quick heating with 100 watts
- Handy controller
- Machine washable at 40°C
- Auto-off after 90 minutes
- Oeko-Tex® Standard 100



Size: approx. 56 × 52 cm
Weight: approx. 0.5 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 100 watts
Article No.: 61156
EAN: 40 15588 61156 8
Warranty: 3 years



HP 630

Head pad for the shoulders and back

Cosy warmth

- Especially fluffy, soft and pleasant material
- Innovative 4D heat flow technology – for an effective heat development
- Flexible heating wire and matrix sensor system for constant heat distribution
- 4 temperature settings

Product features

- APS technology – Active protection system: visual and acoustic temperature control
- Turbo heat with powerful 100 watts
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Machine washable at 40°C
- Flexible shape allows individual adjustment to any shoulder size
- Auto-off after 90 minutes
- Oeko-Tex® Standard 100



Size: approx. 55 × 65 cm
Weight: approx. 0.72 kg
Power supply: 230 V~ 50 Hz
Power output: 100 watts
Article No.: 61158
EAN: 40 15588 61158 2
Warranty: 3 years



HB 675

Cosy heating blanket XXL

XXL version

- Cosy heating blanket with fluffy soft material
- Extra big: 200 × 150 cm
- Reversible blanket in 2-tone design
- 4 temperature settings

Product features

- Safety system with overheating protection
- Handy controller
- Quick heating with 120 watts
- Washable thanks to detachable controller
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 200 × 150 cm
Weight: approx. 2.2 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 120 watts
Article No.: 60231
EAN: 40 15588 60231 3
Warranty: 3 years

Heating blanket



HB 677
3-in-1 heating blanket

3-in-1: Heating blanket, poncho and throw

- Soft and pleasant fleece material
- Even switched off, a perfect companion for cold days
- 3 temperature settings

Product features

- Reversible blanket in 2-tone design
- Easy-care and washable thanks to detachable controller
- Quick heating with 120 watts
- Safety system with overheating protection
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 160 × 130 cm
Weight: approx. 1.56 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 120 watts
Article No.: 61170
EAN: 40 15588 61170 4
Warranty: 3 years



HB 680
Knitted heating blanket

Heating blanket in knitted design

- Cosy material due to Microplush side
- Reversible blanket with 2 designs
- 3 temperature settings

Product features

- Quick heating with 120 watts
- Handy controller
- Safety system with overheating protection
- Washable thanks to detachable controller
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 160 × 120 cm
Weight: approx. 2.2 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 120 watts
Article No.: 60233
EAN: 40 15588 60233 7
Warranty: 3 years



HDW
Heating blanket

Extra soft

- Cosy heating blanket with fluffy soft upper material
- Safety system with overheating protection
- 4 temperature settings

Product features

- Handy controller
- Reversible blanket in 2-tone design
- Quick heating with 120 watts
- Washable thanks to detachable controller
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 180 × 130 cm
Weight: approx. 1.7 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 120 watts
Article No.: 60228
EAN: 40 15588 60228 3
Warranty: 3 years



HU 662

Heated underblanket

With 2 warming zones

- Safety system with overheating protection
- 2 warming zones for body and feet
- 6 temperature settings

Product features

- Pleasantly soft upper material
- Ergonomically designed controller with LED control display on both sides for right and left-handed persons
- Quick heating with 100 watts
- Extra warm foot area
- Washable thanks to detachable controller
- Suitable for mattress sizes 200 x 90 cm and 200 x 100 cm
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 150 × 80 cm
Weight: approx. 0.9 kg
Power supply: 220 – 240 V~ 50Hz
Power output: 100 watts
Article No.: 61220
EAN: 40 15588 61220 6
Warranty: 3 years

Heated underblanket

NEW



HU 666

Heated underblanket

New comfortable design

- Safety system with overheating protection
- Fluffy soft upper material
- 3 temperature settings

Product features

- Washable thanks to detachable controller
- Handy controller
- Suitable for mattress sizes 200 x 90 cm and 200 x 100 cm
- Auto-off after 180 minutes
- Oeko-Tex® standard 100



Size: approx. 150 × 80 cm
Weight: approx. 0.9 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 55 – 60 watts
Article No.: 60214
EAN: 40 15588 60214 6
Warranty: 3 years

Heated underblanket



HU 672

Maxi fleece heated underblanket

With super soft fleece material

- Extra warm foot area
- Safety system with overheating protection
- APS technology – Active protection system: visual and acoustic temperature control
- 4 temperature settings

Product features

- Approx. 2 cm thickness for cosiness
- Machine washable at 30°C
- Two-sided detachable controller
- Auto-off after 120 minutes
- Oeko-Tex® Standard 100



Size: approx. 150 × 80 cm
Weight: approx. 1.1 kg
Power supply: 230 V~ 50 Hz
Power output: 100 watts
Article No.: 60219
EAN: 40 15588 60219 1
Warranty: 3 years



HU 674

Heated underblanket

With Auto-off

- Safety system with overheating protection
- Washable thanks to detachable controller
- Pleasant soft upper material
- 4 temperature settings

Product features

- Handy controller
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 150 × 80 cm
Weight: approx. 0.82 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 100 watts
Article No.: 61232
EAN: 40 15588 61232 9
Warranty: 3 years



HU 676

King - size heated underblanket

Extra large size

- Comfy for two people
- 2 controllers for individual comfort
- Pleasantly soft upper material
- 4 temperature settings

Product features

- Safety system with overheating protection
- Handy controller
- Washable thanks to detachable controller
- Auto-off after 180 minutes
- Oeko-Tex® Standard 100



Size: approx. 150 × 160 cm
Weight: approx. 1.5 kg
Power supply: 230 V~ 50 Hz
Power output: 2 × 100 watts
Article No.: 61230
EAN: 40 15588 61230 5
Warranty: 3 years



FWS **Foot warmer**

Suitable up to shoe size 46

- Removable fluffy soft inner lining
- Large foot space (suitable up to shoe size 46)
- 6 temperature settings and 6 timer settings

Product features

- Safety system with overheating protection
- Extra-soft breathable core
- Auto-off after selected timer
- Ergonomically designed controller with LED control display
- Inner lining machine washable at 30°C
- Quick heating with 100 watts
- Oeko-Tex® Standard 100



Size: approx. 30 × 30 × 23 cm
Weight: approx 0.58 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 100 watts
Article No.: 60258
EAN: 40 15588 60258 0
Warranty: 3 years



FW 150 **Knitted foot warmer**

Attractive design with knitted surface

- For tired feet and cold winter days
- Large foot space (suitable up to shoe size 46)
- 3 temperature settings

Product features

- Quick heating for instantly soothing heat
- Removable and fluffy soft inner lining
- Inner lining washable at 30°C
- Safety system with overheating protection
- Auto-off after 90 minutes
- Oeko-Tex® Standard 100



Size: approx. 30 × 30 × 24 cm
Weight: approx 0.75 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 100 watts
Article No.: 60259
EAN: 40 15588 60259 7
Warranty: 3 years



NEW



OL 100 Heated cape short

Smart business look

- Short warming cape in smart business look
- Mobile use thanks to the powerbank
- Perfectly suited for discreet wearing in the office
- Soothing warmth for shoulders and upper back area

Product features

- 3 temperature levels and 3 timer presets
- Hidden pocket for storing the powerbank and cable
- Can be used indoors and outdoors
- Pleasant and high-quality material
- Magnetic closure at the front
- Safety system with overheating protection
- Automatic shut-off after selected timer (20, 40 or 60 minutes)



Size: approx. 45 × 55 × 0.2 cm
Weight: approx 0.3 kg
Powerbank: 7.4 V | 2,600 mAh
Power output: 19.34 watts
Article No.: 60270
EAN: 40 15588 60270 2
Warranty: 3 years

NEW



OL 150 Heated cape long

Smart business look

- Long warming cape in smart business look
- Mobile use thanks to the powerbank
- Perfectly suited for discreet wearing in the office
- Soothing warmth for shoulders, back and hands
- 2 selectable heating zones in the shoulder and back area and inside the pockets

Product features

- 3 temperature levels and 3 timer presets
- Hidden pocket for storing the powerbank and cable
- Can be used indoors and outdoors
- Pleasant and high-quality material
- Magnetic closure at the front
- Safety system with overheating protection
- Automatic shut-off after selected timer (20, 40 or 60 minutes)



Size: approx. 45 × 88 × 0.2 cm
Weight: approx. 0.7 kg
Powerbank: 7.4 V | 2,600 mAh
Power output: 19.34 watts
Article No.: 60271
EAN: 40 15588 60271 9
Warranty: 3 years

NEW



OL 200 Heating blanket

Large heating zone

- Cosy warming blanket for lap and legs
- Mobile use thanks to the powerbank
- Perfectly suited for discreet use in the office
- 2 large heating zones each 28 x 20.5 cm

Product features

- 3 temperature levels and 3 timer presets
- Hidden pocket for storing the powerbank and cable
- Can be used indoors and outdoors
- Pleasant and high-quality material
- Integrated control panel
- Safety system with overheating protection
- Automatic shut-off after selected timer (20, 40 or 60 minutes)



Size: approx. 130 × 90 cm
Weight: approx 0.8 kg
Powerbank: 7.4 V | 2,600 mAh
Power output: 19.34 watts
Article No.: 60272
EAN: 40 15588 60272 6
Warranty: 3 years



OL 300
Heated wedge cushion

Perfect for the office

- Wedge-shaped heating cushion to improve the sitting posture
- 2-in-1 use as seat cushion or vertically as lumbar support
- Actively stimulates the back muscles, provides blood circulation and prevents incorrect posture of the spine

Product features

- Additional heat function with 3 selectable temperature levels
- Integrated control panel
- Foam core for extra comfort
- Handy and portable thanks to side handle
- Removable and washable cover
- Safety system with overheating protection
- Automatic shut-off



Size: approx. 45 × 37 × 7 cm
Weight: approx. 0.73 kg
Power supply: 100 – 240 V~ 50/60 Hz
Power output: 12 watts
Article No.: 60273
EAN: 40 15588 60273 3
Warranty: 3 years

Warmth for the office




OL 350
Heated wedge cushion

Perfect for the office

- Wedge-shaped heating cushion to improve the sitting posture
- Additional detachable and height-adjustable cushion as lumbar support
- Actively stimulates the back muscles, provides blood circulation and prevents incorrect posture of the spine


Product features

- Additional heat function with 3 selectable temperature levels
- Integrated control panel
- Memory foam core for extra comfort
- Handy and portable thanks to side handle
- Removable and washable cover
- Safety system with overheating protection
- Automatic shut-off



Size: approx. 36 × 44 × 29 cm
Weight: approx. 1 kg
Power supply: 100 – 240 V~ 50/60 Hz
Power output: 24 watts
Article No.: 60274
EAN: 40 15588 60274 0
Warranty: 3 years

Outdoor heat pad



NEW
Available
from autumn
2022

OL 700

Outdoor heated seat cushion

Heat pad to go


- Wireless use thanks to powerbank - perfectly suitable for indoor and outdoor
- Modern design with soft felt surface
- Flat heated seat cushion in compact size (40 x 40 cm)

Product features

- 3 selectable temperature settings
- Convenient "OneTouch" controller
- Controller also usable as carrying handle
- Automatic shutdown after 20 min.
- Easy care, washable pillowcase with zip
- Integrated USB charging port for powerbank (Powerbank included)
- Can be used in combination with outdoor heat cushion
- OL 750 thanks to modular push-button system



Size: approx. 40 × 40 cm
Weight: approx. 0.5 kg
Power supply: Powerbank
Power output: 10 watts
Article No.: 60277
EAN: 40 15588 60277 1
Warranty: 3 years



NEW
Available
from autumn
2022

OL 750

Outdoor heat cushion

Heating pad to go with large heating surface

- Wireless use thanks to powerbank - perfectly suitable for indoor and outdoor
- Modern design with soft felt surface
- Flat heat pad in compact size (40 x 50 cm)

Product features

- 3 selectable temperature settings
- Convenient "OneTouch" controller
- Controller also usable as carrying handle
- Automatic shutdown after 20 min.
- Easy care, washable pillowcase with zip
- Integrated USB charging port for powerbank (Powerbank included)
- Can be used in combination with outdoor heated seat cushion
- OL 700 thanks to modular push-button system



Size: approx. 40 × 50 cm
Weight: approx. 0.6 kg
Power supply: Powerbank
Output: 10 watts
Article No.: 60278
EAN: 40 15588 60278 8
Warranty: 3 years



**FS 350****Foot spa****With reflexology massage**

- Intensive foot reflexology massage
- Function for keeping the water warm

Product features

- Relaxing spa and vibration massage
- Including removable massage rollers in the middle of the foot spa to massage the soles of the feet
- Comfortable size - suitable up to shoe size 12 (UK)



Size: approx. 41 × 34 × 17 cm
Weight: approx. 1.2 kg
Power supply: 230 V~ 50 Hz
Power output: 60 watts
Article No.: 88364
EAN: 40 15588 88364 4
Warranty: 3 years

**FS 881****Foot spa****Our compact bestseller**

- Relaxing spa and vibration massage
- Soothing foot spa
- 3-in-1: Water heating, bubble and vibrating massage

Product features

- Non-slip rubber feet
- Automatic heat function
- Spout-opening
- Comfortable size - suitable up to shoe size 9.5 (UK)



Size: approx. 39 × 33 × 13 cm
Weight: approx. 1.22 kg
Power supply: 220 - 240 V~ 50 Hz
Power output: 60 watts
Article No.: 88380
EAN: 40 15588 88380 4
Warranty: 3 years

**FS 883****Foot spa****With pleasant red light**

- Pleasant red light area
- 3-in-1: Water heating, bubble and vibrating massage

Product features

- Massage rollers for foot reflexology massage
- Massage-aiding footbed
- Dry massage function
- Non-slip rubber feet
- Removable splash protection
- Comfortable size - suitable up to shoe size 11 (UK)



Including 3 pedicure attachments



Size: approx. 41 × 35 × 18 cm
Weight: approx. 1.5 kg
Power supply: 220 - 240 V~ 50/60 Hz
Power output: 80 watts
Article No.: 88387
EAN: 40 15588 88387 3
Warranty: 3 years

Foot spa



FS 886
Foldable foot spa

With pleasant red light

- Foldable foot spa for convenient and space-saving storage
- Relaxing bubble and vibration massage
- Fast heating of the water up to 45°C
- With safety bars to prevent folding during use

Product features

- Pleasant red light area
- Intensive foot reflexology massage
- Massage rollers in the middle of the foot spa to massage the bottom side of the foot
- 3 exchangeable pedicure attachments (callus remover, brush, massage attachment)
- Spout-opening
- Comfortable size - suitable up to shoe size 11 (UK)

blue



grey



Size: approx. 42.3 × 38.5 × 20.1 cm
Weight: approx. 2.1 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 500 watts
Article No.: 88392 (blue)
 88381 (grey)
EAN: 40 15588 88392 7 (blue)
 40 15588 88381 1 (grey)
Warranty: 3 years



FS 888
Foot Spa

Wellness-light in 4 colours

- Powerful bubble foot massage
- Promotes blood circulation and invigorates tired feet
- Fast heating of the water between 35 to 45°C (adjustable in 2°C steps)
- Massage rollers in the middle of the foot spa to massage the bottom side of the foot
- Wellness light in 4 colours: red, green, yellow, blue

Product features

- Integrated aroma dispenser
- Integrated display
- Optional timer: 10 to 60 minutes
- Double-walled housing to keep the water warm
- Non-slip rubber feet
- Spout-opening
- Water level mark
- Comfortable size - suitable up to shoe size 11 (UK)



Size: approx. 36.6 × 26 × 45.9 cm
Weight: approx. 3.0 kg
Power supply: 220 – 240 V
Power output: 430 Watt
Article No.: 88393
EAN: 40 15588 88393 4
Warranty: 3 years





BBS Bath spa

Extra quiet

- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels and combination program
- Heat function selectable

Product features

- Large Soft-Mat (approx. 120 x 36 cm) for lying comfort
- Integrated aroma dispenser for scented oils
- Automatic switch-off
- Timer for programming the operating time
- Extra quiet due to noise-optimised construction
- Suitable for any bath tub
- Suitable for any bath tub



Size: Soft-Mat: approx. 120 x 36 cm
Air-Tube: approx. 2.4 m
Motor: approx. 31 x 21 x 31 cm
Weight: approx. 4.3 kg (total)
Power supply: 220 – 240 V~ 50 Hz
Remote control: 3 V CR2032
Power output: 570 watts
Article No.: 88389
EAN: 40 15588 88389 7
Warranty: 3 years



MBH Bath spa

Relaxes and loosens the muscles

- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels and combination program
- Auto-off after 20 minutes

Product features

- Intensity levels adjustable on the device or via remote control
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated aroma dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma pads
- Incl. remote control



Size: Soft-Mat: approx. 120 x 36 cm
Air-Tube: approx. 2.4 m
Motor: approx. 25 x 16 x 17 cm
Weight: approx. 2.7 kg (total)
Power supply: 220 – 240 V~ 50 Hz
Remote control: 3 V CR2032
Power output: 450 watts
Article No.: 88379
EAN: 40 15588 88379 8
Warranty: 3 years



RS 650 Lounge Chair + Ottomane

Stylish interior Massage Chair

- High seating comfort
- Integrated USB charging socket
- 6 massage types in an automatic sequence program
- 3 intensity levels selectable

Product features

- Kneading, knocking, shiatsu, rolling, clapping and Swedish massage
- Massage zone settings: full, upper or lower back
- Automatic switch-off after 15 minutes
- Appealing, modern design



Size: approx. 105 x 60 x 90 cm
Weight: approx. 28 kg
Power supply: 220 – 240 V~ 50/60 Hz
Power output: 100 watts
Article No.: 88414
EAN: 40 15588 88414 6
Ottomane Article No.: 88415
Ottomane EAN: 40 15588 88415 3
Warranty: 3 years

Massage chair



RS 800
Relax massage chair

Particularly soft massage

- Particularly gentle massage thanks to 2D Massage Quad Roller mechanism
- 5 different massage types: kneading, knocking, shiatsu, clapping and Swedish massage
- 3 automatic program sequences

Product features

- Stepless adjustable backrest and footrest
- 180° chair rotatability
- Relaxing vibration massage in the seat and foot area
- Maximum massage time: 30 minutes
- Integrated overheating protection and surge detector
- Integrated body scanner



Size: approx. 105 × 75 × 110 cm
Colour: beige
Weight: approx. 55 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 120 watts
Cable length: 180 cm
Weight capacity: approx. 120 kg
Article No.: 87703
EAN: 40 15588 87703 2
Warranty: 3 years



RS 810
Relax massage chair

Particularly soft massage

- Particularly gentle massage thanks to 2D Massage Quad Roller mechanism
- 5 different massage types: kneading, knocking, shiatsu, clapping and Swedish massage
- 3 automatic program sequences

Product features

- Stepless adjustable backrest and footrest
- 180° chair rotatability
- Relaxing vibration massage in the seat and foot area
- Maximum massage time: 30 minutes
- Integrated overheating protection and surge detector
- Integrated body scanner



Size: approx. 105 × 75 × 110 cm
Colour: brown
Weight: approx. 55 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 120 watts
Cable length: 180 cm
Weight capacity: approx. 120 kg
Article No.: 87704
EAN: 40 15588 87704 9
Warranty: 3 years



RS 820
Relax massage chair

Particularly soft massage

- Particularly gentle massage thanks to 2D Massage Quad Roller mechanism
- 5 different massage types: kneading, knocking, shiatsu, clapping and Swedish massage
- 3 automatic program sequences

Product features

- Stepless adjustable backrest and footrest
- 180° chair rotatability
- Relaxing vibration massage in the seat and foot area
- Maximum massage time: 30 minutes
- Integrated overheating protection and surge detector
- Integrated body scanner



Size: approx. 105 × 75 × 110 cm
Colour: black
Weight: approx. 55 kg
Power supply: 220 – 240 V~ 50 Hz
Power output: 120 watts
Cable length: 180 cm
Weight capacity: approx. 120 kg
Article No.: 87705
EAN: 40 15588 87705 6
Warranty: 3 years



MP 810

Manicure and pedicure device

Lightweight and powerful

- Care and treatment of nails, cuticles and calluses
- Easy to use and easy replacement of the attachments
- Remove cuticles, file, shape, smooth and polish nails

Product features

- 5 attachments
- Compact device, requires little space
- Small and practical storage bag



Including 5 attachments:
Felt cone, trapezoidal grinder, cylindrical grinder, grinding cone, flame burr



Size: approx. 16.5×3×3 cm
Weight: approx. 0.8 kg
Batteries: 2×1.5 V (AAA)
Article No.: 85140
EAN: 40 15588 85140 7
Warranty: 3 years



MP 815

Manicure and pedicure device

Providing two speed levels

- Care and treatment of nails, cuticles and slight calluses
- 2 speed levels for the individual requirements of care
- High quality and durable attachments

Product features

- Clockwise and counter-clockwise rotation
- Simple to use, easy replacement of the attachments
- Incl. storage bag



Including 7 attachments: Sapphire disc (fine), Sapphire disc (coarse), sapphire callus remover, sapphire cone, felt cone, cylindrical burr and flame burr



Size: approx. 15×3.8×3.7 cm
Weight: approx. 0.1 kg
Power supply: 100 – 240 V~ 50/60 Hz
30 V DC 250 mA
Rotation speed: 2 levels: 3,800/4,600
Cable length: 150 cm
Article No.: 85153
EAN: 40 15588 85153 7
Warranty: 3 years



MP 840

Manicure and pedicure device

Suitable for diabetic patients

- Care and treatment of nails, cuticles and slight calluses
- With LED light for accurate application
- 3 speed levels for the individual requirements of care
- Stops automatically when the pressure is too high – also suitable for diabetics

Product features

- Clockwise and counter-clockwise rotation
- High quality and durable attachments
- Additional nail polishing attachment with 2 smoothing and polishing rollers
- Natural, durable nail gloss (without additional care products)
- Simple to use, easy replacement of the attachments
- Incl. storage bag



Including 7 attachments: Sapphire cone (coarse), sapphire disc (fine), cylindrical burr, flame burr, felt cone, needle burr as well as additional attachment for nail polishing



Size: approx. 15×3.8×3.7 cm
Weight: approx. 0.1 kg
Power supply: 100 – 240 V~ 50/60 Hz
18V DC 500 mA
Rotation speed: 3 levels: 5.000/6.000/7.000
Power output: 9 watts
Cable length: 150 cm
Article No.: 85155
EAN: 40 15588 85155 1
Warranty: 3 years



CR 870
Electric callus remover

Ergonomic handling

- Care and treatment of cuticles and callus
- Gentle and effective removal of excess calluses
- For smoother hands, feet and elbows

Product features

- 3 removable attachments: extra coarse, coarse and fine abrasive roll
- 2 speed levels
- Ergonomic, non-slip handle
- Incl. cleaning brush and batteries



Including 3 removable attachments



Size: approx. 15.5 × 3.43 × 8.04 cm
Weight: approx. 0.148 kg (incl. batt.)
Batteries: 3 × 1.5 V (AAA)
Rotation speed: 2 levels: 800/1,600 rpm
Article No.: 88576
EAN: 40 15588 88576 1
Warranty: 3 years

NEW



DS 400
Facial sauna

With temperature setting for aroma therapy

- Purifying and relaxing effect through moist heat
- Aroma compartment for pleasant aroma therapy
- Temperature setting from 80-95°C for aroma therapy and 100°C for steam

Product features

- For cosmetic use
- Opens and cleanses the pores, promotes blood circulation and relaxes the facial features
- Diffuser for even steam distribution
- Additional nasal attachment for precise application
- LED display for easy reading of applied settings and safe usage
- Variable timer setting from 15 to 60 mins
- Incl. measuring cup for easy filling
- Safety system with overheating protection
- Automatic switch-off after selected timer



Size: approx. 26.2 × 19.9 × 19.6 cm
Weight: approx. 0.645 kg
Power supply: 220-240 V~ 50/60 Hz
Article No.: 88242
EAN: 40 15588 88242 5
Warranty: 3 years



FB 885
Facial cleansing brush

Cleansing with 4 speed levels

- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products
- Rechargeable lithium-ion battery

Product features

- Cleanses the skin in 60 seconds (Timer)
- 4 speed levels
- Automatic switch-off after 1 minute
- Splash-proof
- Incl. protective cover, storage bag and charge base



Including 4 cleansing attachments (sensitive skin, regular skin, peeling, cosmetics sponge)



Size: approx. 18.6 × 5.0 × 2.7 cm
Weight: approx. 0.27 kg
Batteries: Lithium-ion battery
Power output: 500 mAh
Cable length: 1.30 m
Article No.: 88565
EAN: 40 15588 88565 5
Warranty: 3 years



AC 850
Cellulite massager

Intense rolling massage

- Effective massage for stimulating the connective tissue and visibly firmer skin
- Improves blood circulation
- 2 selectable massage intensities
- Light and handy design for convenient self-massage

Product features

- 3×2 rotating massage rollers (removable)
- Individually adjustable hand strap
- Especially powerful thanks to mains operation
- Easy cleaning of the removable massage rollers



Size: approx. 13×11×8.5 cm
Weight: approx. 0.35 kg
Power supply: 110 – 240 V~ 50/60 Hz
 12V DC 1200 mA
Power output: 14.4 watts
Cable length: 200 cm
Article No.: 88540
EAN: 40 15588 88540 2
Warranty: 3 years

NEW



AC 900
Cellulite massager

Multidirectional massage attachment

- Supports the improvement of cellulite appearance
- Promotes blood circulation and effectively massages different skin layers of the connective tissue
- Cordless use thanks to rechargeable battery

Product features

- Invigorating and stimulating massage for the whole body
- Flexible application on legs, buttocks, abdomen, arms and other body parts
- Multidirectional massage roller attachment for a pleasant and effective massage
- 2 intensity levels and 2 rotation directions for individual adjustment



Size: approx. 8.5×14×13 cm
Weight: approx. 0.345 kg
Battery: Lithium-ion battery
 7.4 V 2,200 mAh
 Level 1: 230 rpm±10
 Level 2: 270 rpm±10
 5 hours (Level 1)
Runtime: 5 hours (Level 1)
Charging time: 3.5 hours
Article No.: 88542
EAN: 40 15588 88542 6
Warranty: 3 years

NEW



AC 950
Cellulite massager

Ergonomic design

- Supports the improvement of cellulite appearance
- Promotes blood circulation and effectively massages different skin layers of the connective tissue
- 3 different attachments (multidirectional massage roller, bristle-silicone brush, silicone brush)
- Cordless use thanks to rechargeable battery

Product features

- Invigorating and stimulating massage for the whole body
- Flexible application on legs, buttocks, abdomen, arms and other body parts
- 2 intensity levels and 2 rotation directions for individual adjustment
- Ergonomic handle for body parts difficult to reach
- Detachable and washable attachment
- Convenient loop to hang the product up
- Splash-proof – usable after showering
- Elegant storage bag



Size: approx. 8.5×14×23 cm
Weight: approx. 0.38 kg
Battery: Lithium-ion battery
 7.4 V 2,200 mAh
 Level 1: 230 rpm±10
 Level 2: 270 rpm±10
 5 hours (Level 1)
Runtime: 5 hours (Level 1)
Charging time: 3.5 hours
Article No.: 88544
EAN: 40 15588 88544 0
Warranty: 3 years

Cosmetic mirror



CM 835
Cosmetic mirror

Compact and stable

- Illuminated cosmetic mirror
- 2 swivelling mirrors: normal and 5× magnification
- Uniform ring illumination with 12 white LED-lights
- Diameter of mirror approx. 12 cm

Product features

- High-quality chrome finish
- Stable base
- Incl. 3 batteries (1.5 V AAA)



Size: approx. 16.9 × 8 × 19 cm
Weight: approx. 0.375 kg
Diameter of mirror: approx. ø 12 cm
Batteries: 3 × 1.5 V (AAA)
Article No.: 88554
EAN: 40 15588 88554 9
Warranty: 3 years



CM 840
2-in-1 Cosmetic mirror

Our bestseller

- Illuminated cosmetic mirror
- 2 swivelling mirrors: normal and 5× magnification
- Uniform ring illumination with 18 white LED lights
- Diameter of mirror approx. 13 cm

Product features

- High-quality chrome finish
- Stable base
- Incl. 4 batteries (1.5 V AAA)



Size: approx. 20.5 × 13.5 × 31.5 cm
Weight: approx. 0.95 kg
Diameter of mirror: approx. ø 13 cm
Batteries: 4 × 1.5 V (AAA)
Article No.: 88550
EAN: 40 15588 88550 1
Warranty: 3 years



CM 848
Cosmetic mirror

Impresses 10x magnification

- Illuminated cosmetic mirror with super bright LED's
- Uniform ring illumination with 33 white LED-lights
- 2 mirrors: Normal and 10× magnification
- With big storage tray for make-up and jewellery
- Micro-USB or battery-powered

Product features

- Incl. removable magnetic magnifying mirror – can be easily mounted either on the back or the front of the mirror
- 360 degrees rotation of the mirror surface
- Touch sensor for dimming the light
- Memory light effect: The light switches on with the same light intensity with which it was switched off
- Auto-stop after 20 minutes



Size : Mirror (big) 23.6 × 16.6 cm
Size : Mirror (small) ø 9.5 cm
Size: Tray: 22.4 × 16.2 cm
Weight: approx. 0.6 kg
Power supply: 3 × 1.5 V (AAA) Micro-USB cable
Article No.: 88556
EAN: 40 15588 88556 3
Warranty: 3 years



CM 850

LED suction mirror

Extra large and bright

- Strong adhesion on all smooth surfaces
- Bright LED illumination
- Diameter of mirror approx. 19 cm
- 5× magnification

Product features

- Suction enhancement thanks to special coating
- Detachable and swivelling
- Incl. 3 batteries, 1.5 V AAA
- Telescope bracket



Size: approx. 20×10×20 cm
Weight: approx. 0.3 kg
Diameter of mirror: approx. ø 19 cm
Batteries: 3×1.5 V (AAA)
Article No.: 88558
EAN: 40 15588 88558 7
Warranty: 3 years





UPPER ARM CUFFS for blood pressure monitors

MTP Pro

Large cuff: (32 – 42 cm)
Article No.: 51037
EAN: 40 15588 51037 3

Normal cuff: (22 – 32 cm)
Article No.: 51038 (included in package)
EAN: 40 15588 51038 0

Small cuff: (17 – 22 cm)
Article No.: 51034
EAN: 40 15588 51034 2

BU 535 | BU 535 Voice

Large cuff: (30 – 42 cm)
Article No.: 51126
EAN: 40 15588 51126 4

Normal cuff: (22 – 30 cm)
Article No.: 51135 (included in package)
EAN: 40 15588 51135 6

BU 540 connect | BU 565 | BU 570 connect

Large cuff: (22 – 42 cm)
Article No.: 51299
EAN: 40 15588 51299 5

BU 510 | BU 512 | BU 514 | BU 516

Large cuff: (32 – 42 cm)
Article No.: 51169
EAN: 40 15588 51169 1

Normal cuff: (22 – 36 cm)
Article No.: 51168 (included in package)
EAN: 40 15588 51168 4



MAINS ADAPTER for blood pressure monitors

Mains adapter for MTP Pro

Article No.: 51036
EAN: 40 15588 51036 6

Mains adapter for BU 510

Article No.: 51125
EAN: 40 15588 51125 7



ACCESSORIES for blood glucose meters

MediTouch 2 | MediTouch 2 connect

Control solution
Article No.: 79039
EAN: 40 15588 79039 3

Lancets
Article No.: 79028
EAN: 40 15588 79028 7

Test strips
Article No.: 79038
EAN: 40 15588 79038 6

Lancing device

Article No.: 79002
EAN: 40 15588 79002 7



ACCESSORIES for bath spas

BBS | MBH

Air hose
Article No.: 30666
EAN: 40 15588 30666 2

Soft mat
Article No.: 88383
EAN: 40 15588 88383 5

Fragrance dispenser with 3 pads
Article No.: 30669
EAN: 40 15588 30669 3

BBS

Remote control
Article No.: 88384
EAN: 40 15588 88384 2

MBH

Remote control
Article No.: 88376
EAN: 40 15588 88376 7



WRIST BANDS

for activity tracker

ViFit connect

Black

Article No.: 79417

EAN: 40 15588 79417 9

Blue

Article No.: 79470

EAN: 40 15588 79470 4

Orange

Article No.: 79418

EAN: 40 15588 79418 6

Green

Article No.: 79471

EAN: 40 15588 79471 1

Red

Article No.: 79472

EAN: 40 15588 79472 8

Pink

Article No.: 79473

EAN: 40 15588 79473 5



SPARE PART SETS

for EMS Pad

Body Trainer

Spare part set (2 pads for the abdomen)

Article No.: 79525

EAN: 40 15588 79525 1

Spare part set (4 pads for arms / legs)

Article No.: 79526

EAN: 40 15588 79526 8



ELECTRODES

for pain therapy devices

TDD | TDP

4 Multiple electrodes (approx. 3×5 cm)

Article No.: 87080

EAN: 40 15588 87080 4

2 Multiple electrodes (approx. 5×9 cm)

Article No.: 87100

EAN: 40 15588 87100 9

TT 200 | TT 205

Replacement electrodes (8 pcs.)

Article No.: 88349

EAN: 40 15588 88349 1

Cable set (2 pcs.)

Article No.: 88359

EAN: 40 15588 88359 0

PT 100

Adhesive pads

Article No.: 88356

EAN: 40 15588 88356 9



REPLACEMENT LAMP

for infrared lamps

IRL

150 watts

Article No.: 88255

EAN: 40 15588 88255 5

IR 100 | IRH

100 watts

Article No.: 30101

EAN: 40 15588 30101 8



ACCESSORIES

for inhaler

IN 510 | IN 520

Nebulizer bottle

Article No.: 54570

EAN: 40 15588 54570 2

Angled mouth piece

Article No.: 54571

EAN: 40 15588 54571 9

Nose adapter

Article No.: 54572

EAN: 40 15588 54572 6

Adult mask

Article No.: 54573

EAN: 40 15588 54573 3

Paediatric mask

Article No.: 54574

EAN: 40 15588 54574 0

Infant mask

Article No.: 54575

EAN: 40 15588 54575 7

Air hose

Article No.: 54576

EAN: 40 15588 54576 4

Replacement filter (5 pcs.)

Article No.: 54577

EAN: 40 15588 54577 1

Nasal washer

Article No.: 54578

EAN: 40 15588 54578 8

Accessories set*

Article No.: 54579

EAN: 40 15588 54579 5

*Incl.: Nebulizer bottle, mouth piece, air hose

IN 525

Nebulizer unit

Article No.: 54116

EAN: 40 15588 54116 2

Mask set (adult and paediatric mask)

with mouth piece

Article No.: 54117

EAN: 40 15588 54117 9

Adaptor EU

Article No.: 54118

EAN: 40 15588 54118 6

Adaptor GB

Article No.: 54119

EAN: 40 15588 54119 3



AROMA ESSENCES

for air humidifiers/Aroma diffusers

Aroma Eucalyptus 10 ml

Article No.: 60031

EAN: 40 15588 60031 9

Aroma Lavender 10 ml

Article No.: 60032

EAN: 40 15588 60032 6

Aroma Ylang-Ylang 10 ml

Article No.: 60035

EAN: 40 15588 60035 7

Aroma Orange 10 ml

Article No.: 60037

EAN: 40 15588 60037 1

Aroma Pinewood 10 ml

Article No.: 60038

EAN: 40 15588 60038 8



ACCESSORIES

for inhaler

IN 600 | IN 605 | IN 530 | IN 535

Nebulizer bottle

Article No.: 54560

EAN: 40 15588 54560 3

Angled mouth piece

Article No.: 54561

EAN: 40 15588 54561 0

Nose adapter

Article No.: 54562

EAN: 40 15588 54562 7

Adult mask

Article No.: 54563

EAN: 40 15588 54563 4

Paediatric mask

Article No.: 54564

EAN: 40 15588 54564 1

Infant mask

Article No.: 54565

EAN: 40 15588 54565 8

Air hose

Article No.: 54566

EAN: 40 15588 54566 5

Replacement filter (5 pcs.)

Article No.: 54567

EAN: 40 15588 54567 2

Accessories set*

Article No.: 54568

EAN: 40 15588 54568 9

Baby motif mask (only for IN 165)

Article No.: 54557

EAN: 4015588 54557 3

Infant motif mask (only for IN 165)

Article No.: 54558

EAN: 4015588 54558 0

*Incl.: Nebulizer bottle, mouth piece, nose adapter, adult mask, paediatric mask, air hose, 5 replacement filters



ACCESSORIES

for inhaler

IN 500 | IN 550

Connecting piece

Article No.: 54524

EAN: 40 15588 54524 5

Air filter and nebulizer kit

Article No.: 54525

EAN: 40 15588 54525 2

Mouthpiece and nose adapter

Article No.: 54526

EAN: 40 15588 54526 9

Air hose

Article No.: 54527

EAN: 40 15588 54527 6

Paediatric mask

Article No.: 54528

EAN: 40 15588 54528 3

Adult mask

Article No.: 54529

EAN: 40 15588 54529 0

Nebulizer attachments (3 pcs.)

Article No.: 54531

EAN: 40 15588 54531 3

Baby motif mask (only for IN 155)

Article No.: 54558

EAN: 4015588 54557 3

Baby motif mask (only for IN 155)

Article No.: 54557

EAN: 4015588 54558 0



REPLACEMENT FILTER for air purifier

MEDISANA AIR®

1pc.
Article No.: 60390
EAN: 40 15588 60390 7



REPLACEMENT FILTER for air humidifier

AH 660 | AH 665

2 pcs.
Article No.: 60069 | 60059
EAN: 40 15588 60069 2

AH 670

2 pc.
Article No.: 60057
EAN: 40 15588 60057 9



REPLACEMENT ROLLS for manicure and pedicure devices

CR 860

2 pcs.
Article No.: 88572
EAN: 40 15588 88572 3

CR 870

Replacement roll coarse
Article No.: 88577
EAN: 40 15588 88577 8

Replacement roll extra coarse
Article No.: 88578
EAN: 40 15588 88578 5

NP 860 | MP 840

4 pcs. (2 smoothing and 2 polishing rolls)
Article No.: 85164
EAN: 40 15588 85164 3



REPLACEMENT BRUSHES for facial cleansing brushes

FB 880

3 pcs.
Article No.: 88561
EAN: 40 15588 88561 7

FB 885

4 pcs.
Article No.: 88566
EAN: 40 15588 88566 2



ACCESSORIES for manicure and pedicure devices

MANILUX L | MANILUX S | MP 820

Set of attachments (10 pcs.)*
Article No.: 85450
EAN: 40 15588 85450 7

*Incl.: Sapphire cone (coarse), sapphire disc (coarse), flame burr, sapphire ball burr, felt cone, sapphire cone (long), sapphire disc (fine), needle burr (long), cylindrical burr, drainage grinder

MEDISTYLE S | MEDISTYLE L | MPS | MP 840

Felt cone
Article No.: 85001
EAN: 40 15588 85001 1

Sapphire cone (fine)
Article No.: 85002
EAN: 40 15588 85002 8

Sapphire cone (coarse)
Article No.: 85003
EAN: 40 15588 85003 5

Sapphire disc (fine)
Article No.: 85004
EAN: 40 15588 85004 2

Sapphire ball burr
Article No.: 85005
EAN: 40 15588 85005 9

Needle burr (long)
Article No.: 85006
EAN: 40 15588 85006 6

Sapphire cone (long)
Article No.: 85007
EAN: 40 15588 85007 3

Cuticle pusher
Article No.: 85008
EAN: 40 15588 85008 0

Cylindrical burr
Article No.: 85009
EAN: 40 15588 85009 7

Sapphire disc (coarse)
Article No.: 85010
EAN: 40 15588 85010 3

Flame burr
Article No.: 85011
EAN: 40 15588 85011 0

Attachment adapter
Article No.: 85015
EAN: 40 15588 85015 8



SUBSIDIARIES

GERMANY

medisana Space Technologies GmbH

Carl-Schurz-Str. 2
D-41460 Neuss
Phone: +49 2131 36 68 0
Email: info@medisanaspace.com
Internet: www.medisana.de

AUSTRIA

medisana Austria

Carl-Schurz-Str. 2
D-41460 Neuss
Phone: +43 664 459 54 87
Email: office@medisana.at
Internet: www.medisana.at
Contact: Christian Eichholzer

medisana Service Austria

Gewerbegebiet 4
Flöcking 235/2
8200 Ludersdorf-Wilfersdorf
Phone: +43 3136 52076
Email: service.at@medisana.com

FRANCE

medisana GmbH

Carl-Schurz-Str. 2
D-41460 Neuss
Phone: +49 2131/3 66 80
Email: medisana.export@medisana.de
Internet: www.medisana.fr
Contact: Daniel Wolter

HEAD OFFICE

medisana GmbH

Carl-Schurz-Str. 2, 41460 Neuss,
Germany
Phone: +49 2131 36 68-0
Fax: +49 2131 36 68-50 95
Email: info@medisana.com
Internet: www.medisana.com

Our head office is open from Monday to Thursday
8:00 am – 4:00 pm, Friday 8:00 am – 2:00 pm.

NETHERLANDS/BELGIUM/LUXEMBOURG

medisana Benelux NV

Euregiopark 30
NL-6467 JE Kerkrade
Phone: +31 45 547 08 60
Fax: +31 45 547 08 79
Email: info@medisana.nl
Internet: www.medisana.nl
Contact: Maurice Jongen

Medisana Service Niederlande

Teknihall Benelux
p/a Vreysen
Hazeltonk 6027
4836 LA BREDA
Phone: +31 45 547 08 60
Email: info@medisana.nl
Internet: www.medisana.nl

medisana Service Belgien/Luxemburg

Teknihall bvba
Brusselstraat 33
2321 BE MEER
Phone: +31 45 547 08 60
Email: info@medisana.be
Internet: www.medisana.be

ASIA

medisana Far East Ltd.

Unit E, F, 31/F, YHC Tower
1 Sheung Yuet Road, Kowloon Bay,
Kowloon, Hong Kong
Phone: +852 3746 5902
Fax: +852 23 63 19 55
Email: info@medisana.com.hk
Contact: Kitty Yeung

UNITED ARAB EMIRATES

medisana MENA FZ-LLC.

Dubai Science Park
North Tower, Office 705
Phone: +971 50 476 8055
Email: w.luka@medisana-global.com
Internet: www.medisana.com
Contact: Walid Luka

SERVICE CENTER

Phone: +49 2131 405 88-0
Email: servicecenter@medisana.com

SPAIN

medisana Healthcare S. L.

Portal de l'Angel, 36
08002 Barcelona
Phone: +34 938 33 50 52
Fax: +34 938 33 36 22
Email: info@medisana.es
Internet: www.medisana.es
Contact: Miguel Méndez

RUSSIA

000 medisana Rus

ul. Nagornaja d.20
RU-117186 Moscow
Phone: +7 495 640 26 47
Email: info@medisana.su
Internet: www.medisana.su
Contact: Andrey Usov

GERMANY

Promed GmbH

Kosmetische Erzeugnisse
Lindenweg 11
D-82490 Farchant
Phone: +49 8821 9621-0
Fax: +49 8821 9621-21
Email: info@promed.de
Internet: www.promed.de
Contact: Andreas Schuster

Imprint

medisana GmbH
Carl-Schurz-Str. 2, 41460 Neuss, Germany

Phone: +49 2131 / 36 68-0
Fax: +49 2131 / 36 68-50 95

Email: info@medisana.com
Internet: www.medisana.com

General Manager: Thomas Hollefeld
Kwok Shing Ng

Registration court: Amtsgericht Neuss
Register number: HRB 16348
VAT ID number: DE 123375019

Modifications / Terms and Conditions

Subject to technical modifications and further developments, also without prior notice. We accept no liability for printing errors, typesetting errors or other mistakes. Our current Terms and Conditions can be found at any time on our website.

Copyright

Copyright over the information published in the catalogue (text, pictures, graphics) and their structure is owned by medisana GmbH, and in some cases by third parties. Duplication, distribution, reproduction, dissemination, storage and other use of this information without our written consent is expressly prohibited.

© Copyright 2022, medisana GmbH, Germany



medisana GmbH
Carl-Schurz-Str. 2
41460 Neuss
Germany

Phone: +49 2131 36 68-0
Fax: +49 2131 36 68-50 95

info@medisana.com
www.medisana.com



JULY 2022