

# medisana®

## PRODUCT RANGE 2020



YOUR HEALTH IN  
GOOD HANDS

MOBILE HEALTH HEALTH CONTROL THERAPY & HEALTHY HOME  
WELLNESS SPORT BODY CARE

**medisana**®

# CONTENT

## HEALTH CONTROL

Blood pressure monitors	4
Body analysis scales	10
Personal scales	14
Thermometer	18
Blood glucose meter	20
Pulse oximeter	21

## THERAPY & HEALTHY HOME

Light therapy	22
Inhalers	25
Air & Aroma	28
Lice combs	32
Sleep therapy	33
Pain therapy	33

## WELLNESS

Relax chairs	34
Massage seat covers	35
Neck massager	40
Massage cushion	42
Massage devices	43
Foot massager	44
Bath spa	46
Heat pad	48
Heating blanket	50
Foot warmer	52
Heated underblanket	54

## SPORT

Fascia rolls	56
Body trainer	58
Activity tracker	59

## BODY CARE

Hand & Foot care	60
Cosmetic mirror	63
Body care	64
Accessories	66
Contact	71
Imprint	72





NEW

**BPW 100 CONNECT**  
Blood pressure watch

**2in1 smartwatch and blood pressure watch**

- Blood pressure measurement watch with Bluetooth and VitaDock+ App
- Faster measurement results thanks to inflation technology
- High quality E-ink display with low power consumption



**Technical data**

- Reliable use and accurate blood pressure measurement
- Measurement function adjustable via VitaDock+ App
- Adjustable time and number of blood pressure measurements per day
- Wrist circumferences: 135 – 215 mm
- Display size: 27.6 × 27.6 mm
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- Certified medical device\*\*\*

**Size:** 4.7×4.1×13.5 cm  
**Weight:** approx. 0.065 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-299 mmHg  
Pulse: 40-180 beats/minute  
**Batteries:** 3.7V 210mAh rechargeable Li-ion battery  
**Article No.:** 51215  
**EAN:** 40 15588 51215 5  
**Warranty:** 3 Years

\*iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0



NEW

**BU 570 CONNECT**  
Upper arm blood pressure monitor

**With arrhythmia indicator**

- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- Extra large and easy-to-read display
- Accurate blood pressure measurement on the upper arm with WHO classification



**Technical data**

- Big LCD display: 7.3 x 6 cm with backlight
- 60 memory slots for each of 2 users
- Calculates average value of data of the last 3 results
- Arrhythmia indicator
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- Certified medical device\*\*\*

**Size:** approx. 9.3×13×3.4 cm  
**Weight:** approx. 0.9 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-299 mmHg  
Pulse: 40-180 Beats/minute  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 51206  
**EAN:** 40 15588 51206 3  
**Warranty:** 3 Years

\*\*World Health Organisation  
\*\*\*Certification according to MDD (Medical Device Directive)



NEW

**BU 565**  
Upper arm blood pressure monitor

**Particularly flat design**

- Faster measurement results thanks to inflation technology – The blood pressure is measured during the inflation of the cuff
- Extra large and easy-to-read display
- Accurate blood pressure measurement on the upper arm with WHO classification



**Technical data**

- Big LCD display: 7.3 x 6 cm with backlight
- 60 memory slots for each of 2 users
- Calculates average value of data of the last 3 results
- Arrhythmia indicator
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Certified medical device\*\*\*

**Size:** approx. 9.3×13×3.25 cm  
**Weight:** approx. 0.19 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-299 mmHg  
Pulse: 40-180 beats/minute  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 51202  
**EAN:** 40 15588 51202 5  
**Warranty:** 3 Years



NEW

**BU 560 CONNECT**  
Upper arm blood pressure monitor

**With integrated cuff**

- All in one device – Blood pressure monitor with integrated cuff without hose or cable
- Cuff tightness detection
- 120 memory slots for one user



**Technical data**

- Calculates average value of data of the last 3 results
- Arrhythmia indicator
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Body movement indication
- Cuff size 22 – 42 cm
- With rechargeable Lithium battery 500mHA and USB port
- LCD size: 4.0×3.5 cm
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- Certified medical device\*\*\*

**Size:** 12×5.5×3.4 cm (without cuff)  
**Weight:** approx. 0.23 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-280 mmHg  
Pulse: 40-199 beats/minute  
**Power supply:** Lithium Battery 500 mHA rechargeable  
**Article No.:** 51196  
**EAN:** 40 15588 51196 7  
**Warranty:** 3 Years

**Blood pressure monitors**



NEW

**BU 555**  
Upper arm blood pressure monitor

**Very slim design**

- All in one device – Blood pressure monitor with integrated cuff without hose or cable
- Cuff tightness detection
- 120 memory slots for one user



**Technical data**

- Calculates average value of data of the last 3 results
- Arrhythmia indicator
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Body movement indication
- Cuff size 22 – 42 cm
- With rechargeable Lithium battery 500mHA and USB port
- LCD size: 4.0×3.5 cm
- Certified medical device\*\*\*

**Size:** 12×5.5×3.4 cm (without cuff)  
**Weight:** approx. 0.23 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-280 mmHg  
Pulse: 40-199 beats/minute  
**Power supply:** Lithium Battery 500 mHA rechargeable  
**Article No.:** 51192  
**EAN:** 40 15588 51192 9  
**Warranty:** 3 Years



Bluetooth®

**BU 540 CONNECT**  
Upper arm blood pressure monitor

**Top price-quality ratio!**

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 250 memory slots for each of 2 users
- Very large digits for easy reading
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large upper arm cuff (arm circumferences 22-42 cm)
- Incl. storage bag
- Certified medical device\*\*\*



**Size:** approx. 15×9.5×4.1 cm  
**Weight:** approx. 0.23 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 40-260 mmHg  
 Pulse: 40-180 beats/minute  
**Batteries:** 4×1.5V (AA)  
**Article No.:** 51182  
**EAN:** 40 15588 51182 0  
**Warranty:** 3 Years



Bluetooth®

**BU 530 CONNECT**  
Upper arm blood pressure monitor

**Connect device**

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Measurement for guests without saving measured values
- Battery operation (power adapter optionally available)
- Upper arm cuff for 22-36 cm arm circumference
- Incl. storage bag
- Certified medical device\*\*\*

**Size:** approx. 15×10×6 cm  
**Weight:** approx. 0.4 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-299 mmHg  
 Pulse: 40-199 beats/minute  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 51174  
**EAN:** 40 15588 51174 5  
**Warranty:** 3 Years



**BU 535 VOICE**  
Upper arm blood pressure monitor

**With voice output**

- Accurate blood pressure measurement on the upper arm
- Voice output in selectable languages:  
 Art. 51178: DE, GB, NL, FR, IT, TR  
 Art. 51179: GB, FR, RU ES, GR, AR (arabic)
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Upper arm cuff for 22 - 36 cm arm circumference
- Incl. storage bag
- Certified medical device\*\*\*

**Size:** approx. 15×9.5×4.1 cm  
**Weight:** approx. 0.23 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 40-260 mmHg  
 Pulse: 40-180 beats/minute  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 51178/51179  
**EAN:** 40 15588 51178 3  
 40 15588 51179 3  
**Warranty:** 3 Years



**BU 535**  
Upper arm blood pressure monitor

**With XL display**

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Extra large, easy-to-read display
- Ultra slim design
- Upper arm cuff for 22 - 36 cm arm circumference
- Incl. storage bag
- Certified medical device\*\*\*

**Size:** approx. 15×9.5×4.1 cm  
**Weight:** approx. 0.23 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 40-260 mmHg  
 Pulse: 40-180 beats/minute  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 51176  
**EAN:** 40 15588 51176 9  
**Warranty:** 3 Years



**BU 512**  
Upper arm blood pressure monitor

**Handy upper arm cuff**

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 90 memory slots for each of 2 users
- Very large digits for easy reading (up to 22 mm in height)



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\*\* evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- Battery operation
- Upper arm cuff for 22 - 36 cm arm circumference
- Incl. storage bag
- Certified medical device\*\*\*



**Size:** approx. 13×11×6 cm  
**Weight:** approx. 0.3 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 40-299 mmHg  
 Pulse: 40-199 beats/minute  
**Batteries:** 4×1.5V (AA)  
**Article No.:** 51162  
**EAN:** 40 15588 51162 2  
**Warranty:** 3 Years

\* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0

\*\*World Health Organisation  
 \*\*\*Certification according to MDD (Medical Device Directive)







**MTP PRO**  
Upper arm blood pressure monitor

**MTS**  
Upper arm blood pressure monitor

**Our million-fold bestseller**

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 99 memory spaces for each of 2 users
- Higher accuracy of measurements for severely fluctuating blood pressure thanks to patented 3 MAM technology
- 10 years of warranty

**Easy handling**

- Accurate blood pressure measurement on the upper arm
- Arrhythmia indicator
- 60 memory spaces for each of 2 users
- Very large digits for easy reading



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Extra large, easy-to-read display
- Measurements are classified according to the WHO\* evaluation system using a traffic-light colour scale
- Calculates average value of all stored data
- 1-button operation
- Battery operation (power adapter optionally available)
- Upper arm cuff for 22 - 32 cm arm circumference
- Incl. storage bag
- Incl. blood pressure pass
- Certified medical device\*\*

**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\* evaluation system using a traffic-light colour scale
- Automatic switch-off after 1 min.
- Battery operation (power adapter optionally available)
- Upper arm cuff for 22 - 30 cm arm circumference
- Incl. storage bag
- Certified medical device\*\*



**Size:** approx. 12×16.3×6.6 cm  
**Weight:** approx. 0.5 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 30-280 mmHg  
 Pulse: 40-200 beats/minute  
**Batteries:** 4×1.5V (AA)  
**Article No.:** 51090  
**EAN:** 40 15588 51090 8  
**Warranty:** 10 Years

**Size:** approx. 19.5×22.5×26.5 cm  
**Weight:** approx. 0.3 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-300 mmHg  
 Pulse: 30-180 beats/minute  
**Batteries:** 4×1.5V (AA)  
**Article No.:** 51152  
**EAN:** 40 15588 51152 3  
**Warranty:** 3 Years



**BW 335**  
Wrist blood pressure monitor

**BW 320**  
Wrist blood pressure monitor

**BW 315**  
Wrist blood pressure monitor

**Flat design**

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Very large digits for easy reading
- Ultra slim design

**Elegant design for a top price-quality ratio**

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Extra large, easy-to-read display

**Our bestseller**

- Accurate blood pressure measurement on the wrist
- Arrhythmia indicator
- 120 memory spaces for each of 2 users
- Extra large, easy-to-read display



**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Cuff for wrist circumferences 14 - 19.5 cm
- Incl. storage box
- Certified medical device\*\*

**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Battery level indicator
- Cuff for wrist circumferences 12.5 - 21.5 cm
- Certified medical device\*\*

**Technical data**

- Display items: systolic, diastolic, pulse, date, time
- Measurements are classified according to the WHO\* evaluation system using a traffic-light colour scale
- Calculates average value of data of the last 3 results
- Battery level indicator
- Cuff for wrist circumferences 12.5 - 21.5 cm
- Certified medical device\*\*

**Size:** approx. 8.9×6×3.1 cm  
**Weight:** approx. 0.07 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 0-299 mmHg  
 Pulse: 40-180 beats/minute  
**Batteries:** 2×1.5V (AAA)  
**Article No.:** 51077  
**EAN:** 40 15588 51077 9  
**Warranty:** 3 Years

**Size:** approx. 7.5×6.4×2.9 cm  
**Weight:** approx. 0.117 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 30-280 mmHg  
 Pulse: 40-200 beats/minute  
**Batteries:** 2×1.5V (AAA)  
**Article No.:** 51074  
**EAN:** 40 15588 51074 8  
**Warranty:** 3 Years

**Size:** approx. 6.3×6.3×2.9 cm  
**Weight:** approx. 0.11 kg  
**Measuring method:** Oscillometric  
**Measuring range:** Blood pressure: 30-280 mmHg  
 Pulse: 40-200 beats/minute  
**Batteries:** 2×1.5V (AAA)  
**Article No.:** 51072  
**EAN:** 40 15588 51072 4  
**Warranty:** 3 Years

\*World Health Organisation

\*\*Certification according to MDD (Medical Device Directive)





NEW

**BS 500 CONNECT**  
WiFi body analysis scale

**With WiFi data transfer**

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Heart rate measurement with every weight measurement
- Smart LED lights indicate weight change:  
Red: weight gain  
Green: weight loss  
Blue: no weight change



**Technical data**

- Convenient WiFi data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- High-quality and invisible ITO electrodes
- Automatic recognition up to 8 users on the scale
- Convenient „step-on“ switch
- Invisible LED display (only appears during the weighing)
- Switchable between kg, lb and st
- Body analysis not applicable for persons with implants (e.g. pacemakers)



Bluetooth®

**BS 450 CONNECT**  
Body analysis scale

**Connect device with LED display**

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- High-quality ITO electrodes
- Convenient „step-on“ activation
- Invisible LED display that only lights up when in use
- Switchable between kg, lb and st
- 4 high-precision strain gauge sensors for accurate results
- 30 memory slots for each of the 8 users
- Sensor touch buttons
- Closed underbody for better a stability
- Body analysis not applicable for persons with implants (e.g. pacemakers)



Bluetooth®

**BS 440 CONNECT**  
Body analysis scale

**Connect device with ITO electrodes**

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- High-quality ITO electrodes
- Convenient „step-on“ activation
- Large easy-to-read LCD display with backlight
- Switchable between kg, lb and st
- 4 high-precision strain gauge sensors for accurate results
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)



Bluetooth®

**BS 430 CONNECT**  
Body analysis scale

**With backlight**

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- High-quality stainless steel electrodes
- Convenient „step-on“ activation
- Large easy-to-read LCD display with backlight
- Switchable between kg, lb and st
- 4 high-precision strain gauge sensors for accurate results
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)



**Size:** approx. 30.7×30.7×24 cm  
**Weight:** approx. 1.8 kg  
**Display size:** 6.2×6.2 cm  
**Weight display:** 5 to 180 kg  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40443 (white), 40442 (black)  
**EAN:** 40 15588 40443 6 (white)  
40 15588 40442 9 (black)  
**Warranty:** 3 Years

**Size:** approx. 31.9×31.9×2.3 cm  
**Weight:** approx. 1.7 kg  
**Display size:** 7.4×2.9 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 40511 (white), 40510 (black)  
**EAN:** 40 15588 40511 2 (white)  
40 15588 40510 5 (black)  
**Warranty:** 3 Years

**Size:** approx. 32×32×2.7 cm  
**Weight:** approx. 2 kg  
**Display size:** 7×4.9 cm  
**Weight display:** 5 to 180 kg  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 40423  
**EAN:** 40 15588 40423 8  
**Warranty:** 3 Years



**Size:** approx. 32.5×32.5×2.5 cm  
**Weight:** approx. 2 kg  
**Display size:** 7×4.9 cm  
**Weight display:** 5 to 180 kg  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 40422  
**EAN:** 40 15588 40422 1  
**Warranty:** 3 Years

\* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0







Bluetooth®

**BS 444 CONNECT**  
Body analysis scale

**Our Bestseller**

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- Determination of the BMI value
- Integrated calorie requirement analysis (BMR)
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- High-quality stainless steel electrodes
- Convenient „step-on“ activation
- Large easy-to-read LCD display
- Switchable between kg, lb and st
- 4 high-precision strain gauge sensors for accurate results
- 30 memory spaces for each of 8 users
- Body analysis not applicable for persons with implants (e.g. pacemakers)



Bluetooth®

**TARGETSCALE® 3**  
Body analysis scale

**Athlete mode for athletes**

- Precise determination of weight (up to 180 kg) Body Mass Index (BMI), kcal per day, calculates percentage of body fat and water as well as muscle and bone mass
- Automatic recognition of up to 8 users
- Data transfer via Bluetooth® to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- Lighted rings depict the difference between current and target weight
- Automatic recognition and individual definition of target weight
- Large easy-to-read LCD display with white luminous digits
- Athlete mode for athletes
- Stores up to 240 measurements in the scale memory (automatic data transfer to smartphone)
- Switchable between kg, lb and st
- Body analysis not applicable for persons with implants (e.g. pacemakers)



**BS 465**  
Body analysis scale

**With color indication**

- Measurement of weight, body fat, body water, muscle percentage and bone weight
- LCD display shows by color the difference between the current weight and the target weight
- Green: under target weight  
Blue: target weight  
Red: above target weight



**Technical data**

- Automatic recognition of up to 8 users
- 9 memory slots for 9 person
- Individual adjustment of the target weight possible
- Four high-precision strain gauge sensors for accurate readings
- 6 mm tempered safety glass platform
- Athlete mode for athletes
- Step on & auto-off
- Automatic switch-off after 8 seconds
- Integrated calorie requirement analysis
- Switchable between kg, lb and st
- Body analysis not applicable for persons with implants (e.g. pacemakers)



**Size:** approx. 30×30×2.3 cm  
**Weight:** approx. 1.6 kg  
**Display size:** 7.3×4.8 cm  
**Weight display:** 5 to 180 kg  
**Graduation:** 100g, 0.2lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40444  
**EAN:** 40 15588 40444 3  
**Warranty:** 3 Years

**Size:** approx. ø 36×4.75 cm  
**Weight:** approx. 2 kg  
**Display size:** 7×5 cm  
**Weight display:** 5 to 180 kg  
**Graduation:** 100g, 0.2lb, 1/4 lb  
**Batteries:** 4×1.5V (AA)  
**Article No.:** 40414  
**EAN:** 40 15588 40414 6  
**Warranty:** 3 Years

**Size:** approx. 30×30×3 cm  
**Weight:** approx. 1.7 kg  
**Display size:** 4.4×9.1 cm  
**Weight display:** 5 to 180 kg  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 40483 (white), 40484 (black)  
**EAN:** 40 15588 40483 2 (white)  
40 15588 40484 9 (black)  
**Warranty:** 3 Years

**Body analysis scales**



**BS 460**  
Body analysis scale

**With removeable display**

- Measurement of weight, body fat, body water, muscle and bone mass
- Integrated calorie analysis (BMR)
- Maximum capacity: 180 kg
- Removable display with stand and wall mount
- Convenient „step-on“ activation



**Technical data**

- Modern blue illumination when in use (shiny cross magic light)
- Sensor touch buttons
- Tread area of 6 mm safety glass with rounded corners
- 4 high-precision strain gauge sensors for accurate results
- Displays time
- 10 memory slots
- Energy saving automatic switch-off after 10 seconds
- Overload indication
- Low battery warning
- Switchable between kg, lb and st
- Body analysis not applicable for persons with implants (e.g. pacemakers)

**Size:** approx. 31.9×31.9×2.8 cm  
**Weight:** approx. 2 kg  
**Display size:** 9×4.5 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 7×1.5V (AAA)  
**Article No.:** 40482  
**EAN:** 40 15588 40482 5  
**Warranty:** 3 Years



\* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0





NEW

**PS 436**  
Design glass scale

**Design „Golden Sprinkle“**

- Invisible LED display (only appears during the weighing)
- Large easy-to-read LED display (28 mm)
- Convenient „step-on“ switch
- Energy-saving automatic cut-off after 10 seconds



**Technical data**

- Treads of the standing area made of 6 mm safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Extra slim design
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 30×30×2.1 cm  
**Weight:** approx. 1.5 kg  
**Display size:** 7×5 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40520  
**EAN:** 40 15588 40520 4  
**Warranty:** 3 Years



NEW

**PS 437**  
Design glass scale

**Design „Silver Lining“**

- Invisible LED display (only appears during the weighing)
- Large easy-to-read LED display (28 mm)
- Convenient „step-on“ switch
- Energy-saving automatic cut-off after 10 seconds



**Technical data**

- Treads of the standing area made of 6 mm safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Extra slim design
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 30×30×2.1 cm  
**Weight:** approx. 1.5 kg  
**Display size:** 7×5 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40522  
**EAN:** 40 15588 40522 8  
**Warranty:** 3 Years



NEW

**PS 438**  
Design glass scale

**Design „Enjoy“**

- Invisible LED display (only appears during the weighing)
- Large easy-to-read LED display (28 mm)
- Convenient „step-on“ switch
- Energy-saving automatic cut-off after 10 seconds



**Technical data**

- Treads of the standing area made of 6 mm safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Extra slim design
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 30×30×2.1 cm  
**Weight:** approx. 1.5 kg  
**Display size:** 7×5 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40524  
**EAN:** 40 15588 40524 2  
**Warranty:** 3 Years



NEW

**PS 439**  
Design glass scale

**Design „Pineapple“**

- Invisible LED display (only appears during the weighing)
- Large easy-to-read LED display (28 mm)
- Convenient „step-on“ switch
- Energy-saving automatic cut-off after 10 seconds



**Technical data**

- Treads of the standing area made of 6 mm safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Extra slim design
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 30×30×2.1 cm  
**Weight:** approx. 1.5 kg  
**Display size:** 7×5 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40526  
**EAN:** 40 15588 40526 6  
**Warranty:** 3 Years



**PS 435**  
Design glass scale

**Noble design**

- Invisible LED display (only appears during the weighing)
- Large easy-to-read LED display (28 mm)
- Convenient „step-on“ switch
- Energy-saving automatic cut-off after 10 seconds



**Technical data**

- Treads of the standing area made of 6 mm safety glass with rounded corners
- Four high-precision strain gauge sensors for accurate readings
- Extra slim design
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 30×30×2.1 cm  
**Weight:** approx. 1.5 kg  
**Display size:** 7×5 cm  
**Weight display:** bis 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40434  
**EAN:** 40 15588 40434 4  
**Warranty:** 3 Years



NEW

**PS 500**  
Glass personal scale

**Trendy glitter look**

- Glitter and glossy speckled surface
- Standing area made of 6 mm safety glass
- High-precision strain gauge sensors for accurate readings



**Technical data**

- Convenient „step-on“ switch
- Low battery indication
- Overload indication
- Switchable between kg, lb and st

**Size:** approx. 30×30×2.1 cm  
**Weight:** approx. 1.5 kg  
**Display size:** 2.9×7.3 cm  
**Weight display:** up to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 1×CR2032  
**Article No.:** 40490 (gold), 40491 (grey)  
**EAN:** 40 15588 40490 0 (gold)  
 40 15588 40491 7 (grey)  
**Warranty:** 3 Years



**PS 440**  
Bamboo personal scale

**Magic LED display**

- Decorative personal scale made of bamboo
- Invisible LED display with red luminous numbers
- Energy-saving automatic switch-off



**Technical data**

- 4 high-precision strain gauge sensors for accurate results
- Convenient „step-on“ activation
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 30×30×3.2 cm  
**Weight:** approx. 1 kg  
**Character height:** 2.8 cm  
**Weight display:** 5 to 180 kg, 396 lb, 28 st 4 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40544  
**EAN:** 40 15588 40544 0  
**Warranty:** 3 Years



**PS 425**  
Personal scale with voice output

**Our talking scale**

- Voice output in 4 languages: DE, GB, FR, ES
- Easy-to-read LCD display
- Convenient „step-on“ activation
- Energy-saving automatic switch-off



**Technical data**

- Tread area made of 6 mm safety glass with rounded corners
- 4 high-precision strain gauge sensors for accurate results
- Switchable between kg and lb
- Overload indication
- Low battery display

**Size:** approx. 31×31×2.4 cm  
**Weight:** approx. 1.8 kg  
**Character height:** 2.8 cm  
**Weight display:** 4 to 180 kg  
**Graduation:** 100 g, 0.2 lb  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 40448  
**EAN:** 40 15588 40448 1  
**Warranty:** 3 Years



**PS 470**  
XL glass personal scale

**Particularly resilient**

- Extra large, easy-to-read display (71 mm high)
- Tread area made of extra white safety glass with rounded corners
- Energy-saving automatic switch-off



**Technical data**

- 4 high-precision strain gauge sensors for accurate results
- Convenient „step-on“ activation
- Switchable between kg, lb and st
- Overload indication
- Low battery display

**Size:** approx. 34×34×2.5 cm  
**Weight:** approx. 2 kg  
**Display size:** 9.3×7.7 cm  
**Weight display:** up to 250 kg, 550 lb, 39 st 6 lb  
**Graduation:** 100 g, 0.2 lb, 1/4 lb  
**Batteries:** 1×CR2032  
**Article No.:** 40547  
**EAN:** 40 15588 40547 1  
**Warranty:** 3 Years



NEW

**PS 100**  
Analogue personal scale

**Retro-Style 2.0**

- Trendy retro design
- Large full-view analogue dial
- Non-slip tread area



**Technical data**

- Battery-free, mechanical operation
- Solid metal casing

**Size:** approx. 29×43×8.5 cm  
**Weight:** approx. 3.2 kg  
**Weight display:** up to 150 kg  
**Graduation:** 500 g  
**Article No.:** 40498 (red), 40499 (black)  
**EAN:** 40 15588 40498 6 (red)  
 40 15588 40499 3 (black)  
**Warranty:** 3 Years



**PSD**  
Analogue personal scale

**Trendy retro design**

- Trendy retro design
- Battery-free, mechanical operation
- Large analogue full view scale



**Technical data**

- Non-slip tread area
- Solid metal casing

**Size:** approx. 28.5×42.5×8.5 cm  
**Weight:** approx. 3.6 kg  
**Weight display:** up to 150 kg  
**Graduation:** 500 g  
**Article No.:** 40461  
**EAN:** 40 15588 40461 0  
**Warranty:** 3 Years





<p><b>TM 780</b> Infrared multifunctional thermometer</p> <p><b>Non-contact measurement</b></p> <ul style="list-style-type: none"> <li>Precise non-contact measurement of the body temperature</li> <li>Additional function: measurement of ambient temperature, liquids and surfaces (2-100 °C)</li> <li>Visual and acoustic fever alarm with colour change – green: normal temperature, red: raised temperature</li> <li>Distance detection to indicate right distance for measurement (error signal at &gt; 8cm)</li> </ul>	<p><b>TM 760</b> Infrared multifunctional thermometer</p> <p><b>With LED tracking light</b></p> <ul style="list-style-type: none"> <li>Precise non-contact measurement of the body temperature</li> <li>Additional function: measurement of ambient temperature and surfaces (0-100 °C)</li> <li>Visual and acoustic fever alarm with colour change – green: normal temperature, red: raised temperature</li> <li>Extra large, easy-to-read display</li> </ul>	<p><b>TM 750 CONNECT</b> Multifunctional thermometer</p> <p><b>Connect device</b></p> <ul style="list-style-type: none"> <li>Precise measurement of the body temperature in the ear and on the forehead in seconds</li> <li>Additional function: measurement of ambient temperature, liquids and surfaces (0-100 °C)</li> <li>Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online**</li> </ul>
--	--	---




<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>Acoustic signal at the end of temperature reading</li> <li>Automatic display of the last measured temperature</li> <li>Automatic power off after 60 seconds</li> <li>Low battery indication</li> <li>Conversion from °C to °F</li> <li>Incl. storage case</li> <li>Certified medical device*</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>LED Tracking light to indicate the correct point of measurement</li> <li>Acoustic signal at the beginning and end of temperature reading (can be switched off to avoid disturbance during sleep)</li> <li>Automatic storage of the last 30 measurement results</li> <li>Display of date and time</li> <li>Conversion from °C to °F</li> <li>Ergonomic and light-weight design</li> <li>Certified medical device*</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>Automatic storage of the last 30 measurement results</li> <li>Visual fever alarm with colour change – green: normal temperature, red: raised temperature</li> <li>Acoustic signal at the beginning and end of temperature reading</li> <li>Display of date and time</li> <li>Incl. storage box</li> <li>Certified medical device*</li> </ul>
---	---	--

<p><b>Size:</b> 14.2×3.7×3.5 cm <b>Weight:</b> approx. 0.08 kg <b>Measuring range:</b> 2-100 °C <b>Measuring accuracy:</b> Forehead: ± 0.3 °C (34-43 °C) Object: 2-100 °C ± 1 °C <b>Batteries:</b> 2×1.5V (AAA) <b>Article No.:</b> 76175 <b>EAN:</b> 40 15588 76175 1 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> 15.8×4.2×4.3 cm <b>Weight:</b> approx. 0.067 kg <b>Measuring range:</b> Body: 34-43 °C Object: 0-100 °C <b>Measuring accuracy:</b> Forehead: ± 0.3 °C (34-42 °C) Object: ± 1 °C (0-100 °C) <b>Batteries:</b> 1×CR2032 <b>Article No.:</b> 76121 <b>EAN:</b> 40 15588 76121 8 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 13.4×4.2×4.7 cm <b>Weight:</b> approx. 0.065 kg <b>Measuring range:</b> Body: 34-43 °C Object: 0-100 °C <b>Measuring accuracy:</b> Forehead: ± 0.3 °C (34-43 °C) Ear: ± 0.3 °C (35.5-42 °C) Object: ≤30 ±1.5 °C / &gt;30 ±5% <b>Batteries:</b> 2×1.5V (AAA) <b>Article No.:</b> 76145 <b>EAN:</b> 40 15588 76145 4 <b>Warranty:</b> 3 Years</p>
---	---	---



<p><b>TM 750</b> Multifunctional thermometer</p> <p><b>Extended memory function</b></p> <ul style="list-style-type: none"> <li>Precise measurement of the body temperature in the ear and on the forehead in seconds</li> <li>Additional function: measurement of ambient temperature, liquids and surfaces (10-100 °C)</li> </ul>	<p><b>TM 700</b> Thermometer</p> <p><b>Our gentle classic</b></p> <ul style="list-style-type: none"> <li>Precise measurement of the body temperature: oral, axillary, rectal</li> <li>Flexible tip</li> <li>Measurement in 10 seconds</li> <li>Water resistant</li> </ul>	<p><b>FTC</b> Thermometer</p> <p><b>Top price-quality ratio</b></p> <ul style="list-style-type: none"> <li>Precise measurement of the body temperature: oral, axillary, rectal</li> <li>Acoustic signal at the beginning and end of temperature reading</li> <li>Water resistant</li> </ul>
--	---	---



<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>Automatic storage of the last 30 measurement results</li> <li>Visual fever alarm with colour change – green: normal temperature, red: raised temperature</li> <li>Acoustic signal at the beginning and end of temperature reading</li> <li>Display of date and time</li> <li>Incl. storage box</li> <li>Certified medical device*</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>Acoustic fever alarm</li> <li>Acoustic signal at the beginning and end of temperature reading</li> <li>Automatic storage of the last measurement result</li> <li>Automatic switch-off after 10 minutes</li> <li>Digital LCD display</li> <li>Conversion from °C to °F</li> <li>Easy battery replacement</li> <li>Incl. storage box</li> <li>Certified medical device*</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>Acoustic fever alarm</li> <li>Automatic storage of the last measurement result</li> <li>Automatic switch-off after 8 minutes</li> <li>Digital LCD display</li> <li>Easy battery replacement</li> <li>Incl. storage box</li> <li>Certified medical device*</li> </ul>
--	--	--

<p><b>Size:</b> approx. 13.4×4.2×4.7 cm <b>Weight:</b> approx. 0.065 kg <b>Measuring range:</b> Body: 34-43 °C Object: 0-100 °C <b>Measuring accuracy:</b> Forehead: ± 0.3 °C (34-43 °C) Ear: ± 0.3 °C (35.5-42 °C) Object: ≤30 ±1.5 °C / &gt;30 ±5% <b>Batteries:</b> 2×1.5V (AAA) <b>Article No.:</b> 76140 <b>EAN:</b> 40 15588 76140 9 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 13.9×2.2×1.3 cm <b>Weight:</b> approx. 0.011 kg <b>Measuring range:</b> 32-42.9 °C <b>Scaling:</b> 0.1 °C <b>Batteries:</b> 1×1.5V (AG3/LR41) <b>Article No.:</b> 77040 <b>EAN:</b> 40 15588 77040 1 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 12.8×1×1.9 cm <b>Weight:</b> approx. 0.01 kg <b>Measuring range:</b> 32-42 °C <b>Scaling:</b> 0.1 °C <b>Batteries:</b> 1×1.5V (AG3/LR41) <b>Article No.:</b> 77030 <b>EAN:</b> 40 15588 77030 2 <b>Warranty:</b> 3 Years</p>
---	---	--

\*\* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0

\* Certification according to MDD (Medical Device Directive)



NEW

Bluetooth®

**GLUCOWATCH CONNECT**  
Blood glucose meter watch

**Blood glucose meter watch with Bluetooth® and VitaDock+ App**

- Healthcare management function: step counter and sleep monitor
- Glucose level warning by color system indication
- Stainless steel watch housing with detachable display
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- Measuring unit individually adjustable: mg/dL or mmol/L
- Standard smartwatch features: displays incoming calls and messages
- Reliable use and blood pressure measurement
- Water-proof test strip slot design to fully protect blood glucose measurement quality, and give the most accurate results
- Precise measurement accuracy according to specialized temperature calibration algorithm: avoid interference from body and environment temperature



Mini-to-go set incl. 20 sealed test strips, lancing device und 10 lancets

**Size:** approx. 4.5×1.2 cm  
**Weight:** approx. 0.084 kg  
**Batteries:** 435mAh/h rechargeable lithium battery  
**Power output:** 5V/1A adaptor  
**Cable length:** 120 cm  
**Article No.:** 79150  
**EAN:** 40 15588 79150 5  
**Warranty:** 3 Years



Bluetooth®

**MEDITOUCH 2 CONNECT**  
Blood glucose meter

**Connect device**

- Precise measurement accuracy according to the new ISO 15197: 2015 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- USB port for data transfer to VitaDock Online
- Pre and post meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0.6 µl amount of blood
- Extra large, easy-to-read display
- Incl. diabetic diary
- Certified medical device\*\*



Starter set incl. 10 test strips, lancing device, 10 lancets, AST cap, control solution and bag

**Size:** approx. 5×9.8×1.1 cm  
**Weight:** approx. 0.047 kg  
**Batteries:** 2×3V (CR 2032)  
**dual Article No.:** 79048  
**EAN:** 40 15588 79048 5  
**Warranty:** 3 Years



Bluetooth®

**MEDITOUCH 2 CONNECT CARBON**  
Blood glucose meter

**Connect device in a carbon look**

- Precise measurement accuracy according to the new ISO 15197: 2015 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- USB port for data transfer to VitaDock Online
- Measuring unit individually adjustable: mg/dL and mmol/L
- Pre and post meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0.6 µl amount of blood
- Extra large, easy-to-read display
- Certified medical device\*\*



Starter set incl. 10 test strips, lancing device, 10 lancets, AST cap, control solution and bag

**Size:** approx. 5×9.8×1.1 cm  
**Weight:** approx. 0.047 kg  
**Batteries:** 2×3V (CR 2032)  
**dual Article No.:** 79062  
**EAN:** 40 15588 79062 1  
**Warranty:** 3 Years



**MEDITOUCH 2**  
Blood glucose meter

**Our classic**

- Precise measurement accuracy according to the new ISO 15197: 2015 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- Acoustic test reminders
- 480 memory slots



**Technical data**

- USB port for data transfer to VitaDock Online
- Calculates average values for 7/14/30/90 days
- Pre and post meal markers
- Automatic ketone warning
- Measuring time within 5 seconds
- 0.6 µl amount of blood
- Extra large, easy-to-read display
- Certified medical device\*\*



Starter set incl. 10 test strips, lancing device, 10 lancets, AST cap, control solution and bag

**Size:** approx. 5×9.8×1.1 cm  
**Weight:** approx. 0.03 kg  
**Batteries:** 2×3V (CR 2032)  
**mg/dL Article No.:** 79030  
**EAN:** 40 15588 79030 0  
**mmol/L Article No.:** 79034  
**EAN:** 40 15588 79034 8  
**Warranty:** 3 Years



Bluetooth®

**PM 100 CONNECT**  
Pulse oximeter

**Now as connect device**

- Pulse oximeter to measure the saturation of oxygen in blood (SpO2%) and the heart rate (pulse)
- View measurements on OLED display
- Display brightness adjustable in 10 levels
- Simple „one touch“ operation – ideal for on the go
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- Oxygen saturation (SpO2) and pulse rate can be displayed as bars or waves
- 6 different views
- Automatic switch-off after 8 seconds
- Low battery display
- Incl. carrying strap
- Certified medical device\*\*



**PM 100**  
Pulse oximeter

**Simple „one touch“ operation**

- Pulse oximeter to measure the saturation of oxygen in blood (SpO2 %) and the heart rate (pulse)
- View measurements on OLED display
- Display brightness adjustable in 10 levels
- Simple „one touch“ operation – ideal for on the go



**Technical data**

- Oxygen saturation (SpO2) and pulse rate can be displayed as bars or waves
- 6 different views
- Automatic switch-off after 8 seconds
- Low battery display
- Incl. carrying strap
- Certified medical device\*\*

**Size:** approx. 5.8×3.4×3.5 cm  
**Weight:** approx. 0.03 kg  
**Measuring range:** SpO2: 70% ~ 100%  
 Pulse: 30-235 beats/minute  
**Batteries:** 2×1.5V (AAA)  
**Article No.:** 79455  
**EAN:** 40 15588 79455 1  
**Warranty:** 3 Years

\* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0

\*\*Certification according to MDD (Medical Device Directive)






**IR 885**  
Infrared lamp

**With timer function**

- Intense infrared light to relax muscles or treat colds
- Stimulates blood circulation
- Relieves pain
- Durable lamp (halogen)



**Technical data**

- Adjustable tilt angle (35°)
- On/Off button
- Timer (1 to 15 min)
- Fan caster for cooling
- Cable rewind
- Carrying handle
- Large-surface application
- Incl. lamp (halogen)
- Certified medical device\*

**Size:** approx. 19.8×18.6×29.5 cm  
**Weight:** approx. 1.5 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 300 wattss  
**Article No.:** 88257  
**EAN:** 40 15588 88257 9  
**Warranty:** 3 Years



**IR 850**  
Infrared lamp

**With illuminated LED button**

- Intense infrared light to relax muscles or treat colds
- Three-stage electronic timer (5, 10 or 15 minutes) with LED indicator
- Cool down function with auto off



**Technical data**

- Ceramic glass to block UV light
- Adjustable tilt angle (-15° to +35°)
- Recommended treatment distance: 40 cm
- Illuminated LED On/Off button
- Hot warning indicator
- Incl. lamp (halogen)
- Certified medical device\*

**Size:** approx. 20×18×28 cm  
**Weight:** approx. 1.2 kg  
**Power supply:** 220V~ 50 Hz  
**Power output:** 300 watts  
**Article No.:** 88303  
**EAN:** 40 15588 88303 3  
**Warranty:** 3 Years



**IR 350**  
Infrared lamp

**With timer function in 3 stages**

- Intense infrared light to relax muscles
- Three-stage electronic timer (5, 10 or 15 minutes) with LED indicator
- Stimulates blood circulation



**Technical data**

- Adjustable angle of inclination
- Durable lamp (Philips)
- Automatic switch-off

**Size:** approx. 26.4×18.3×18.3 cm  
**Weight:** approx. 1.17 kg  
**Power supply:** 220V~ 50Hz, 2.5A  
**Power output:** 150 watts  
**Article No.:** 88226  
**EAN:** 40 15588 88226 5  
**Warranty:** 3 Years



**IR 200**  
Infrared lamp

**Individually adjustable angle of inclination**

- Intense infrared light to relax muscles
- Three-stage electronic timer (5, 10 or 15 minutes) with LED indicator
- Stimulates blood circulation



**Technical data**

- Adjustable angle of inclination
- Durable lamp (Philips)
- Automatic switch-off

**Size:** approx. 23.3×20.2×17.8 cm  
**Weight:** approx. 0.55 kg  
**Power supply:** 220V~ 50Hz, 2.5A  
**Power output:** 100 watts  
**Article No.:** 88220  
**EAN:** 40 15588 88220 3  
**Warranty:** 3 Years



**IR 100**  
Infrared lamp

**Helps to relieve colds**

- Infrared lamp to relax the muscles
- Soothing if you have a cold
- Stimulates blood circulation
- New with safety handle




**Technical data**

- Adjustable angle of inclination
- Durable lamp (Philips)
- Made in EU

**Size:** approx. 16.5×13.5×20.5 cm  
**Weight:** approx. 0.36 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 100 watts  
**Article No.:** 88232  
**EAN:** 40 15588 88232 6  
**Warranty:** 3 Years

\*Certification according to MDD (Medical Device Directive)





**LT 460**  
Daylight lamp

**LT 500**  
Daylight lamp

**As flat as a tablet**

- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Uniform light distribution with LED technology
- 2 light intensities can be selected

**Particularly silent**

- Energizing daylight lamp for increasing the energy balance and avoiding symptoms of winter depression
- Uniform light distribution through LED technology
- Colour-changing wellness light with 4 colours



**Technical data**

- Light intensity: up to 10,000 Lux
- Comfortable one-button operation
- Can be used horizontally and vertically
- Ideal for home and on the road
- LED technology
- Incl. storage bag
- Certified medical device\*

**Technical data**

- Flicker-free
- Long-lasting, energy-saving LED lighting
- Light intensity: up to 10,000 Lux
- 2 light intensities can be selected
- Particularly silent
- Comfortable two-button operation
- Certified medical device\*

**Size:** approx. 26×15×1.7 cm  
**Weight:** approx. 0.5 kg  
**Light area:** 20.5×12 cm  
**Power supply:** Input: 100-240V~ 50/60 Hz, 500mA; Output: 12V-1 A  
**Power output:** 12 watts  
**Cable length:** approx. 1.8 m  
**Article No.:** 45220  
**EAN:** 40 15588 45220 8  
**Warranty:** 3 Years

**Size:** approx. 25×25×8 cm  
**Weight:** approx. 0.6 kg  
**Light area:** ø 25 cm  
**Power supply:** 220-240V~ 50 Hz 0.2A  
**Power output:** 23 watts  
**Cable length:** approx. 1.8 m  
**Article No.:** 45226  
**EAN:** 40 15588 45226 0  
**Warranty:** 3 Years



**IN 600**  
Inhaler



**IN 510**  
Inhaler



**IN 550 PRO**  
Inhaler

**With rechargeable battery**

- High efficiency through innovative micro compressor technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Powered by rechargeable battery and Micro USB: ideal for on the go

**With comprehensive accessories**

- Nebulisation with compressed air technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories

**3 attachments for a more targeted application**

- Nebulisation with compressed air technology
- Including 3 attachments to produce different particle sizes:
  - 3.0 – 4.0 µm for the deeper areas of the lung
  - 3.5 – 4.5 µm for respiratory tract and lungs front
  - 4.5 – 6.0 µm for mouth and throat
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)



**Technical data**

- 4 LED battery charge indicator
- Comprehensive set of accessories
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulisation
- Compact, light-weighted and easy to use
- Automatic switch-off with an empty solution tank
- Incl. string pouch
- Certified medical device\*

**Technical data**

- Durable and modern design
- Easy to use
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Incl. string pouch
- Certified medical device\*

**Technical data**

- Comprehensive set of accessories: Inhalation with mask, mouthpiece or nose adapter – ideal for children
- Highly efficient due to very precise nebulisation
- Particularly effective thanks to high proportion of respirable aerosols
- Extra long hose – 2 m
- Integrated accessory tray
- Certified medical device\*

Nebuliser bottle, air tube, adult mask, pediatric mask, angled mouth piece, nose adapter, 5 replacement filters

Nebulizer bottle, air tube, adult mask, pediatric mask, angled mouth piece, nose adapter, 5 replacement filters

Incl. mouthpiece, masks and nose adapters for children and adults, connecting piece, 3 attachments to produce different particle sizes, 5 replacement filters, 2 replacement fuses

**Size:** approx. 10.8×7.2×4.4 cm  
**Weight:** approx. 0.22 kg  
**Power supply:** 100-240V~ 50/60Hz  
**Power output:** 10 watts  
**Filling volume:** max. 6 ml  
**Nebulizing output:** 0.25 ml/min  
**Particle size:** approx. 2.9 µm  
**Hose length:** approx. 1.8 m  
**Article No.:** 54542  
**EAN:** 40 15588 54542 9  
**Warranty:** 3 Years

**Size:** approx. 18×14×9.4 cm  
**Weight:** approx. 1.45 kg  
**Power supply:** 230V~ 50 Hz  
**Power output:** 100 watts  
**Nebulizing output:** approx. 0.30 ml/min  
**Particle size:** approx. 3 µm  
**Hose length:** approx. 1.8 m  
**Article No.:** 54547  
**EAN:** 40 15588 54547 4  
**Warranty:** 3 Years

**Size:** approx. 33×20×13.8 cm  
**Weight:** approx. 2 kg  
**Power supply:** 230V~ 50 Hz  
**Power output:** 60 watts  
**Filling volume:** max. 10 ml  
**Nebulizing output:** ≤ 0.2 ml/min  
**Particle size:** >60% ≤ 5 µm  
**Hose length:** approx. 2 m  
**Article No.:** 54530  
**EAN:** 40 15588 54530 6  
**Warranty:** 3 Years

\*Certification according to MDD (Medical Device Directive)





**IN 500 COMPACT**  
Inhaler

**Our bestseller**

- Nebulisation with compressed air technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)



**Technical data**

- Comprehensive set of accessories: Inhalation with mask, mouthpiece or nose adapter – ideal for children
- Highly efficient due to very precise nebulisation
- Particularly effective thanks to high proportion of respirable aerosols
- Extra long hose – 2 m
- Certified medical device\*

Incl. mouthpiece, masks and nose adapters for children and adults, connecting piece, 5 replacement filters

**Size:** approx. 19.5×13.5×9.2 cm  
**Weight:** approx. 1.3 kg  
**Power supply:** 230V~ 50 Hz  
**Power output:** 60 watts  
**Filling volume:** max. 10 ml  
**Nebulizing output:** ≤ 0.2 ml/min  
**Particle size:** >60% ≤ 5 µm  
**Hose length:** approx. 2 m  
**Article No.:** 54520  
**EAN:** 40 15588 54520 7  
**Warranty:** 3 Years



**IN 530**  
Inhaler

**Compact and lightweight**

- High efficiency through innovative micro compressor technology
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Powered by rechargeable battery: ideal for on the go
- Comprehensive set of accessories



**Technical data**

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulisation
- Compact, light-weighted and easy to use
- Automatic switch-off with an empty solution tank
- Incl. storage case
- Certified medical device\*

Nebuliser bottle, adult mask, pediatric mask, angled mouth piece, nose adapter, 5 replacement filters

**Size:** approx. 4.5×4.3×11.5 cm  
**Weight:** approx. 0.19 kg  
**Batteries:** Lithium-ion battery  
**Power output:** 10 watts  
**Filling volume:** max. 6 ml  
**Nebulizing output:** 0.25 ml/min  
**Particle size:** approx. 2.9 µm  
**Article No.:** 54540  
**EAN:** 40 15588 54540 5  
**Warranty:** 3 Years



**IN 525**  
Ultrasonic inhaler

**Small inhaler – ideal for on the go**

- High efficiency through innovative micromembrane nebulisation (mesh technology) with ultrafine droplets
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories
- Use for a short time in the supine position is also possible



**Technical data**

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulisation
- Low power consumption (no more than 0.9 watts)
- Automatic switch-off with an empty solution tank
- Incl. storage case
- Certified medical device\*

Incl. mouthpiece and masks for children and adults

**Size:** 7.2×4×10.8 cm  
**Weight:** approx. 0.1 kg (without batt.)  
**Batteries:** 2×1.5V (AA)  
**Filling volume:** max. 8 ml  
**Nebulizing output:** ≥ 0.2 ml/min  
**Particle size:** Average <5 µm  
**Article No.:** 54115  
**EAN:** 40 15588 54115 5  
**Warranty:** 3 Years



**USC**  
Ultrasonic inhaler

**Our classic**

- High efficiency through innovative micro-membrane nebulisation (mesh technology) with ultrafine droplets
- For targeted treatment of diseases of the upper and lower respiratory tract (e.g. asthma and cold)
- Comprehensive set of accessories



**Technical data**

- Particularly effective thanks to high proportion of respirable aerosols
- Inhalation with mask or mouthpiece – also eminently suitable for children
- Low-noise nebulisation
- Incl. storage case
- 4 batteries (1.5V AA)
- Certified medical device\*

Incl. mouthpiece, masks for children and adults and connecting piece

**Size:** approx. 14.2×7.5×7.5 cm  
**Weight:** approx. 0.3 kg (without batt.)  
**Batteries:** 4×1.5V (AA)  
**Filling volume:** max. 8 ml  
**Nebulizing output:** ≥ 0.25 ml/min  
**Particle size:** Average <4.8 µm  
**Article No.:** 54105  
**EAN:** 40 15588 54105 6  
**Warranty:** 3 Years

\*Certification according to MDD (Medical Device Directive)



**AH 670**  
Air humidifier

**With wellness light**

- Air humidifier with ultrasound technology to improve the breathing air quality
- Wellness light
- Integrated aroma dispenser for scented oils
- Smart control by remote control



**Technical data**

- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Two 360° adjustable nozzles
- Water tank 1,800 ml
- Timer function (1/2/4 hrs)
- High efficiency through micro-fine atomisation
- 2 intensity levels to choose from
- Nebulizing output: 200 ml/h
- Incl. water filter
- Operating time: approx. 9 h



**AH 661**  
Air humidifier

**Germ-reducing pre-heating**

- Air humidifier with ultrasound technology to improve the breathing air quality
- Hygienically clean through an integrated pre-heating function to reduce germs and bacteria
- Integrated aroma dispenser for scented oils



**Technical data**

- 3 fogging intensities to choose from
- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Water tank 3,500 ml
- Touch control buttons
- Low-noise and energy-saving
- Timer for programming the operating time: 2/4/6 hours



**AH 662**  
Air humidifier

**Aroma dispenser for scented oils**

- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Integrated aroma dispenser for scented oils



**Technical data**

- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Water level indicator (full=blue, empty=red)
- Water tank 900 ml
- Incl. cleaning brush



**AH 660**  
Air humidifier

**Including water filter**

- Air humidifier with ultrasound technology to improve the breathing air quality
- High efficiency through micro-fine atomisation
- Steplessly adjustable misting levels



**Technical data**

- Removable tank
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving: 30 watts
- Water tank 4,500 ml
- For rooms up to approx. 30 m<sup>2</sup>
- Incl. water filter



**Size:** approx. 24×17×36 cm  
**Weight:** approx. 1.3 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 30 watts  
**Nebulizing output:** max. 300 ml/h  
**Operating time:** approx. 15 Std.  
**Article No.:** 60067  
**EAN:** 40 15588 60067 8  
**Warranty:** 3 Years



**HG 100**  
Hygrometer

**Including magnetic mount**

- Displays humidity and room temperature for a healthy climate indoors
- Saves the highest and lowest values
- For hanging on the wall (magnetic mount) or standing upright



**Technical data**

- Temperature measurement range: -10 °C to +50 °C (+14° to +122 F°)
- Humidity measurement range: 20% - 90%
- Comfort zone indicator (smiley display)
- Time display
- Alarm function
- Incl. battery (1.5V AAA)

**Size:** approx. 8.2×7×2 cm  
**Weight:** approx. 0.07 kg  
**Batteries:** 1×1.5V (AAA)  
**Article No.:** 60079  
**EAN:** 40 15588 60079 1  
**Warranty:** 3 Years

**Size:** approx. 18×18×30.3 cm  
**Weight:** approx. 1.3 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 35 watts  
**Nebulizing output:** max. 200 ml/h  
**Operating time:** approx. 9 Stunden  
**Article No.:** 60054  
**EAN:** 40 15588 60054 8  
**Warranty:** 3 Years

**Size:** approx. 18.5×18.5×28.1 cm  
**Weight:** approx. 1.1 kg  
**Power supply:** AC 220-240V~ 50/60 Hz  
**Power output:** 75 watts  
**Nebulizing output:** cold 250ml/h, warm 350ml/h  
**Operating time:** approx. 10-14 Std.  
**Article No.:** 60052  
**EAN:** 40 15588 60052 4  
**Warranty:** 3 Years

**Size:** approx. 15×15×22 cm  
**Weight:** approx. 0.45 kg  
**Power supply:** 100-240V~ 50/60 Hz  
**Power output:** 30 watts  
**Nebulizing output:** max. 60 ml/h  
**Operating time:** approx. 15 Std.  
**Article No.:** 60077  
**EAN:** 40 15588 60077 7  
**Warranty:** 3 Years





**AD 625**  
Aroma diffuser

**Stylish bamboo diffuser**

- Decorative bamboo aroma diffuser for scented oils
- By using ultrasound technology, fine mist of fragrance is created and is spread throughout your entire room
- Wellness light with a colour change in 6 colours
- Stylish living accessory without visible control elements



**Technical data**

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving

**Size:** approx.  $\varnothing 16.6 \times 7.2$  cm  
**Weight:** approx. 0.68 kg  
**Power supply:** Mains adaptor  
**Water tank:** 100 ml  
**Nebulizing output:** max. 16 ml/h  
**Operating time:** approx. 6 hours  
**Article No.:** 60083  
**EAN:** 40 15588 60083 8  
**Warranty:** 3 Years



**AD 620**  
Aroma diffuser

**Our bestseller**

- Aroma diffuser for scented oils
- By using ultrasound technology, fine mist of fragrance is created and is spread throughout your entire room
- Wellness light with a colour change in 5 colours
- Timer for programming operation: 2 and 4 hours



**Technical data**

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving
- Incl. refill cup

**Size:** approx.  $13.2 \times 13.2 \times 14.5$  cm  
**Weight:** approx. 0.25 kg  
**Power supply:** 100-240V~ 50/60 Hz, DC 24V/500mA; 12W  
**Water tank:** 300 ml  
**Nebulizing output:** max. 30 ml/h  
**Laufzeit:** approx. 10 hours  
**Article No.:** 60082  
**EAN:** 40 15588 60082 1  
**Warranty:** 3 Years



**AD 635**  
Aroma diffuser

**Made of frosted glass**

- Aroma diffuser for scented oils with high-quality housing made of frosted glass
- By using ultrasound technology, fine mist of fragrance is created and is spread throughout your entire room
- Wellness light with a colour-change in 6 colours
- Stylish accessory for the living room without visible switches



**Technical data**

- High efficiency through micro-fine atomisation
- Automatic shut-off when the tank is empty
- Low-noise and energy-saving

**Size:**  $\varnothing 15 \times 17.1$  cm  
**Weight:** approx. 0.6 kg  
**Power supply:** Mains adaptor  
**Water tank:** 100 ml  
**Nebulizing output:** max. 16 ml/h  
**Laufzeit:** approx. 6 hours  
**Article No.:** 60085  
**EAN:** 4015588 60085 2  
**Warranty:** 3 Years





**LC 870**  
Electric lice comb

**With LED light**

- Integrated LED light for an easier detection and removal of lice
- Simple combing frees hair from head lice and nits
- No side effects, no use of chemical agents
- Also suitable for dogs and cats



**Technical data**

- The lice are destroyed by an electrical current upon contact with the teeth of the comb

**Size:** approx. 20.9×6.6×3.1 cm  
**Weight:** approx. 0.1 kg  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 41019  
**EAN:** 40 15588 41019 2  
**Warranty:** 3 Years



**LC 860**  
Electric lice comb

**Hair-gentle treatment**

- Simple combing frees hair from head lice and nits
- No side effects, no use of chemical agents
- Also suitable for dogs and cats



**Technical data**

- The lice are destroyed by an electrical current upon contact with the teeth of the comb

**Size:** approx. 19.3×6.6×2.9 cm  
**Weight:** approx. 0.095 kg  
**Batteries:** 1×1.5V (AA)  
**Article No.:** 41017  
**EAN:** 40 15588 41017 8  
**Warranty:** 3 Years



**NEW**

**SP 100**  
SleepWell Pillow

**Promotes relaxation, sleep and stress reduction**

- Foam pillow with integrated stereo sound and pre-installed relaxation music
- Binaural beats and sleep music
- Includes Sleep Genius Music program that helps you to fall asleep faster, to sleep deeper and longer



**Technical data**

- Large, orthopaedic pillow suitable for all sleep types
- High quality loudspeakers
- Bluetooth® for streaming music
- High density memory foam core
- Mini USB port for charging and uploading MP3 files
- Audio connection for headphones

**Size:** approx. 57.5×38×12.5 cm  
**Weight:** approx. 1.5 kg  
**Batteries:** Lithium-ion battery  
**Power supply:** -5V/110V  
**Power output:** 0.6 watts  
**Article No.:** 61185  
**EAN:** 40 15588 61185 8  
**Warranty:** 3 Years



**NEW**

**PT 100**  
PainShield Ultrasonic Therapy

**Promotes blood circulation**

- Therapeutic, portable pain therapy device for soft tissue and musculoskeletal system
- Improves blood circulation and supports the treatment of pain
- Easy to use anywhere and anytime – at home, at work or even while you sleep (no clinical setting necessary)



**Technical data**

- Low intensity, low frequency ultrasound
- Treatment of large (body) surfaces
- Allows the treatment of inaccessible areas
- No side effects
- No ultrasound gel required

**Size:** approx. 24×16×7 cm  
**Weight:** approx. 0.5 kg  
**Batteries:** Lithium-ion battery  
**Power supply:** 100-240V~ 50/60Hz  
**Power output:** 6 watts  
**Article No.:** 88346  
**EAN:** 40 15588 88346 0  
**Warranty:** 3 Years





**RC 410**  
Shiatsu relax chair

**With cantilever function**

- 2 massage types for the back area: Shiatsu and spot massage
- 3 selectable massage zones: total, upper or lower back
- Additional heat function
- Soft cushioning for sitting comfort



**Technical data**

- High quality design, easy to clean material
- Timer function: 5/10/15 minutes
- High backrest
- Handy remote control
- Loadable up to 150 kg

**Size:** approx. 67×90×107 cm  
**Weight:** approx. 13 kg  
**Power supply:** 100-240V~ 50/60 Hz  
 12V DC, 3000mA  
**Power output:** 36 watts  
**Article No.:** 88410  
**EAN:** 40 15588 88410 8  
**Warranty:** 3 Years



**RS 650**  
Lounge Chair

**Stylish interior accessory**

- High seating comfort
- USB charging socket
- 6 massage types in an automatic sequence program
- 3 intensity levels selectable



**Technical data**

- Kneading-, knocking-, Shiatsu, rolling-, clapping and Swedish massage
- Massage zone settings: full, upper or lower back
- Automatic switch-off after 15 min utes
- Appealing, modern design

**Size:** approx. 105×60×90 cm  
**Weight:** approx. 28 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 100 watts  
**Article No.:** 88414  
**EAN:** 40 15588 88414 6  
**Warranty:** 3 Years



**MCG 820**  
COMFORT Shiatsu acupressure  
massage seat cover

**Soft-Touch technology**

- Gel massage heads provide a soft and natural massage experience
- Soothing Shiatsu massage for the back
- Soft sprung massage heads gently massage the neck area



**Technical data**

- 3 selectable massage zones for Shiatsu massage: total, upper or lower back or point-specific spot massage
- Acupressure massage in the waist area and thighs with 3 selectable intensity levels
- Separate, height-adjustable neck massage
- Additional red light and heat function
- Handy cable remote control
- Removable and washable neck cover
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 124.5×19×43.5 cm  
**Weight:** approx. 7 kg  
**Power supply:** 220-240~50/60Hz  
**Power output:** 40 watts  
**Article No.:** 88991 (silver)  
 88995 (black)  
**EAN:** 40 15588 88991 2 (silver)  
 40 15588 88995 0 (black)  
**Warranty:** 3 Years



**MCG 810**  
COMFORT Shiatsu massage seat cover

**Soft-Touch technology**

- Gel massage heads provide a soft and natural massage experience
- Soothing Shiatsu massage for the back
- Soft sprung massage heads gently massage the neck area



**Technical data**

- 3 selectable massage zones for Shiatsu massage: total, upper or lower back or point-specific spot massage
- Separate, height-adjustable neck massage
- Separate vibration massage in the seat with 3 selectable intensity levels
- Additional red light and heat function
- Handy cable remote control
- Removable and washable neck cover
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 124.5×19×43.5 cm  
**Weight:** approx. 5.5 kg  
**Power supply:** 220-240~50/60Hz  
**Power output:** 48 watts  
**Article No.:** 88990  
**EAN:** 40 15588 88990 5  
**Warranty:** 3 Years

**Massage seat covers**



**MC 828**  
Shiatsu massage seat cover

**Heating and cooling function**

- Targeted, intensive massage for tangible sense of well-being
- Air stream through the mesh of the back area provides a pleasant heating or cooling function
- Intensive Shiatsu and tapping massage in the back area



**Technical data**

- 3 massage types for the back area: Shiatsu, tapping and spot massage
- 3 intensity levels of tapping massage selectable
- Separate, height-adjustable neck massage with two directions
- 5 automatic massage programs selectable
- Separate vibration massage in the seat
- Automatic power off after 15 minutes
- Handy cable remote control
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 124×18.3×44 cm  
**Weight:** approx. 8.8 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 140 watts  
**Article No.:** 88928  
**EAN:** EAN 40 15588 88928 8  
**Warranty:** 3 Years



**MC 826**  
Shiatsu massage seat cover

**Holistic relaxation with tapping and acupressure massage**

- Targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist area with 3 selectable intensity levels
- 2 massage programs for the back selectable (shiatsu and tapping massage)



**Technical data**

- Relaxing massage for the entire upper or lower back
- Neck and shoulder massage (height-adjustable)
- Separate vibration massage in the seat with 3 selectable intensity levels
- Targeted spot massage for the tapping massage
- For tapping massage 3 selectable programs as well as 3 intensity levels selectable
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 120×26×46 cm  
**Weight:** approx. 8.5 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** approx. 60 watts  
**Article No.:** 88926  
**EAN:** 40 15588 88926 4  
**Warranty:** 3 Years



**MC 825 PLUS**  
Shiatsu acupressure massage seat cover

**With spot massage**

- Targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist and thigh area
- Separate, height-adjustable neck massage



**Technical data**

- 3 selectable massage zones: total, upper or lower back
- 3 intensity levels of acupressure massage selectable
- Intensive shiatsu and spot massage in the back area
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 110×42.5×18 cm  
**Weight:** approx. 5.9 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 40 watts  
**Article No.:** 88956  
**EAN:** 40 15588 88956 1  
**Warranty:** 3 Years



**MC 825**  
Shiatsu acupressure massage seat cover

**Acupressure for waist and thigh area**

- Targeted, intensive massage for tangible sense of well-being
- Acupressure massage in the waist and thigh area
- Separate, height-adjustable neck massage



**Technical data**

- 3 selectable massage zones: total, upper or lower back
- 3 intensity levels of acupressure massage selectable
- Intensive shiatsu massage in the back area
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 110×4.5×18 cm  
**Weight:** approx. 5.9 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 40 watts  
**Article No.:** 88939  
**EAN:** 40 15588 88939 4  
**Warranty:** 3 Years



**MC 824**  
Shiatsu massage seat cover

**With body-scan function**

- Pleasant massage along the spine with new rocking technique
- 3 massage types for the back area: Shiatsu, rolling and spot massage
- Body-scan function detects when shoulders are reached and automatically reverses back down
- Shiatsu and rolling massage in 2 speeds



**Technical data**

- 3 massage zones for shiatsu massage selectable: total, lower or upper back or shoulder area
- Separate vibration massage in the seat
- Additional heat function
- Handy cable remote control
- Auto-stop after 15 minutes
- Removable and washable protective cover and neck cover
- Easy to clean material
- Convenient velcro for fastening at the back of the chair
- 2 Demo-programs

**Size:** approx. 135×48×21 cm  
**Weight:** approx. 7.8 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** approx. 60 watts  
**Article No.:** 88921  
**EAN:** 40 15588 88921 9  
**Warranty:** 3 Years



**MC 822**  
Shiatsu massage seat cover

**Electrical and height-adjustable neck massage**

- Targeted, intensive massage for tangible sense of well-being
- 3 types of massage in the back area: Shiatsu, rolling and spot massage



**Technical data**

- 3 selectable massage zones: total, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 118×40×9 cm  
**Weight:** approx. 6.5 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 40 watts  
**Article No.:** 88922  
**EAN:** 40 15588 88922 6  
**Warranty:** 3 Years



**MC 815**  
Shiatsu massage seat cover

**With 4 rotating massage heads**

- Targeted, intensive massage for tangible sense of well-being
- 3 selectable massage zones: total, upper or lower back



**Technical data**

- 4 rotating massage heads
- Separate vibration massage in the seat with 3 selectable intensity levels
- Additional red light and heat function
- Timer function: 5/10/15 minutes
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair
- Demo function

**Size:** approx. 105×45×9 cm  
**Weight:** approx. 3.4 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 24 watts  
**Article No.:** 88916  
**EAN:** 40 15588 88916 5  
**Warranty:** 3 Years





**MC 830**  
Shiatsu massage seat cover  
with gel heads

**Soft back massage**

- Targeted, intensive massage for tangible sense of well-being
- Massage heads made of high-quality gel for a softer and more natural massage
- More gentle to hips, spine and shoulder blades
- 3 selectable massage zones: total, upper or lower back



**Technical data**

- Separate vibration massage in the seat with 3 selectable intensity levels
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 105×43×8.5 cm  
**Weight:** approx. 3.4 kg  
**Power supply:** 100-240V~ 50/60 Hz  
**Power output:** 18 watts  
**Article No.:** 88944  
**EAN:** 40 15588 88944 8  
**Warranty:** 3 Years



**MCN PRO**  
Shiatsu massage seat cover

**Our bestseller with height-adjustable neck massage**

- Targeted, intensive massage for tangible sense of well-being
- Separate adjustable neck massage
- Separate vibration massage in the seat with 3 selectable intensity levels



**Technical data**

- 3 selectable massage zones: total, upper or lower back
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 125×43×36 cm  
**Weight:** approx. 6 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 48 watts  
**Article No.:** 88970  
**EAN:** 40 15588 88970 7  
**Warranty:** 3 Years



**MCN**  
Shiatsu massage seat cover

**Our bestseller**

- Targeted, intensive massage for tangible sense of well-being
- Separate neck massage
- Separate vibration massage in the seat with 3 selectable intensity levels



**Technical data**

- 3 selectable massage zones: total, upper or lower back
- Additional red light and heat function
- Auto-stop after 15 minutes
- Removable and washable neck cover
- Handy cable remote control
- Easy to clean material
- Convenient velcro for fastening at the back of the chair

**Size:** approx. 118×40×9 cm  
**Weight:** approx. 5.5 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 48 watts  
**Article No.:** 88930  
**EAN:** 40 15588 88930 1  
**Warranty:** 3 Years



**MC 810**  
Vibration massage seat cover

**With 9 selectable massage programs**

- 4 massage zones in the back and seat area
- 3 selectable massage intensities
- 9 selectable massage programs
- Intensive heat function with or without massage



**Technical data**

- 6 vibration motors
- Auto-stop after 15 minutes
- Handy cable remote control
- Easy to clean material
- Very handy and easy to store
- Also ideal for office use or in the car
- Incl. car power adapter

**Size:** approx. 107×46×4 cm  
**Weight:** approx. 1.4 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12 V DC 2500mA  
**Power output:** 30 watts  
**Article No.:** 88937  
**EAN:** 40 15588 88937 0  
**Warranty:** 3 Years

**Massage seat covers**



**MCH**  
Vibration massage seat cover

**With additional heat function**


- 4 massage zones in the back and seat area
- 3 selectable massage intensities
- 5 selectable massage programs
- Additional heat function



**Technical data**

- 5 vibration motors
- Auto-stop after 15 minutes
- Handy cable remote control
- Easy to clean material
- Very handy and easy to store
- Also ideal for office use or in the car
- Incl. car power adapter

**Size:** approx. 107×47×3.5 cm  
**Weight:** approx. 1.1 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 800mA  
**Power output:** 9.6 watts  
**Article No.:** 88935  
**EAN:** 40 15588 88935 6  
**Warranty:** 3 Years



NEW



**NM 890**  
Shiatsu neck massager

**Like a real kneading massage**

- Intensive shiatsu massage for the neck
- Specifically adopted to the ergonomics of the neck area
- 2 massage modes selectable clockwise and anti-clockwise rotation
- 3 levels of rotation speed selectable



**Technical data**

- Intensity can be individually adjustable by pulling the holding straps
- Also applicable on shoulder, back, abdomen, upper and lower leg
- Easy wireless operation through integrated control panel
- Additional heat function
- Removable and washable neck cover

**Size:** approx. 38×18.5×29.5 cm  
**Weight:** approx. 2.1 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 2.500 mA  
**Power output:** 30 watts  
**Article No.:** 88975  
**EAN:** 40 15588 88975 2  
**Warranty:** 3 Years

NEW



**NMG 850**  
COMFORT Shiatsu neck massager

**Soft-Touch technology**

- Gel massage heads provide a soft and natural massage experience
- Relaxing Shiatsu massage for the neck and shoulders
- 3 massage speed levels selectable



**Technical data**

- Massage intensity can be individually adjusted by pulling on the holding straps
- Multi functional usability - also applicable on shoulders, back, abdomen, upper and lower legs
- Easy wireless operation through integrated control panel
- Easy to clean material
- Additional red light and heat function
- Removable and washable neck cover

**Size:** approx. 39.5×18×11 cm  
**Weight:** approx. 1.69 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 2.000 mA  
**Power output:** 24 watts  
**Article No.:** 88993 (silver), 88997 (black)  
**EAN:** 40 15588 88993 6 (silver), 40 15588 88997 4 (black)  
**Warranty:** 3 Years



**NM 880**  
Shiatsu neck massager

**With rechargeable battery**

- Relaxing shiatsu massage for the neck and shoulders
- Wireless operation through rechargeable battery
- Additional red light and heat function



**Technical data**

- Ergonomic fit
- Easy care polyester cover and a pleasant foam padding
- Integrated control panel
- Intensity can be individually adjusted by pulling the holding straps
- Also applicable on shoulder, back, abdomen, upper and lower leg

**Size:** approx. 38×33.5×17 cm  
**Weight:** approx. 1.6 kg  
**Power supply:** Li-ion battery: 11.1V 2.200 mAh  
**Power output:** 18 watts  
**Article No.:** 88948  
**EAN:** 40 15588 88948 6  
**Warranty:** 3 Years



**Neck massager**



**NM 860 SILVER**  
Shiatsu neck massager

**For an optimal neck massage**

- Intensive shiatsu massage in the neck area
- Also applicable on shoulder, back, abdomen, upper and lower leg
- 4 rotating massage heads
- Massage intensity individually adjustable through pulling the holding straps



**Technical data**

- 2 levels of rotation speed selectable
- Additional light and heat function in the upper massage heads
- Easy wireless operation through integrated control panel
- Removable and washable neck cover

**Size:** approx. 35.5×28.5×20 cm  
**Weight:** approx. 1.55 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 1.500 mA  
**Power output:** 18 watts  
**Article No.:** 88961  
**EAN:** 40 15588 88961 5  
**Warranty:** 3 Years



**NM 866**  
Vibration neck massager

**With intensive heat function**

- 2-levelled vibration massage for the neck and shoulders
- Ideal for home and on the go
- Optional intensive heat function (with AC power supply)



**Technical data**

- Functional design
- Intensity can be individually adjusted by pulling the holding straps
- Optional mains or battery operation
- Easy wireless operation through integrated control panel

**Size:** approx. 38×33.5×17 cm  
**Weight:** approx. 1.6 kg  
**Power supply:** 100-240V~ 50/60 Hz 6V DC 800 mA  
**Batteries:** 4×1.5V (AA)  
**Power output:** 24 watts  
**Article No.:** 88972  
**EAN:** 40 15588 88972 1  
**Warranty:** 3 Years





NEW

**NM 870**  
**Shiatsu neck massager**

**Relaxed travelling companion**

- Pleasant vibration massage for the neck and shoulders
- Ideal for travelling and at home
- Extra soft and fluffy fleece cover

**MCG 800**  
**COMFORT Shiatsu massage cushion**

**Soft-Touch technology**

- Gel massage heads provide a soft and natural massage experience
- Soothing Shiatsu massage for neck, shoulder, back and legs
- Multi functional usability

**MC 840**  
**Shiatsu massage cushion**

**Particularly flat design**

- Soothing shiatsu massage for neck, shoulder, back and legs
- 4 rotating massage heads
- Additional red light and heat function
- Particularly flat design for flexible, convenient operation



**Technical data**

- Neck massager made of comfortable molded foam – fits perfectly to the contours of the head and neck
- With integrated „One Touch“- Control panel
- With integrated storage bag – space-saving storage or on the go

**Technical data**

- 4 rotating gel massage heads with selectable rotating direction
- Additional red light and heat function
- Particularly flat design for flexible, convenient operation
- Easy to clean material
- Convenient velcro for fastening at the back of the chair
- Handy cable remote control

**Technical data**

- Handy cable remote control
- Convenient velcro for fastening at the back of the chair
- Easy to clean material

**Size:** approx. 28 × 26 × 9 cm  
**Weight:** approx. 0.425 kg  
**Batteries:** 2 x 1.5V (AA)  
**Power output:** 3.6 watts  
**Article No.:** 88946  
**EAN:** 40 15588 88946 2  
**Warranty:** 3 Years

**Size:** approx. 34.5×22×125 cm  
**Weight:** approx. 1.8 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 1500 mA  
**Power output:** 18 watts  
**Article No.:** 88992 (silver), 88996 (black)  
**EAN:** 40 15588 88992 9 (silver), 40 15588 88996 7 (black)  
**Warranty:** 3 Years

**Size:** approx. 31×23×11 cm  
**Weight:** approx. 1.3 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 1000 mA  
**Power output:** 12 watts  
**Article No.:** 88949  
**EAN:** 40 15588 88949 3  
**Warranty:** 3 Years



**HM 886**  
**Hand massager**

**Including 5 massage attachments**

- Red light and heat function
- 2 selectable intensity levels
- Different massage angles adjustable
- Including 5 massage attachments

**HM 858**  
**Hand massager**


**Vibration, acupressure and spot massage**

- Pleasant vibration massage
- Continuously adjustment of massage intensity
- Additional red light and heat function
- Various massage attachments: acupressure and spot massage

**ITM**  
**Hand massager**

**Tapping, acupressure and spot massage**

- Vigorous tapping massage
- Continuously adjustment of massage intensity
- Additional red light and heat function
- Various massage attachments: acupressure and spot massage



**Technical data**




- Handy design for a convenient self-massage
- Vibration massage

**Technical data**

- Easy to use thanks to its ergonomic shape
- Non-slip grip
- Extra-long cable

**Technical data**

- Easy to use thanks to its ergonomic shape
- Non-slip grip
- Extra-long cable



**Size:** approx. 41.1×10.2×12.2 cm  
**Weight:** approx. 0.3 kg  
**Power supply:** 220-240V~ 50/60 Hz, 6.5 watts  
**Article No.:** 88268  
**EAN:** 40 15588 88268 5  
**Warranty:** 3 Years

**Size:** approx. 32.5×17×14.3 cm  
**Weight:** approx. 1.8 kg  
**Power supply:** 220-240V~ 50/60 Hz, 28 watts  
**Article No.:** 88292  
**EAN:** 40 15588 88292 0  
**Warranty:** 3 Years

**Size:** approx. 44×14×13 cm  
**Weight:** approx. 1.7 kg  
**Power supply:** 220-240V~ 50 Hz, 35 watts  
**Article No.:** 88275  
**EAN:** 40 15588 88955 4  
**Warranty:** 3 Years



**MM 825**  
Massage mat

**Massage for the entire body**

- 4 selectable massage zones: upper and lower back, upper and lower leg incl. 2 intensity levels for each massage zone
- Comfortable vibration massage for the whole body
- Soothing heat function in the back area



**Technical data**

- 2 intensity levels for each massage zone
- 5 selectable massage programs
- Additional heat function
- Auto-stop after 15 minutes
- Handy cable remote control
- Extra soft and cosy fleece cover
- Cushion in the head area of the mat for lying comfortably

**Size:** approx. 180 × 62 × 10 cm  
**Weight:** approx. 1.8 kg  
**Power supply:** 100-240V~ 50/60 Hz, 12V DC 1200mA  
**Power output:** 14.4 watts  
**Article No.:** 88955  
**EAN:** 40 15588 88955 4  
**Warranty:** 3 Years



**NEW**

**FMG 880**  
COMFORT Shiatsu foot massager

**Soft-Touch technology**

- Gel massage heads provide a soft and natural massage experience
- Deep kneading Shiatsu massage of the foot and heels
- Upper foot compression massage with 3 selectable intensity levels



**Technical data**

- Promotes blood circulation, relaxes tense muscles and invigorates tired feet
- Additional red light and heat function
- Comfortable to use with the foot
- Easy to clean material

**Size:** approx. 34.5 × 32 × 15 cm  
**Weight:** approx. 2.4 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 30 watts  
**Article No.:** 88994  
**EAN:** 40 15588 88994 3  
**Warranty:** 3 Years



**FM 888**  
Shiatsu foot massager

**Intensive foot reflexology massage**

- 2 different massage programs selectable
- Shiatsu massage of the foot soles and heels
- Upper foot compression massage
- Comfortable size – suitable up to shoe size 11 (UK)



**Technical data**

- Additional red light and heat function
- 2 intensity levels for the compression massage selectable
- Promotes the blood circulation and invigorates tired feet
- Removable and washable inner lining
- Easy to clean material

**Size:** approx. 41.7 × 35.7 × 22.4 cm  
**Weight:** approx. 4.1 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 30 watts  
**Article No.:** 88398  
**EAN:** 40 15588 88398 9  
**Warranty:** 3 Years



**FM 885**  
Shiatsu foot massager

**Shiatsu massage of the foot soles and heels**

- Compression massage of the foot topsides in 3 intensity levels
- Switchable heat in the shiatsu massage heads



**Technical data**

- Space-saving softcover
- Improves blood circulation and acts invigorating
- Comfortable to use with the foot

**Size:** approx. 40.2 × 37.5 × 17.7 cm  
**Weight:** approx. 3.5 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 28 watts  
**Article No.:** 88397  
**EAN:** 40 15588 88397 2  
**Warranty:** 3 Years



**FM 883**  
Shiatsu foot massager

**Additional red light and heat function**

- Soothing shiatsu massage
- Suitable for feet and back\*
- 2 selectable massage intensities



**Technical data**

- 6 × 3 rotating massage heads
- Two height adjustable feet to adjust the massage position
- Handy cable rewind under the device
- Simple one-button operation

**Size:** approx. 40.5 × 36.5 × 12 cm  
**Weight:** approx. 3.3 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 30 watts  
**Article No.:** 88396  
**EAN:** 40 15588 88396 5  
**Warranty:** 3 Years

\* Device is attachable on chair backrests with adjustable feet





**FS 885 COMFORT**  
Foot spa

**5 in 1 comfort functions**

- Relaxing spa and vibration massage
- 5in1: Bubble- and vibrating massage, red light and heat function and massage rollers in the middle of the foot bubble bath
- Water can be heated to the preferred temperature (20 °C to 48 °C)
- Timer function for 20/30/40/50/60 min.



**Technical data**

- Massage rollers in the middle of the foot bubble bath to massage the inner side of the foot
- Digital display for timer setting and water temperature
- Red light and heat function
- Ergonomically arranged air bubble nozzles
- Removable massage rollers in the area of the soles
- Good heat insulation through double-walled casing
- Intensive foot reflexology massage
- Comfortable size - suitable up to shoe size 11 (UK)



**Size:** approx. 38.6×24.3×45.9 cm  
**Weight:** approx. 2.9 kg  
**Power supply:** 220-240V~ 50 Hz  
**Power output:** 390 watts  
**Article No.:** 88378  
**EAN:** 40 15588 88378 1  
**Warranty:** 3 Years



**FS 883**  
Foot spa

**With pleasant red light**

- Relaxing spa and vibration massage
- Pleasant red light area
- 3in1: Water heating, bubble- and vibrating massage
- Function for keeping the water warm
- 3 exchangeable pedicure attachments (callus remover, brush, massage attachment)



**Technical data**

- Massage rollers for foot reflexology massage
- Massage-aiding footbed
- Dry massage function
- Non-slip rubber feet
- Removable splash protection
- Comfortable size - suitable up to shoe size 11 (UK)



3 exchangeable pedicure attachments

**Size:** approx. 41×35×18 cm  
**Weight:** approx. 1.5 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 80 watts  
**Article No.:** 88387  
**EAN:** 40 15588 88387 3  
**Warranty:** 3 Years



**FS 881**  
Foot spa

**Intensive foot reflexology massage**

- Relaxing spa and vibration massage
- Soothing foot spa
- 3in1: Water heating, bubble- and vibrating massage



**Technical data**

- Non-slip rubber feet
- Automatic heat function
- Spout-opening
- Comfortable size - suitable up to shoe size 9.5 (UK)

**Size:** approx. 39×33×13 cm  
**Weight:** approx. 1.22 kg  
**Power supply:** 220-240V~ 50 Hz  
**Power output:** 60 watts  
**Article No.:** 88380  
**EAN:** 40 15588 88380 4  
**Warranty:** 3 Years



**MBH**  
Bath spa

**Relaxes and loosens the muscles**

- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels and combination program
- Auto-stop after 20 minutes



**Technical data**

- Automatic heat function
- Intensity levels adjustable on the device or via remote control
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated aroma dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma pads
- Incl. remote control

**Size:** Soft mat: approx. 120×36 cm  
 Air hose: approx. 2.4 m  
 Motor: approx. 25×16×17 cm  
 approx. 2.7 kg (in total)  
**Weight:** approx. 2.7 kg (in total)  
**Power supply:** 220-240V~ 50 Hz  
 Remote control: 3 V CR2032  
**Power output:** 450 watts  
**Article No.:** 88375  
**EAN:** 40 15588 88375 0  
**Warranty:** 3 Years



**BBS**  
Bath spa

**Extra quiet**

- Relaxing air bubble bath to stimulate blood circulation and to relieve muscles
- Massage/bubble function with 3 intensity levels and combination program
- Intensive heat function selectable


















**Technical data**

- Automatic heat function
- Intensity levels adjustable on the device or via remote control
- Timer function: 10/20/30 minutes
- Non-slip soft mat for lying comfort
- Suitable for any bath tub
- Integrated aroma dispenser for scented oils
- Double safety valve for air hose and device
- Incl. 3 aroma pads
- Incl. remote control

**Size:** Soft mat: approx. 120×36 cm  
 Air hose: approx. 2.4 m  
 Motor: approx. 31×21×31 cm  
 approx. 4.3 kg (in total)  
**Weight:** approx. 4.3 kg (in total)  
**Power supply:** 220-240V~ 50 Hz  
 Remote control: 3 V CR2032  
**Power output:** 570 watts  
**Article No.:** 88386  
**EAN:** 40 15588 88386 6  
**Warranty:** 3 Years



 <p><b>NEW</b></p>	 <p><b>NEW</b></p>	 <p><b>NEW</b></p>	 <p><b>NEW</b></p>	
<p><b>HP 650</b> Heat pad XL</p>	<p><b>HP 605</b> Heat pad</p>	<p><b>HP 630</b> Heat pad for the shoulders and back</p>	<p><b>HP 622</b> Heat pad for the shoulders</p>	<p><b>HP 626</b> Heat cape for neck and shoulder</p>
<p><b>New XL Size</b></p> <ul style="list-style-type: none"> <li>• Heat pad in comfortable size (45×35 cm)</li> <li>• 4 temperature settings</li> <li>• Automatic shut-off after 90 minutes</li> <li>• Washable thanks to detachable controller</li> <li>• Quick heating with 100 watts</li> </ul>	<p><b>With removable cover</b></p> <ul style="list-style-type: none"> <li>• Removable cover of pleasantly soft upper material</li> <li>• Hand washable thanks to detachable controller</li> <li>• Quick heating with 100 watts</li> <li>• Detachable cover and cable</li> <li>• 4 temperature settings</li> </ul>	<p><b>Fast and comfortable warmth</b></p> <ul style="list-style-type: none"> <li>• Ultra-soft breathable cushion</li> <li>• Innovative 4D heat flow technology – for an effective heat development</li> <li>• Flexible heating wires and a matrix sensor system for even heat distribution</li> <li>• 4 temperature settings</li> </ul>	<p><b>With 6 temperature settings</b></p> <ul style="list-style-type: none"> <li>• Flexible heating wires and a matrix sensor system for even heat distribution</li> <li>• 6 temperature settings</li> <li>• Ergonomic shape for pleasant warmth at shoulder and neck</li> </ul>	<p><b>With rechargeable battery</b></p> <ul style="list-style-type: none"> <li>• For indoor and outdoor use</li> <li>• 4 temperature levels</li> <li>• Battery running time: approx. 5 hours</li> <li>• Individually adjustable fit</li> </ul>
				
<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Particularly soft, pleasant material</li> <li>• Safety system with overheating protection</li> <li>• Quick heating within 2 minutes for low level and appr. 4 minutes for highest level.</li> <li>• Oeko-Tex® Standard 100</li> <li>• Handy controller</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Soft upper material</li> <li>• Safety system with overheating protection</li> <li>• Automatic shut-off after 90 minutes</li> <li>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</li> <li>• Cover machine-washable at 30 °C</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• APS technology – Active protection system: visual and acoustic temperature control</li> <li>• Turbo heat with powerful 100 watts</li> <li>• Automatic shut-off after 90 minutes</li> <li>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</li> <li>• Machine washable at 40 °C</li> <li>• Flexible shape allows individual adjustment to any shoulder size</li> <li>• Fluffy soft uppermaterial</li> <li>• Oeko-Tex® Standard 100</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Turbo heat with powerful 100 watts</li> <li>• Automatic shut-off after 90 minutes</li> <li>• Handy controller</li> <li>• Machine washable at 40 °C</li> <li>• Oeko-Tex® Standard 100</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Fast direct warming</li> <li>• Charging power supply for a fast recharging (approx. 2-3 hours)</li> <li>• Automatic switch-off after 90 minutes</li> <li>• Overheating protection</li> </ul>
				
<p><b>Size:</b> approx. 45×35 cm <b>Weight:</b> approx. 0.33 kg <b>Power supply:</b> 230V~ 50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 61150 <b>EAN:</b> 40 15588 61150 6 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 40×30 cm <b>Weight:</b> approx. 0.35 kg <b>Power supply:</b> 230V~ 50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 61148 <b>EAN:</b> 40 15588 61148 3 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 55×65 cm <b>Weight:</b> approx. 0.72 kg <b>Power supply:</b> 230V~ 50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 61158 <b>EAN:</b> 40 15588 61158 2 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 56×52 cm <b>Weight:</b> approx. 0.5 kg <b>Power supply:</b> 220-240V~ 50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 61156 <b>EAN:</b> 40 15588 61156 8 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 49×42 cm <b>Weight:</b> approx. 0.2 kg <b>Power supply:</b> 8.4 V DC 1.2 A <b>Article No.:</b> 61165 <b>EAN:</b> 40 15588 61165 0 <b>Warranty:</b> 3 Years</p>





NEW

**HB 674**  
Mobile heating blanket 3in1

**Mobile Warmth**

- Cosy, versatile 3in1 heating blanket: Poncho, throw, blanket
- Reversible blanket with 2 design
- Flexible heat „on the Go“
- Pleasantly soft upper material with great cosyness inside
- Even without heat a perfect companion
- 2 Temperature settings



**Technical data**

- Comfortable size: 162×62 cm
- Handy controller
- Heat on the go due to powerbank
- Quick heating in 3 minutes
- Safety system with overheating protection
- Machine washable
- Automatic switch-off after 30 minutes

**Size:** approx. 162×62 cm  
**Weight:** approx. 0.85 kg  
**Powerbank:** 5.300mAh,5V 2A  
**Power output:** 7 watts  
**Article No.:** 61153  
**EAN:** 40 15588 61153 7  
**Warranty:** 3 Years



NEW

**HB 677**  
3in1 heating blanket

**3in1: Heating blanket, poncho and throw**

- Cosy, versatile 3in1 heating blanket
- Can be used as blanket, poncho or throw
- Reversible blanket in 2-tone design
- Even when switched off, it's a perfect companion for cold days
- Soft and pleasant material
- 3 temperature settings



**Technical data**

- Quick heating with 120 watts
- Easy-care and washable thanks to detachable heating element
- Includes safety system: Switches off auto-matically after approx. 180 minutes to prevent overheating
- Oeko-Tex® Standard 100

**Size:** approx. 160×130 cm  
**Weight:** approx. 1.56 kg  
**Power supply:** 220-240 V ~ 50 Hz  
**Power output:** 120 watts  
**Article No.:** 61170  
**EAN:** 40 15588 61170 4  
**Warranty:** 3 Years



NEW

**HB 680**  
Knitted heating blanket

**Heating blanket in knitted design**

- Reversible blanket with 2 designs
- Quick heating with 120 watts
- Automatic switch-off after 180 minutes
- Pleasantly soft material
- 3 temperature settings



**Technical data**

- Cosy due to Microplush side
- Handy controller
- Safety system with overheating protection
- Washable thanks to detachable controller

**Size:** approx. 160×120 cm  
**Weight:** approx. 2.2 kg  
**Power supply:** 220-240 V ~ 50 Hz  
**Power output:** 120 watts  
**Article No.:** 60233  
**EAN:** 40 15588 60233 7  
**Warranty:** 3 Years



NEW

**HB 675**  
Cosy heating blanket XXL

**XXL version**

- Cosy heating blanket with 4 temperature settings
- Quick heating with 120 watts
- Extra big: 200×150 cm
- Automatic shut-off after 180 minutes



**Technical data**

- Safety system with overheating protection
- Handy controller
- Reversible blanket in 2-tone design
- Fluffy soft uppermaterial
- Hand washable thanks to detachable controller
- Material certified to Oeko-Tex® Standard 100

**Size:** approx. 200×150 cm  
**Weight:** approx. 2.2 kg  
**Power supply:** 220-240V ~ 50 Hz  
**Power output:** 120 watts  
**Article No.:** 60231  
**EAN:** 40 15588 60231 3  
**Warranty:** 3 Years







NEW

**HDW**  
Heating blanket

**Extra soft**

- Cosy heating blanket with 4 temperatures
- Quick heating with 120 watts
- Safety system with overheating protection
- Automatic shut-off after 180 minutes



**Technical data**

- Handy controller
- Reversible blanket in 2-tone design
- Fluffy soft upper material
- Material certified to Oeko-Tex® Standard 100
- Hand washable thanks to detachable controller



**Size:** approx. 180 × 130 cm  
**Weight:** approx. 1.7 kg  
**Power supply:** 220-240V~ 50 Hz  
**Power output:** 120 watts  
**Article No.:** 60228  
**EAN:** 40 15588 60228 3  
**Warranty:** 3 Years



NEW

**FW 100**  
Massage foot warmer

**With massage function**

- For tired feet and for cold winter days
- Super comforting massage
- Instantly smoothing heat
- 2 temperature settings
- Large foot space (suitable up to shoe size 46)



**Technical data**

- Removable fluffy soft lining – Machine washable at 30 °C
- Super comforting massage
- Quick heating
- Safety system with overheating protection
- Oeko-Tex® Standard 100



**Size:** approx. 30 × 30 × 24 cm  
**Weight:** approx. 0.68 kg  
**Power supply:** 100-240V  
**Power output:** 7.2 watts  
**Article No.:** 60260  
**EAN:** 40 15588 60260 3  
**Warranty:** 3 Years



NEW

**FWS**  
Foot warmer

**Large foot space**

- Quick heating with 100 watts
- Electronic foot warmer with 3 temperature settings
- Removable fluffy soft lining
- Large foot space (suitable up to shoe size 46)



**Technical data**

- Safety system with overheating protection
- Extra-soft breathable core
- Automatic shut-off after 90 minutes
- Ergonomically designed controller with LED controll display
- Interior Plush fleece removable and washable by 30 °C
- Oeko-Tex® Standard 100



**Size:** approx. 30 × 30 × 23 cm  
**Weight:** approx. 0.6 kg  
**Power supply:** 220-240V~ 50/60 Hz  
**Power output:** 100 watts  
**Article No.:** 60258  
**EAN:** 40 15588 60258 0  
**Warranty:** 3 Years



NEW

**FW 150**  
Knitted foot warmer

**Foot warmer in knitted design**

- For tired feet and for cold winter days
- 3 temperature settings
- Automatic switch off after 3 hours
- Large foot space (suitable up to shoe size 46)



**Technical data**

- Instantly smoothing heat
- Removable and washable inlay
- Safety system with overheating protection
- Nice design in knitted surface
- Interior Plush fleece removable and washable by 30 °C

**Size:** approx. 30 × 30 × 24 cm  
**Weight:** approx. 0.6 kg  
**Power supply:** 220 - 240 V  
**Power output:** 100 watts  
**Article No.:** 60259  
**EAN:** 40 15588 60259 7  
**Warranty:** 3 Years








<p><b>HU 676</b> Kingsize heated underblanket</p> <p><b>Kingsize</b></p> <ul style="list-style-type: none"> <li>• Extra Large Size – comfy for two people</li> <li>• 2 controller for individual comfort</li> <li>• Safety system with overheating protection</li> <li>• Washable thanks to detachable controller</li> <li>• 4 temperature settings</li> </ul>	<p><b>HU 672</b> Maxi fleece heated underblanket</p> <p><b>Sleep like a baby</b></p> <ul style="list-style-type: none"> <li>• Safety system with overheating protection</li> <li>• Heated underblanket with super soft fleece material</li> <li>• Extra warm foot area</li> <li>• 4 temperature settings</li> </ul>	<p><b>HU 674</b> Heated underblanket</p> <p><b>With automatic shut-off</b></p> <ul style="list-style-type: none"> <li>• Safety system with overheating protection</li> <li>• Washable thanks to detachable controller</li> <li>• Pleasant soft upper material</li> <li>• 4 temperature settings</li> </ul>
--	---	--

<p><b>HU 670</b> Fleece heated underblanket</p> <p><b>Soft fleece material</b></p> <ul style="list-style-type: none"> <li>• Heated underblanket with super soft fleece material</li> <li>• 4 temperature settings</li> <li>• Safety system with overheating protection</li> <li>• Extra warm foot area</li> </ul>	<p><b>HU 665</b> Heated underblanket</p> <p><b>Our bestseller</b></p> <ul style="list-style-type: none"> <li>• Safety system with overheating protection</li> <li>• 3 temperature levels</li> <li>• Washable thanks to detachable controller</li> </ul>	<p><b>HU 662</b> Heated underblanket</p> <p><b>Extra warm foot zone</b></p> <ul style="list-style-type: none"> <li>• Safety system with overheating protection</li> <li>• 6 temperature levels</li> <li>• Quick heating with 100 watts</li> </ul>
---	---	---



<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Pleasantly soft upper material</li> <li>• Handy controller</li> <li>• Automatic shut-off after 180 minutes</li> <li>• Oeko-Tex® Standard 100</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Automatic shut-off after 120 minutes</li> <li>• Approx. 2 cm thickness for cozyness</li> <li>• APS technology – Active protection system: visual and acoustic temperature control</li> <li>• Machine washable at 30 °C</li> <li>• 2 side detachable controller</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Handy controller</li> <li>• Automatic shut-off after 180 minutes</li> <li>• Oeko-Tex® Standard 100</li> </ul>
---	---	---

<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Automatic shut-off after 120 minutes</li> <li>• Approx. 1 cm thickness for cozyness</li> <li>• APS technology – Active protection system: visual and acoustic temperature control</li> <li>• Machine washable at 30 °C</li> <li>• 2 side detachable controller</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• Pleasantly soft upper material</li> <li>• Automatic switch-off after 180 minutes</li> <li>• Handy controller</li> <li>• Suitable for mattress sizes 200×90-100 cm</li> <li>• Oeko-Tex® Standard 100</li> </ul>	<p><b>Technical data</b></p> <ul style="list-style-type: none"> <li>• 2 warming zones for body and feet</li> <li>• Automatic switch-off after 180 minutes</li> <li>• Ergonomically designed controller with LED control display on both sides for right and left-handed persons</li> <li>• Washable thanks to detachable controller</li> <li>• Pleasantly soft upper material</li> <li>• Suitable for mattress sizes 200×90-100 cm</li> <li>• Oeko-Tex® Standard 100</li> </ul>
---	--	---



<p><b>Size:</b> approx. 150×160 cm <b>Weight:</b> approx. 1.5 kg <b>Power supply:</b> 230 V/50 Hz <b>Power output:</b> 2 ×100 watts <b>Article No.:</b> 61230 <b>EAN:</b> 40 15588 61230 5 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 80×150 cm <b>Weight:</b> approx. 1.1 kg <b>Power supply:</b> 230 V/50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 60219 <b>EAN:</b> 40 15588 60219 1 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 150×80 cm <b>Weight:</b> approx. 0.82 kg <b>Power supply:</b> 230 V/50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 61232 <b>EAN:</b> 40 15588 61232 9 <b>Warranty:</b> 3 Years</p>
---	---	--

<p><b>Size:</b> approx. 80×150 cm <b>Weight:</b> approx. 1 kg <b>Power supply:</b> 230V 50Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 60218 <b>EAN:</b> 40 15588 60218 4 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 150×80 cm <b>Weight:</b> approx. 1 kg <b>Power supply:</b> 230V~ 50 Hz <b>Power output:</b> 60 watts <b>Article No.:</b> 60217 <b>EAN:</b> 40 15588 60217 7 <b>Warranty:</b> 3 Years</p>	<p><b>Size:</b> approx. 150×80 cm <b>Weight:</b> approx. 1 kg <b>Power supply:</b> 230V~ 50 Hz <b>Power output:</b> 100 watts <b>Article No.:</b> 61220 <b>EAN:</b> 40 15588 61220 6 <b>Warranty:</b> 3 Years</p>
---	--	---



**POWERROLL CLASSIC**  
Fascia roll

**Unique surface structure**

- Vibrating effect to reach deep-seated muscle tissues
- New surface structure to support punctual self-massage and increase blood circulation
- For activation and regeneration
- 8 intensity levels
- Rechargeable Lithium-Ion battery (incl. charger)



**Technical data**

- LED indication of the selected intensity level
- Surface is resistant to moisture
- Capacity information
- Max. operating time with fully charged battery: approx. 180 min. (depending on the selected intensity level)
- Auto-stop after 10 minutes
- Hardness grade: standard



Includes training plan

**Size:** approx. 31 cm long, ø 15 cm  
**Weight:** approx. 0.75 kg  
**Batteries:** Lithium-ion battery  
**Article No.:** 79530  
**EAN:** 40 15588 79530 5  
**Warranty:** 3 Years



**POWERROLL XT**  
Fascia roll

**Extra long**

- Vibrating effect to reach deep-seated muscle tissues
- New surface structure to support punctual self-massage and increase blood circulation
- For activation and regeneration
- 8 intensity levels
- Rechargeable Lithium-Ion battery (incl. charger)
- Length: 45 cm



**Technical data**

- LED indication of the selected intensity level
- Surface structure supports the massage effect
- Surface is resistant to moisture
- Capacity information
- Max. operating time with fully charged battery: approx. 180 min. (depending on the selected intensity level)
- Auto-stop after 10 minutes
- Hardness grade: standard



Includes training plan

**Size:** approx. 45 cm long, ø 15 cm  
**Weight:** approx. 0.9 kg  
**Batteries:** Lithium-ion battery  
**Article No.:** 79532  
**EAN:** 40 15588 79532 9  
**Warranty:** 3 Years



**POWERROLL XT ULTRASOFT**  
Fascia roll

**Hardness grade ultrasoft**

- Vibrating effect to reach deep-seated muscle tissues
- For activation and regeneration
- Surface structure supports punctual self-massage and increase blood circulation
- 8 intensity levels
- Length: 45 cm
- Hardness grade: Ultrasoft



**Technical data**

- LED indication of the selected intensity level
- Surface structure supports the massage effect
- Surface is resistant to moisture
- Capacity information
- Max. operating time with fully charged battery: approx. 180 min. (depending on the selected intensity level)
- Auto-stop after 10 minutes



Includes training plan

**Size:** approx. 45 cm long, ø 15 cm  
**Weight:** approx. 0.9 kg  
**Batteries:** Lithium-ion battery  
**Article No.:** 99078  
**EAN:** 40 15588 99078 6  
**Warranty:** 3 Years



**VARIOROLL**  
Foam roller


**Ultra flexible for any body position**

- Particular suitable for punctual self-massage – for example on the back, neck, etc.
- For the targeted use on fascia, trigger points and muscles
- Specially soft on the spine, thanks to 3-step size adjustment – distance between the balls (20 - 24 cm)
- For extended applications such as the arms, legs etc.



**Technical data**

- Resistant to moisture



Includes training plan

**Size:** approx. 20×10×10 cm, ø 15 cm  
**Weight:** approx. 0.13 kg  
**Article No.:** 79516  
**EAN:** 40 15588 79516 9  
**Warranty:** 3 Years





**EMS**  
EMS body trainer

**Targeted muscle building**

- Modern, sporty design
- Helps to build muscles
- Auto programs: Massage, acupuncture, tapping, combination
- Set including 3 controllers with 1 pad for the abdomen and 2 pads for arms or legs
- Automatic switch-off after 23 min.



**Technical data**

- Intensity adjustable in 15 steps
- PU leather pads
- Incl. 3 batteries (CR2032)



Includes 2 pads for arms and legs

**Size:** 17×18.5 cm (abdomen pad)  
5×20 cm (arm pad)  
**Weight:** 30 g (abdomen pad)  
13 g (arm pad)  
**Batteries:** 3×CR2032  
**Article No.:** 79524  
**EAN:** 40 15588 79524 4  
**Warranty:** 3 Years



**AM 880**  
Abdominal muscle stimulator

**50 intensity levels**

- Flexible neoprene waist belt with incorporated electrodes
- Choice of 6 abdominal muscle training programs (3 à 25 min, 3 à 30 min)
- Intensity can be set over 50 levels
- Can be worn under clothing thanks to the flat design



**Technical data**

- Automatic switch-off at the end of the program
- Simple operation with a large LCD display
- Device is integrated into the belt
- 1 channel with 4 electrodes
- Waist belt with Velcro closer
- Incl. extension belt for waist circumference up to 140 cm
- Water contact electrodes, no need for contact gel

**Size:** approx. 110×16×3 cm  
**Weight:** approx. 0.37 kg  
**Batteries:** 3×1.5V (AAA)  
**Curve form:** Bi-phase square wave pulse  
**Intensity:** 0–70mA  
**Frequency:** 30–90 Hz  
**Pulse width:** 200–300 µs  
**Article No.:** 88322  
**EAN:** 40 15588 88322 4  
**Warranty:** 3 Years



**Body trainer | Activity tracker**



**VIFIT RUN**  
Activity tracker

**With heart rate measurement**

- Tracks steps, distance, burned calories, activity time, sleep and heart rate function
- Individual daily target can be set - with reminder function
- Displays incoming calls or messages - SMS, WhatsApp, Facebook Messenger
- Real time heart rate dynamic measurement



**Technical data**

- Five heart rate training zones: Warm Up, Fat Burn, Aerobics, Cardio, High
- Sleep mode: tracks activity during sleep and hours slept
- Waterproof up to 50 metres
- Displays Date, time and battery status
- Alarm clock function with vibration alarm
- Display of the current day's weather
- Stores measurements of the last 30 days
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*
- Incl. USB cable

**Size:** approx. 4.15×2×1.18 cm  
**Weight:** approx. 0.032 kg  
**Batteries:** Lithium-ion battery  
**Article No.:** 79492  
**EAN:** 40 15588 79492 6  
**Warranty:** 3 Years



**VIFIT TOUCH**  
Activity tracker

**Activity and sleep tracker in one**

- High-quality OLED display
- Displays incoming calls or messages
- Stores the results of the last 15 days
- Individual daily goals and movement reminder
- Bluetooth® data transfer to the VitaDock+ app for iOS and Android and to VitaDock Online\*



**Technical data**

- Activity mode: tracks steps, distance, burned calories, activity time in hours and daily target in %
- Sleep mode: tracks movement activity during sleep and duration of sleep
- Motion sensor for activity measurement: counts actual steps
- Splash proof
- Vibration alarm clock
- Battery life: 5 - 7 days
- Rechargeable Lithium-Ion battery
- Display of time, date and battery indicator
- Incl. USB cable

**Size:** approx. 2.3×2×0.9 cm  
**Weight:** approx. 0.022 kg  
**Batteries:** Lithium-ion battery  
**Article No.:** 79486  
**EAN:** 40 15588 79486 5  
**Warranty:** 3 Years



\* iOS: iPhone 6S and newer, iPad Air 2, iPad Mini 4 and newer. Android: devices, which support Google Android version 6.x and the Bluetooth Smart 4.0





**MP 840**  
Manicure and pedicure device

**With 3 speed levels**

- Care and treatment of nails, cuticles and slight calluses
- With LED light for accurate application
- 3 speed levels for the individual requirements of care
- Stops automatically when the pressure is too high – also suitable for diabetics



**Technical data**

- Clockwise and counter-clockwise rotation
- High quality and durable attachments
- Stops automatically when the pressure is too high – also suitable for diabetics
- Natural, durable nail gloss (without additional care products)
- Simple to use, easy replacement of the attachments
- Incl. storage bag



Including 7 attachments: Sapphire cone (coarse), Sapphire disc (fine), Cylindrical burr, Flame burr, Felt cone, Needle burr as well as additional attachment for nail polishing

**Size:** approx. 15×3.8×3.7 cm  
**Weight:** approx. 0.1 kg  
**Power supply:** 100-240V~ 50/60 Hz, DC18V/500mA  
**Rotation speed:** 3 levels: 5.000/6.000/7.000  
**Power output:** 9 watts  
**Cable length:** 1.5 m  
**Article No.:** 85155  
**EAN:** 40 15588 85155 1  
**Warranty:** 3 Years



**MPS**  
Manicure and pedicure device

**With high-quality sapphire attachments**

- Care and treatment of nails, cuticles and slight calluses
- Continuous speed adjustment



**Technical data**

- Stops automatically when the pressure is too high – also suitable for diabetics
- Clockwise and counter-clockwise rotation
- High quality and durable sapphire attachments
- Incl. storage bag



Including 4 attachments: Sapphire cone (coarse), Sapphire disc (fine), Flame burr, Felt cone



**Size:** approx. 15×4×4 cm  
**Weight:** approx. 0.14 kg  
**Power supply:** Netzadapter: 100-240V~ 50/60 Hz, 18V DC 500 mA  
**Rotation speed:** 1600-5000 U/min  
**Article No.:** 85126  
**EAN:** 40 15588 85126 1  
**Warranty:** 3 Years



**MP 815**  
Manicure and pedicure device

**For individual needs of care**

- Care and treatment of nails, cuticles and slight calluses
- 2 speed levels
- High quality and durable attachments



**Technical data**

- Clockwise and counter-clockwise rotation
- Simple to use, easy replacement of the attachments
- Incl. storage bage



Including 7 attachments: Sapphire disc (fine), Sapphire disc (coarse), Sapphire callus remover, Sapphire cone, Felt cone, Cylindrical burr and Flame burr

**Size:** approx. 15×3.8×3.7 cm  
**Weight:** approx. 0.1 kg  
**Power supply:** 100-240V~ 50/60 Hz, 30V DC 250mA  
**Rotation speed:** 2 levels: 3.800/4.600  
**Cable length:** 1.5 m  
**Article No.:** 85153  
**EAN:** 40 15588 85153 7  
**Warranty:** 3 Years



**MP 810**  
Manicure and pedicure device

**Light weight and powerful**

- Care and treatment of nails, cuticles and calluses
- Easy to use and easy replacement of the attachments
- Remove cuticles, file, shape, smooth and polish nails




**Technical data**

- 5 attachments
- Compact device, requires little space
- Small and practical storage bag



Incl. 5 attachments: Felt cone, Trapezoidal grinder, Cylindrical grinder, Grinding cone, Flame burr

**Size:** approx. 16.5×3×3 cm  
**Weight:** approx. 0.8 kg  
**Batteries:** 2×1.5V (AAA)  
**Article No.:** 85140  
**EAN:** 40 15588 85140 7  
**Warranty:** 3 Years



**Hand & Foot care**



**NP 860**  
Electric nail polisher

**Long lasting shine without chemicals**

- For glossy, natural-looking nails
- Professional nail care at home
- Durable nail gloss (without additional care products)



**Technical data**

- Simple to use on fingernails and toenails
- Small and handy – ideal for on the go
- 2 speed levels for individual needs of nail care
- Ergonomic design
- Smoothing and polishing – in two steps to nails with “wow” effect



Including 2 smoothing and polishing rollers

**Size:** approx. 15.5×3.6×3.3 cm  
**Weight:** approx. 0.15 kg (incl. batteries)  
**Batteries:** 2×1.5V (AAA)  
**Rotation speed:** 2 levels: 2.000/2.400 U/min  
**Article No.:** 85160  
**EAN:** 40 15588 85160 5  
**Warranty:** 3 Years





**ND 870**  
UV nail hardener

**Professional nail design**

- Professional nail design for home use
- Ideal for the hardening of gel and liquid acrylic
- 360 degree reflective coating for an even curing
- Easy operation thanks to 4-step timer (90sec / 120sec / 30min / 60min)



**Technical data**

- Detachable floor section for an easy lamp replacement
- 4 non-slip rubber feet for a firm standing


**Size:** approx. 25.3×22.1×11.3 cm  
**Weight:** approx. 1.036 kg  
**Power supply:** 220-240V~ 50/60 Hz; 36W  
**Bulb wave length:** λ=365 nm  
**Article No.:** 85470  
**EAN:** 40 15588 85470 5  
**Warranty:** 3 Years



**CR 870**  
Electric callus remover

**Ergonomic handling**

- Care and treatment of nails, cuticles and callus
- Gentle and effective remove excess calluses
- For smoother hands, feet and elbows



**Technical data**

- 3 removable attachments: extra coarse, coarse and fine abrasive roll
- 2 speed levels
- Apply to dry skin
- Ergonomic, non-slip handle
- Incl. cleaning brush and batteries



Incl. 3 removable attachments

**Size:** approx. 15.5×3.43×8.04 cm  
**Weight:** approx. 0.148 kg  
**Batteries:** 3×1.5V AAA  
**Rotation speed:** 2 levels: 1.000/1.600  
**Article No.:** 88576  
**EAN:** 40 15588 88576 1  
**Warranty:** 3 Years



**CM 850**  
LED suction mirror

**Extra large and bright**

- Strong adhesion on all smooth surfaces
- Bright LED illumination
- Diameter of mirror approx. 19 cm
- 5× magnification



**Technical data**

- Suction enhancement thanks to special coating
- Removable and swivelling
- Battery powered
- Telescope bracket

**Size:** approx. 20×10×20 cm  
**Weight:** approx. 0.3 kg  
**Diameter of mirror:** approx. ø 19 cm  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 88558  
**EAN:** 40 15588 88558 7  
**Warranty:** 3 Years



**CM 840**  
Zin1 Cosmetic mirror

**Our bestseller**

- Illuminated cosmetic mirror
- 2 swivelling mirrors: normal and 5× magnification
- Uniform ring illumination with 18 LED-lights
- Diameter of mirror approx. 13 cm



**Technical data**

- High-quality chrome finish
- Stable base

**Size:** approx. 20.5×13.5×31.5 cm  
**Weight:** approx. 0.95 kg  
**Diameter of mirror:** approx. ø 13 cm  
**Batteries:** 4×1.5V (AAA)  
**Article No.:** 88550  
**EAN:** 40 15588 88550 1  
**Warranty:** 3 Years



**CM 835**  
Cosmetic mirror

**Compact and stable**

- Illuminated cosmetic mirror
- 2 swivelling mirrors: Normal and 5× magnification
- Uniform ring illumination with 12 white LED-lights
- Diameter of mirror approx. 12 cm



**Technical data**

- High-quality chrome finish
- Stable base

**Size:** approx. 16.9×8×19 cm  
**Weight:** approx. 0.375 kg  
**Diameter of mirror:** approx. ø 12 cm  
**Batteries:** 3×1.5V (AAA)  
**Article No.:** 88554  
**EAN:** 40 15588 88554 9  
**Warranty:** 3 Years



**DC 300**  
LED light therapy

**2 light therapies in one device**

- Handy light therapy
- Blue light to improve acne
- Red light for skin inflammation



**Technical data**

- 2×5 min daily are sufficient
- Auto switch-off after 5 min.
- Red light: 630 nm
- Blue light: 415 nm
- Light and handy - also for travelling
- Skin compatibility dermatologically tested\*

**Size:** approx. 2.6×2.6×14.1 cm  
**Weight:** approx. 0.23 kg  
**Batteries:** 1× AA  
**Article No.:** 85180  
**EAN:** 40 15588 85180 3  
**Warranty:** 3 Years



**FB 885**  
Facial cleansing brush

**Cleansing with 4 speed levels**

- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products
- Rechargeable lithium-ion battery



**Technical data**

- Cleanses the skin in 60 seconds (Timer)
- 4 speed settings
- Automatic switch-off after 1 min.
- Including 4 cleansing attachments (Sensitive skin, regular skin, peeling, cosmetics sponge)
- Splash-proof
- Incl. protective cover, storage bag and charge base



Including 4 cleansing attachments (sensitive skin, regular skin, peeling, cosmetics sponge)

**Size:** approx. 18.6×5.0×2.7 cm  
**Weight:** approx. 0.27 kg  
**Batteries:** Lithium-ion battery  
**Power output:** 500 mAh  
**Cable length:** 1.30 m  
**Article No.:** 88565  
**EAN:** 40 15588 88565 5  
**Warranty:** 3 Years



**FB 880**  
Facial cleansing brush

**Soft cleansing with vibration**

- Facial cleansing brush for a gentle and thorough cleaning
- Frees the skin from dead skin particles
- Gentle vibration stimulates the blood circulation and promotes the absorption of care products



**Technical data**

- Cleanses the skin in 60 seconds (Timer)
- 2 speed settings: normal or intensive cleansing
- Automatic switch-off after 1 min.
- Including 3 cleansing attachments (sensitive skin, regular skin, peeling)
- Battery running time lasts up to 70 hours
- Splash-proof
- Incl. protective cover, storage bag and batteries (2 x 1.5V AA)



Including 3 cleansing attachments (sensitive skin, regular skin, peeling)

**Size:** approx. 15.7×5.7×7.9 cm  
**Weight:** approx. 0.22 kg  
**Batteries:** 2×1.5V (AA)  
**Article No.:** 88560  
**EAN:** 40 15588 88560 0  
**Warranty:** 3 Years



**FSS**  
Facial sauna

**Soothing with aroma attachment**

- Opens the pores, cleans the skin and relaxes facial lines with humid heat
- Additional attachment for soothing aroma therapies
- Easy cleaning and handling



**Technical data**

- For cosmetic use
- Diffuser for even steam distribution
- Additional nose piece
- High-quality heating element
- Indicator light
- Measuring cup for easy filling



Including nose piece and measuring cup

**Size:** approx. 30×20×19.5 cm  
**Weight:** approx. 0.56 kg  
**Power supply:** 220-240V~ 50/60 Hz, max. 100 watts  
**Power output:** 88245  
**Article No.:** 88245  
**EAN:** 40 15588 88245 6  
**Warranty:** 3 Years



**AC 850**  
Cellulite massager

**Intensive rolling massage**

- Effective massage for stimulating the connective tissue for tangibly firmer skin
- Improves blood circulation
- 2 selectable massage intensities
- Light and handy design for convenient self-massage



**Technical data**

- 3×2 rotating massage rollers (removable)
- Individually adjustable hand strap
- Especially powerful thanks to mains operation
- Easy cleaning of the removable massage rollers

**Size:** approx. 13×11×8.5 cm  
**Weight:** approx. 0.35 kg  
**Power supply:** 110-240V~ 50/60 Hz; 12V DC 1200 mA  
**Power output:** 14.4 watts  
**Cable length:** 2 m  
**Article No.:** 88540  
**EAN:** 40 15588 88540 2  
**Warranty:** 3 Years

\*Dermatological clinical study on 20 participants





**UPPER ARM CUFFS**  
for blood pressure monitors

**MTM | MTP**

Large Cuff (32 – 42 cm)  
Article No.: 51037  
EAN: 40 15588 51037 3

Normal cuff (22 – 32 cm)  
Article No.: 51038 (included in package)  
EAN: 40 15588 51038 0

Small cuff (17 – 22 cm)  
Article No.: 51034  
EAN: 40 15588 51034 2

**BU 535 | BU 535 Voice | MTV | MTC | MTD | MTS | MTR | CardioCompact**

Large cuff (30 – 42 cm)  
Article No.: 51126  
EAN: 40 15588 51126 4

Normal cuff (22 – 30 cm)  
Article No.: 51135 (included in package)  
EAN: 40 15588 51135 6

**MTX**

Large „Easy-Form“ cuff (22 – 42 cm)  
Article No.: 51087  
EAN: 40 15588 51087 8

Normal „Easy-Form“ cuff (22 – 32 cm)  
Article No.: 51088  
EAN: 40 15588 51088 5

**BU 550 Connect | BU 575 connect**

Normal cuff (22 – 42 cm)  
Article No.: 51299  
EAN: 40 15588 51299 5

**BU 510 | BU 520 | BU 530 connect**

Large cuff (32 – 42 cm)  
Article No.: 51169  
EAN: 40 15588 51169 1

Normal cuff (22 – 36 cm)  
Article No.: 51168 (included in package)  
EAN: 40 15588 51168 4

**CardioDock | CardioDock 2**

Large „Easy-Form“ cuff (30 – 42 cm)  
Article No.: 51284  
EAN: 40 15588 51284 1

Normal „Easy-Form“ cuff (22 – 32 cm)  
Article No.: 51283  
EAN: 40 15588 51283 4



**MAINS ADAPTER**  
for blood pressure monitors

**Mains adapter for BU 550 connect**

Article No.: 51013  
EAN: 40 15588 51013 7

**Mains adapter for BU 575 connect | BW 300 connect**

Article No.: 51014  
EAN: 40 15588 51014 7

**Mains adapter for MTM | MTP | MTX**

Article No.: 51036  
EAN: 40 15588 51036 6

**Mains adapter for CardioCompact**

Article No.: 51095  
EAN: 40 15588 51095 3

**Mains adapter for MTV | MTC | MTD | MTS | BU 510 connect | BU 520 | BU 530 Connect**

Article No.: 51125  
EAN: 40 15588 51125 7



**PROTECTIVE SLEEVES FOR**  
for thermometer

**FTF**

50 pieces  
Article No.: 76111  
EAN: 40 15588 76111 9



**ACCESSORIES**  
for blood glucose meters

**Control solution for MediTouch 2 | MediTouch 2 connect | MediTouch 2 Carbon**

Article No.: 79039  
EAN: 40 15588 79039 3

**Lancets for MediTouch 2 | MediTouch 2 connect | MediTouch 2 Carbon**

Article No.: 79028  
EAN: 40 15588 79028 7

**Test strips for MediTouch 2 | MediTouch 2 connect | MediTouch 2 Carbon**

Article No.: 79038  
EAN: 40 15588 79038 6

**Test strips for Gluowatch connect**

Test strips 50 and 10 lancets  
Article No.: 79152  
EAN: 40 15588 79152 9

**Lancing device**

Article No.: 79002  
EAN: 40 15588 79002 7

**USB cable for Meditouch 2**

Article No.: 79049  
EAN: 40 15588 79049 2



**WRISTBANDS**  
for activity tracker

**ViFit connect**

**Black**  
Article No.: 79417  
EAN: 40 15588 79417 9

**Orange**  
Article No.: 79418  
EAN: 40 15588 79418 6

**Blue**  
Article No.: 79470  
EAN: 40 15588 79470 4

**Green**  
Article No.: 79471  
EAN: 40 15588 79471 1

**Red**  
Article No.: 79472  
EAN: 40 15588 79472 8

**Pink**  
Article No.: 79473  
EAN: 40 15588 79473 5



**SPARE PART SETS**  
for EMS Pad

**Body Trainer**

Spare part sets (2 pads for the abdomen)  
Article No.: 79525  
EAN: 40 15588 79525 1

Spare part sets (4 pads for arms / legs)  
Article No.: 79526  
EAN: 40 15588 79526 8



**ACCESSORIES**  
for bath spas

**BBS | MBH**

**Air hose**  
Article No.: 30666  
EAN: 40 15588 30666 2

**Soft mat**  
Article No.: 88383  
EAN: 40 15588 88383 5

**Fragrance dispenser with 3 pads**  
Article No.: 30669  
EAN: 40 15588 30669 3

**BBS**

**Remote control**  
Article No.: 88384  
EAN: 40 15588 88384 2

**MBH**

**Remote control**  
Article No.: 88376  
EAN: 40 15588 88376 7



**LOOFAH MASSAGE PAD**  
for mini handmassager

**HM 840**

Article No.: 88526  
EAN: 40 15588 88526 6



**ELECTRODES**

for pain therapy devices

**TDD | TDP | TED**  
 4 Multiple electrodes (approx. 3 x 5 cm)  
 Article No.: 87080  
 EAN: 40 15588 87080 4

2 Multiple electrodes (approx. 5 x 9 cm)  
 Article No.: 87100  
 EAN: 40 15588 87100 9

Note: Not suitable for hair removal system EPD (Article No.: 85550)

Permanent electrodes  
 Article No.: 13586  
 EAN: 40 15588 13586 6

**TDP**  
 Cable set  
 Article No.: 88315  
 EAN: 40 15588 88315 6

**BT 850**  
 Adhesive pads (2 pcs.)  
 Article No.: 88334  
 EAN: 40 15588 88334 7

**TDB**  
 Replacement electrodes  
 Article No.: 86359  
 EAN: 40 15588 86359 2



**REPLACEMENT LAMP**

for infrared lamps

**IRH**  
 100 watts  
 Article No.: 30101  
 EAN: 40 15588 30101 8

**IRL**  
 150 watts  
 Article No.: 88255  
 EAN: 40 15588 88255 5

**IR 100**  
 100 watts  
 Article No.: 30101  
 EAN: 4015588301018



**ACCESSORIES**

for Inhalers

**IN 500 | IN 550**  
 Connecting piece  
 Article No.: 54524  
 EAN: 40 15588 54524 5

Air filter and nebulizer kit  
 Article No.: 54525  
 EAN: 40 15588 54525 2

Mouthpiece and nose adapter  
 Article No.: 54526  
 EAN: 40 15588 54526 9

Air hose  
 Article No.: 54527  
 EAN: 40 15588 54527 6

Mask for children  
 Article No.: 54528  
 EAN: 40 15588 54528 3

Mask for adults  
 Article No.: 54529  
 EAN: 40 15588 54529 0

Nebulizer attachments (3 pcs.)  
 Article No.: 54531  
 EAN: 40 15588 54531 3

**IN 525**  
 Nebulizer kit  
 Article No.: 54116  
 EAN: 40 15588 54116 2

Mask set (children + adults) with mouthpiece  
 Article No.: 54117  
 EAN: 40 15588 54117 9

Mains adapter EU  
 Article No.: 54118  
 EAN: 40 15588 54118 6

**USC**  
 Mask set (2 masks + adapter)  
 Article No.: 54101  
 EAN: 40 15588 54101 8

Mouthpiece  
 Article No.: 54102  
 EAN: 40 15588 54102 5

Covering cap  
 Article No.: 54103  
 EAN: 40 15588 54103 2

Nebulizer kit  
 Article No.: 54104  
 EAN: 40 15588 54104 9



**AROMA ESSENCES**

for air humidifiers/Aroma diffusers

Aroma Lemon 10 ml  
 Article No.: 60030  
 EAN: 40 15588 60030 2

Aroma Eukalyptus 10 ml  
 Article No.: 60031  
 EAN: 40 15588 60031 9

Aroma Lavender 10 ml  
 Article No.: 60032  
 EAN: 40 15588 60032 6

Aroma Citronella 10 ml  
 Article No.: 60033  
 EAN: 40 15588 60033 3

Aroma Ylang-Ylang 10 ml  
 Article No.: 60035  
 EAN: 40 15588 60035 7

Aroma Peppermint 10 ml  
 Article No.: 60036  
 EAN: 40 15588 60036 4

Aroma Orange 10 ml  
 Article No.: 60037  
 EAN: 40 15588 60037 1

Aroma Pinewood 10 ml  
 Article No.: 60038  
 EAN: 40 15588 60038 8

Aroma Wildberries 10 ml  
 Article No.: 60039  
 EAN: 40 15588 60039 5

Aroma Sea breeze 10 ml  
 Article No.: 60040  
 EAN: 40 15588 60040 1



**REPLACEMENT FILTER**

for air purifier

**MEDISANA AIR®**  
 1 pc.  
 Article No.: 60390  
 EAN: 40 15588 60390 7

2 pcs.  
 Article No.: 60391  
 EAN: 40 15588 60391 4



**REPLACEMENT WATER FILTER**

for air humidifier

**AH 660 | AH 665**  
 2 pcs.  
 Article No.: 60069 | 60059  
 EAN: 40 15588 60069 2

**AH 670**  
 1 pc.  
 Article No.: 60057  
 EAN: 40 15588 60057 9



**REPLACEMENT BRUSHES**

for facial cleansing brushes

**FB 880**  
 3 pcs.  
 Article No.: 88561  
 EAN: 40 15588 88561 7

**FB 885**  
 4 pcs.  
 Article No.: 88566  
 EAN: 40 15588 88566 2



**REPLACEMENT CARTRIDGE**

for hair removal system

**IPL 800**  
 Article No.: 88581  
 EAN: 40 15588 88581 5

**IPL 805**  
 Article No.: 88584  
 EAN: 40 15588 88584 6

**IPL 836 | IPL 840**  
 Article No.: 88586  
 EAN: 40 15588 88586 0

**IPL 850**  
 Article No.: 88589  
 EAN: 40 15588 88589 1



**REPLACEMENT ROLLS**

for electric callus remover

**CR 860**  
 2 pcs.  
 Article No.: 88572  
 EAN: 40 15588 88572 3



**REPLACEMENT ROLLS**

for manicure and pedicure devices

**NP 860 | MP 840**  
 4 pcs.  
 Article No.: 85164  
 EAN: 40 15588 85164 3



**ATTACHMENTS**

for manicure and pedicure devices

**MANILUX L | MANILUX S | MP 820**  
 Set of attachments (10 pieces)\*  
 Article No.: 85450  
 EAN: 40 15588 85450 7

\*contains: sapphire cone (coarse), sapphire disc (coarse), flame burr, sapphire ball burr, felt cone, sapphire cone (long), sapphire disc (fine), needle burr (long), cylindrical burr, drainage grinder

**MEDISTYLE S | MEDISTYLE L | MPS | MP 840**

Felt cone  
 Article No.: 85001  
 EAN: 40 15588 85001 1

Sapphire cone (coarse)  
 Article No.: 85003  
 EAN: 40 15588 85003 5

Sapphire disc (fine)  
 Article No.: 85004  
 EAN: 40 15588 85004 2

Sapphire ball burr  
 Article No.: 85005  
 EAN: 40 15588 85005 9

Needle burr (long)  
 Article No.: 85006  
 EAN: 40 15588 85006 6

Sapphire cone (long)  
 Article No.: 85007  
 EAN: 40 15588 85007 3

Cuticle pusher  
 Article No.: 85008  
 EAN: 40 15588 85008 0

Cylindrical burr  
 Article No.: 85009  
 EAN: 40 15588 85009 7

Sapphire disc (coarse)  
 Article No.: 85010  
 EAN: 40 15588 85010 3

Flame burr  
 Article No.: 85011  
 EAN: 40 15588 85011 0

Attachment adapter  
 Article No.: 85015  
 EAN: 40 15588 85015 8

Grinding set for artificial nails  
 Article No.: 85460  
 EAN: 40 15588 85460 6

Sanding Bands (30 pieces)  
 Article No.: 85461  
 EAN: 40 15588 85461 3





**HEAD OFFICE**

**medisana GmbH**  
 Jagenbergstr. 19, 41468 Neuss, Germany  
 Phone: +49 2131 36 68-0  
 Fax: +49 2131 36 68-50 95  
 eMail: [info@medisana.de](mailto:info@medisana.de)  
 Internet: [www.medisana.de](http://www.medisana.de)

Our head office is open from Monday to Thursday  
 8:00 am – 4:00 pm, Friday 8:00 am – 2:00 pm.  
 Outside these hours, you can leave a message on  
 our voice mail.

**SERVICE CENTER**

Phone: +49 2131 405 88-0  
 eMail: [servicecenter@medisana.de](mailto:servicecenter@medisana.de)

**SUBSIDIARIES**

**GERMANY**

**medisana Space Technologies GmbH**  
 Jagenbergstr. 19  
 41468 Neuss  
 Phone: +49 2131 36 68 0  
 eMail: [info@medisanaspace.com](mailto:info@medisanaspace.com)  
 Internet: [www.medisana.de](http://www.medisana.de)

**AUSTRIA**

**medisana Austria**  
 Wiener Str. 35a  
 A-8605 Kapfenberg  
 Phone: +43 386 22 10 52  
 Fax: +43 316 23 11 23 52 12  
 eMail: [office@medisana.at](mailto:office@medisana.at)  
 Internet: [www.medisana.at](http://www.medisana.at)  
 Contact: Thorsten Hüning

**medisana Service Austria**

Burgstaller Strasse 2  
 A-8143 Dobl  
 Phone: +43 313 65 20 76  
 eMail: [service.at@medisana.com](mailto:service.at@medisana.com)  
 Contact: Christian Eichholzer

**FRANCE**

**medisana GmbH**  
 Jagenbergstr. 19  
 D-41468 Neuss  
 Phone: +49 2131/3 66 80  
 eMail: [medisana.export@medisana.de](mailto:medisana.export@medisana.de)  
 Internet: [www.medisana.fr](http://www.medisana.fr)  
 Contact: Adem Otlu

**NETHERLANDS/BELGIUM/LUXEMBOURG**

**medisana Benelux NV**  
 Euregiopark 30  
 NL-6467 JE Kerkrade  
 Phone: +31 45 547 08 60  
 Fax: +31 45 547 08 79  
 eMail: [info@medisana.nl](mailto:info@medisana.nl)  
 Internet: [www.medisana.nl](http://www.medisana.nl)  
 Contact: Maurice Jongen

**medisana Service Niederlande**

Teknihall Benelux  
 p/a Vreysen  
 Hazeldonk 6027  
 4836 LA BREDA  
 Phone: +31 45 547 08 60  
 eMail: [info@medisana.nl](mailto:info@medisana.nl)  
 Internet: [www.medisana.nl](http://www.medisana.nl)

**medisana Service Belgien/Luxemburg**

Teknihall bvba  
 Brusselstraat 33  
 2321 BE MEER  
 Phone: +31 45 547 08 60  
 eMail: [info@medisana.be](mailto:info@medisana.be)  
 Internet: [www.medisana.be](http://www.medisana.be)

**ISRAEL**

**medisana Israel**  
 13 Zarchin Street  
 Building A  
 Raanana 43662 Israel  
 Phone: +972 9 7711418  
 Fax: +972 9 7740131  
 eMail: [a.cohen@medisana-global.com](mailto:a.cohen@medisana-global.com)  
 Internet: [www.medisana.com](http://www.medisana.com)  
 Contact: Avi Cohen

**ASIA**

**medisana Far East Ltd.**  
 Unit E, F, 31/F, YHC Tower  
 1 Sheung Yuet Road, Kowloon Bay,  
 Kowloon, Hong Kong  
 Phone: +852 23 63 17 78  
 Fax: +852 23 63 19 55  
 eMail: [info@medisana.com.hk](mailto:info@medisana.com.hk)  
 Contact: Kitty Yeung

**UNITED ARAB EMIRATES**

**medisana MENA FZ-LLC.**  
 Dubai Science Park  
 North Tower, Office 705  
 Phone: +971 50 476 8055  
 eMail: [w.luka@medisana-global.com](mailto:w.luka@medisana-global.com)  
 Internet: [www.medisana.com](http://www.medisana.com)  
 Contact: Walid Luka

**SPAIN**

**medisana Healthcare S. L.**  
 C/Clementina Arderiu, Nave 3  
 Pol. Industrial Les Vives  
 ES-08295 Sant Vicens de Castellet/Barcelona  
 Phone: +34 938 33 50 52  
 Fax: +34 938 33 36 22  
 eMail: [info@medisana.es](mailto:info@medisana.es)  
 Internet: [www.medisana.es](http://www.medisana.es)  
 Contact: Cornelius Hellmund

**GREECE**

**medisana Hellas L.t.d**  
 11 Helmou Str.  
 GR-14122 N. Iraklion – Athens T.K  
 Phone: +30 210 / 275 09 32  
 Fax: +30 210 / 275 00 19  
 eMail: [info@medisana.gr](mailto:info@medisana.gr)  
 Internet: [www.medisana.gr](http://www.medisana.gr)  
 Contact: Petros C. Koutras

**RUSSIA**

**000 medisana Rus**  
 ul. Nagornaja d.20, korp. 1, komn. 6  
 RU-117186 Moskau  
 Phone: +7 495 640 26 47  
 Fax: +7 495 640 26 47  
 eMail: [info@medisana.su](mailto:info@medisana.su)  
 Internet: [www.medisana.su](http://www.medisana.su)  
 Contact: Andrei Usov

**GERMANY**

**Promed GmbH**  
 Kosmetische Erzeugnisse  
 Lindenweg 11  
 D-82490 Farchant  
 Phone: +49 8821 9621-0  
 Fax: +49 8821 9621-21  
 eMail: [info@promed.de](mailto:info@promed.de)  
 Internet: [www.promed.de](http://www.promed.de)  
 Contact: Andreas Schuster

**Imprint**

medisana GmbH  
 Jagenbergstr. 19, 41468 Neuss, Deutschland

Phone: +49 2131 / 36 68-0  
 Fax: +49 2131 / 36 68-50 95

eMail: [info@medisana.de](mailto:info@medisana.de)  
 Internet: [www.medisana.de](http://www.medisana.de)

General manager: Marco Getz  
 Kwok Shing Ng

Registration court: Neuss County Court  
 Register number: HRB 16348  
 VAT ID number: DE 123375019

**Modifications / Terms and Conditions**

Subject to technical modifications and further developments,  
 also without prior notice. We accept no liability for printing  
 errors, typesetting errors or other mistakes. Our current Terms  
 and Conditions can be found at any time on our website.

**Copyright**

Copyright over the information published in the catalogue  
 (text, pictures, graphics) and their structure is owned by  
 medisana GmbH, and in some cases by third parties. Duplication,  
 distribution, reproduction, dissemination, storage and  
 other use of this information without our written consent is  
 expressly prohibited.





# medisana<sup>®</sup>

YOUR HEALTH IN GOOD HANDS

Jagenbergstraße 19  
41468 Neuss  
Germany

Phone: +49 2131 36 68-0  
Fax: +49 2131 36 68-50 95

eMail: [info@medisana.de](mailto:info@medisana.de)  
Internet: [www.medisana.de](http://www.medisana.de)



September 2019